

CHEM 130
Chemical Principles I
Laboratory Manual

Copyright 2016 by Truman State University. All rights reserved.

Table of Contents

Determination of Density.....	2
Preparation and Analysis of Alum.....	11
Hydrolysis of Oil of Wintergreen	29
Vitamin C Analysis.....	43
Enthalpies of Solution.....	49
Selected Enthalpies of Formation	59
Using Enthalpies of Formation to Calculate Enthalpies of Reaction	61
The First Law and the Sign Convention used in Thermodynamics.....	62
Kinetics of Crystal Violet Bleaching	65
Spectrophotometric Determination of an Equilibrium Constant	71
Equipment and Techniques.....	79
Operation of the Vernier LabQuest 2.....	80
Basic Data Collection Using Logger Pro.....	85
Basic Data Analysis Using Logger Pro	88
Operating Ocean Optics and Vernier Spectrometers with LoggerPro.....	90
pH Measurements with Logger Pro	98
Operating the Spectronic 20 Genesys Spectrometer.....	101
Melting Point Determination Apparatus.....	104
Data Analysis, Recordkeeping and Reporting Results	107
What You need to do Before Coming to Laboratory.....	108
The Laboratory Notebook.....	109
Introduction to Statistics in Chemistry	116
Preparing Graphs	124
Guide to Excel.....	128
Advanced Regression with Microsoft Excel.....	137
Propagation of Uncertainty	141
The Laboratory Report.....	147
Example Formal Report.....	161

Determination of Density

B. D. Lamp, D. L. McCurdy, V. M. Pultz and J. M. McCormick*

Introduction

Not so long ago a statistical data analysis of any data set larger than a few points was a time-consuming and tedious procedure.^{1, 2} This was changed first by the introduction of personal computers and then by spreadsheets, which are computer programs that allow the user to enter and manipulate numerical data. Spreadsheets were originally designed for business applications, but have become essential tools for data analysis in all of the sciences, because of the ease with which they can perform complex calculations, and graph the results. Many hand-held calculators can perform similar tasks, but spreadsheets have the advantage because they store data in an easily edited form and produce higher-quality graphs.

In this exercise you will learn the basics of statistical data analysis and of spreadsheet operations using the program Excel. The data that you will manipulate will be measured values of copper's density obtained by first measuring a copper block's dimensions and then by water displacement. Before reading this exercise and preparing your notebook, read [*The Laboratory Notebook*](#), the [*What You need to do Before Coming to Laboratory*](#), [*Introduction to Statistics in Chemistry*](#), [*Preparing Graphs*](#) and [*Guide to Excel*](#) sections of the laboratory manual. You will want to consult these documents as you prepare your notebook for this exercise. Save your work frequently.

Experimental

This laboratory has no particular safety hazards associated with it, as long as normal laboratory protocols are followed.

IMPORTANT! You are required to bring either a laptop or a tablet with Microsoft Excel installed on it to the laboratory this week. Note that the screen shots in this experiment are from Excel 2010, but later versions of Excel are essentially the same. Other similar spreadsheet programs are acceptable, but the screen shots shown in this exercise and the some of the functions that we will use may be different in a different program. The "help" function in Excel or other online resources are useful if you run into problems locating specific features.

Before coming to lab, prepare tables similar to Table 1 and Table 2, shown above, in your laboratory notebook and label them as shown. Leave enough space so that each table has twelve blank rows for data (there may be up to twelve groups in your laboratory). Note how each table organizes the data in an easily read and understood format.

Table 1. Partial view of the table used to record the dimensions of the copper blocks and each block's mass in the laboratory notebook. Your table may have to accommodate data for up to twenty four blocks, depending on your instructor's wishes.

Block Number	Length (cm)	Width (cm)	Height (cm)	Mass (g)	Volume (cm ³)	Uncertainty in the Volume (cm ³)

Table 2. Partial table for recording the class data for copper's density as determined by water displacement. Again, the table in your notebook should accommodate data for as many as 24 blocks.

Block Number	Mass (g)	Volume H ₂ O Displaced (mL)	Density (g/cm ³)	Uncertainty in the Density (g/cm ³)

You will be assigned a copper block; write down the number of your block in your notebook, and make all of your measurements on the same block. Describe the block's color, texture and appearance in your notebook's *Results* section being as descriptive as possible.

Determination of Density using a Ruler to Measure the Volume

Obtain the mass of the block to three decimal places using one of the top-loading balances located in the laboratory. It is **never** good lab practice of set a chemical directly on a balance pan. Therefore, place a piece of weighing paper or a plastic weigh boat on the pan. Zero the balance by pressing the tare button and then place the copper block on the weighing paper or in the weigh boat. Record the block's mass in Table 1 and in Table 2, making sure that all three decimal places are recorded, even if some, or all, of them are 0. If a balance is not displaying three decimal places, or if the number of decimal places changes when you put your block on the balance, bring it to the attention of your instructor and he or she will assist you. We will assume that the uncertainty associated with the mass measurement (Δm) is ± 1 in the last decimal place measured (i. e., $\Delta m = \pm 0.001$ g for balances reading three decimal places).

Measure the length, width and height of your block using a plastic ruler. The plastic rulers are marked off every 0.1 cm, but you can estimate and report the measurements to ± 0.01 cm (see Fig. 1). We will assume that the copper pieces are perfect rectangular blocks (the four lengths are the same, as are the four widths and the four heights, and all sides meet at 90° angles).




Figure 1. How to estimate the length of an object to the nearest 0.01 cm with a ruler that is graduated in increments of 0.1 cm by dividing the distance between the gradations with ten imaginary lines (shown in red). The length of this object is thus 1.22 cm.

Determination of Density by Water Displacement

Add enough water to a 50-mL graduated cylinder so that your copper block will be completely submerged. The 50-mL graduated cylinder is marked in 1-mL increments, but you should be able to estimate the volume to the nearest 0.1 mL, as shown in Fig. 2. Record the starting volume of water in your notebook. Carefully place the copper into the graduated cylinder being careful not to splash any water out of the cylinder. Gently tap the cylinder to dislodge any air bubbles that are clinging to the copper, and record the new volume. Calculate the difference between the final and initial volumes to determine the volume of water displaced; enter this value in Table 2.




Figure 2. How to estimate the volume to the nearest 0.1 mL in a graduated cylinder that is marked with 1 mL gradations by dividing the distances between gradations into imaginary lines (shown in red). The volume in the case shown would be 18.6 mL.

Results and Analysis

Graphical Determination of Density using a Ruler to Measure the Volume

Calculate the block's volume, V , from its dimensions paying particular attention to the significant figures in your calculation. Determine the uncertainty in the volume, ΔV , from the uncertainties in each dimension (Δx , Δy , Δz) using Equation 1 (Eqn. 1). Your instructor will tell you how to share your data with the whole class.

$$\Delta V = \pm V \sqrt{\left(\frac{\Delta x}{x}\right)^2 + \left(\frac{\Delta y}{y}\right)^2 + \left(\frac{\Delta z}{z}\right)^2} \quad (1)$$

Calculate the density, d , of your block, and determine the uncertainty associated with this single measurement of the density, Δd , using Eqn. 2, and your value of ΔV , assuming that Δm is ± 0.001 g. Record these calculations in your laboratory notebook.

$$\Delta d = \pm d \sqrt{\left(\frac{\Delta V}{V}\right)^2 + \left(\frac{\Delta m}{m}\right)^2} \quad (2)$$

Prepare a graph of the copper blocks' volume as a function of their mass in Excel (see the [Guide to Excel](#) section). First, open Excel and set up the first work sheet as shown in Fig. 3.³ Enter the date, your name and your lab partner's in the cells where it says *Date* and *Names*, respectively. Enter the class data for the copper blocks in columns B through E starting with block 1 in row 4

and continuing on to block 12 in row 15. If not all twelve blocks were measured, leave those rows blank.

Be sure that all significant figures are shown! Spreadsheets drop trailing zeros, even if they are significant, and you will need to adjust the significant figures displayed in the cell using the increase and decrease decimal places buttons (helpful hint: select multiple cells before clicking on one of the adjust decimal places buttons, see the [Guide to Excel](#) section for instructions on how to do this, if needed).




Figure 3. Worksheet for the data from the determination of copper's density from the block's dimensions and its mass.

Enter a formula in cell F4 to calculate the volume of the block from its dimensions.⁴ Copy this cell and paste it into cells F5 through F15. Check to see that your calculated volumes are the same as those you and your classmates found. If there are mistakes, locate the errors' sources, and correct them.

Translate Eqn. 1 into Excel format and enter it in cell G4.⁵ Copy and paste cell G4 into cells G5 through G15. Once again, check these values and correct any mistakes.

In cell H4 enter a formula to calculate $\Delta V/V$ (the uncertainty in the volume measurement divided by the volume). Copy and paste cell H4 into cells H5 through H15. Set the number of significant figures in cells H4 through H15 to two (this is not the correct number of significant figures, but it will be useful later on).

Now that the data have been entered into the spreadsheet, we need to find the relationship between the block's mass (m) and its volume (V). We will assume that there is a linear relationship, which we can write as Eqn. 3, where a is the slope of the line and b is the y-intercept. To avoid confusion, the slope is given the symbol " a ", instead of the usual " m ",

$$V = a \cdot m + b \quad (3)$$

because the mass already has that symbol. Prepare a graph in Excel with volume as the dependent variable and mass as the independent variable.⁶ Insert a trend line in the graph so that the line's equation is displayed on the graph (see the [Guide to Excel](#) section for instructions on how to do this). Print out two copies of the graph so that each fills a half-sheet of paper. Attach one copy to an original page of your notebook and its mate to the corresponding duplicate page.⁷ Helpful hint: it is easier to copy the graph and paste it into *Word* before printing (see the [Guide to Excel](#) section).

Using Excel's regression package (see the [Guide to Excel](#) section), calculate the slope, the intercept and the uncertainty in the slope and intercept at the 95% confidence limit for these data. Be sure that the slope and intercept found using the regression package are the same as those determined from the trend line. If they are not exactly the same, there is a problem somewhere that you must correct before continuing. Adjust the width of the columns on the regression-output worksheets so that all of the headings can be read. Use the print set up/print preview options to print these sheets such that each fits legibly onto a half-sheet of paper (again, copying and pasting into *Word* may give the best results). Print out two copies so that you can attach one copy to an original and one to a duplicate page in your notebook. Be sure to save the Excel file and to record the name you saved it under in your laboratory notebook.

Write down in your notebook the final values of the slope and intercept, and give their 95% confidence interval. Watch your significant figures and units! Remember that the uncertainty is telling you the position of the last significant figure (see the [Introduction to Statistics in Chemistry](#) section).

Determine copper's density, d , from the slope of the best-fit line from your Excel graph and from your hand-drawn graph. Calculate the uncertainty in copper's density (Δd) at the 95% confidence limit from the standard deviation in the slope of your Excel graph (Δa , labeled "Standard Error" in the output for the regression package) using Eqn. 4. This must be done because the slope is $1/d$, and so Δa is not Δd . Write this value in your notebook using the proper format (see the [Introduction to Statistics in Chemistry](#) section). Calculate a percent error for the average value using the accepted density of copper (8.96 g/cm^3).⁸

$$\Delta d = \pm d \sqrt{\left(\frac{\Delta a}{a}\right)^2} \quad (4)$$

From this uncertainty, it is possible to determine the 95% confidence interval for our experimentally-determined density. When dealing with two-dimensional data sets, we must use a slightly different approach to calculate the confidence interval, as compared to one-dimensional data. In the two-dimensional case, the confidence limit is calculated as $\Delta = t \cdot \sigma$, where σ is the standard deviation and t is determined based on $n - 2$ degrees of freedom, where n is the number of x, y pairs in the dataset. (Notice the absence of the square root of n term, this is not a typo!) You can use either the TINV function (see [Guide to Excel](#) section or Table 1 in the [Introduction to Statistics in Chemistry](#)) to find the appropriate value for Student's t .




Figure 4. Worksheet for the data from the determination of copper's density by water displacement.

Determination of Density by Water Displacement

Prepare the second work sheet in your Excel spread sheet for the water displacement data so that it looks like that shown in Fig. 4.⁹

Enter the class data starting with block 1 in cell B5. Calculate the density of each block from its mass and volume by entering the correct formula in cells D5 through D16. Remember to adjust the number of decimal places in each to reflect the correct number of significant figures.

In cell D17, calculate the average density using the AVERAGE function (see [Guide to Excel](#) section). Adjust the number of significant figures in the calculated average (remember that the average can be no more precise than the least precise number used to calculate it).

Inspect the data to identify whether any point seems to be out of place. If you find a point that you think is an outlier, first check that there were no computational or other gross errors, then perform the *Q*-test ([Introduction to Statistics in Chemistry](#), Eqn. 6 and Table 2) on the suspect point. Record these calculations in your notebook. If you can eliminate a point, remove it completely from the spreadsheet. Write in your notebook beside Table 2 that this point was “eliminated on the basis of a *Q*-test.”

Calculate the estimated standard deviation, *S*, of the data using Eqn. 5, where *N* is the number of points in the data set, *x_i* is each individual measurement and *x_{avg}* is the average. First, we will do it in a step-wise calculation, and then use a built-in Excel function.

$$S = \sqrt{\frac{1}{N-1} \sum_{i=1}^N (x_i - x_{avg})^2} \quad (5)$$

In cell E5 enter a formula to calculate $(x_i - x_{avg})$ using an absolute reference to the cell which contains the average.¹⁰ This will stop Excel from automatically changing the reference to the cell with the average when we copy and paste cell E5 into cells E6 through E16. Look at the equations in the cells after you paste them in to convince yourself that what you expect to happen is actually happening. Adjust the significant figures as needed.

Enter a formula in cell F5 to calculate $(x_i - x_{avg})^2$, and then copy and paste it into the appropriate cells in column F.¹¹ Be sure that each cell displays the correct number of significant figures.

In cell D18 enter a formula to calculate $\sum (x_i - x_{avg})^2$ using the SUM function (see [Guide to Excel](#) section) and in cell E18 enter the units for the numerical value in cell D18.

Enter a formula in cell D19 to calculate S from cell D18 using the SQRT function (see [Guide to Excel](#) section). Treat the square root as you would treat division to determine the number of significant figures. The number of data points, N , is an exact number, and as such does not affect the number of significant figures in the standard deviation. Note that, as with the average, the standard deviation cannot be more precise than the original data! Enter the units of S in cell E19. In cell D20 employ the STDEV function to calculate S directly from the data in cells D5 through D16 (see [Guide to Excel](#) section). Enter the units for cell D20 in cell E20. Note that from now on, if you need to calculate any statistical quantity (i. e., average or standard deviation) you can use the spreadsheet's built-in functions, instead of the elaborate procedure that you used in this lab.

Enter an equation in cell D21 to calculate the uncertainty at the 95% confidence level from the estimated standard deviation in cell D21. The expression that you will need is Eqn. 5 in the *Introduction to Statistics in Chemistry* section. You can use either the TINV function (see *Guide to Excel*) or Table 1 in the *Introduction to Statistics in Chemistry* to find the appropriate value of Student's t . In cell E21, enter the units of Δ .

Print out two copies of this worksheet and attach one copy to a page and the other copy to the duplicate page in your notebook. In your notebook write down the average value for the density of the copper as measured by water displacement, its standard deviation and the confidence limits at the 95% confidence level. Again take care that your calculation has the proper units and number of significant figures. Calculate the percent error in copper's density as measured by this method. Save the Excel file and to record the name you saved it under in your laboratory notebook.

Conclusions

As discussed in the [The Laboratory Notebook](#) section, there are three types of exercises that are encountered in chemistry. This exercise's focus was on measurement, so use the [Outline for Measurement Experiments](#). You have data on copper's density obtained from two different methods, and each will need to be discussed. But, you do not need to have two separate discussions of the two methods; they can be discussed simultaneously.

As part of your conclusions include a discussion of whether the assumptions made in each method were likely to be valid, and how the results would be different if the assumptions were not valid. This discussion can be qualitative (i. e., no calculations required), but it must address all of the assumptions.

When you determined the density by measuring the copper block's volume using the ruler you calculated $\Delta V/V$, the ratio of the error in the volume to the total volume. Include in your discussion what you notice about $\Delta V/V$ for the small blocks as compared to the large blocks. What does this mean and why might this be a problem? How would you modify the experiment to minimize its effect on the density?

Summary of Results

The table that summarizes your results should look like Table 3; substitute your values for those shown. You do not need to list all of the data.

Table 3. Summary of class results for the determination of copper's density by two different methods. All confidence intervals are at the 95% confidence limit.

	From Volume Determined by Ruler	Determined by Water Displacement
Density	$8.94 \pm 0.05 \text{ g/cm}^3$	$9.0 \pm 0.1 \text{ g/cm}^3$
Number of Points	12	11

References and Notes

1. Schlotter, N. E. *J. Chem. Educ.* **2013**, *90*, 51-55. This article is available as a PDF file at <http://pubs.acs.org/doi/pdf/10.1021/ed300334e> for Truman addresses and *J. Chem. Educ.* subscribers.
2. Jolly, W. L. *Encounters in Experimental Chemistry*, 2nd Ed.; Harcourt Brace Jovanovich: New York, 1985, 52-54.
3. In cell F3 you will actually enter "Volume (cm3)", in cell G3 "DV (cm3)" and in cell H3 "DV/V". Use *Format, Cells*, as required, to subscript or superscript text, or to change text to symbol font.
4. The volume of a rectangular box is its *length* x *width* x *height*. Your block's length, width and height are in cells B4, C4 and D4, respectively. So, to calculate the volume in cell F4 enter "= B4*C4*D4" (without the quotes) in cell F4.
5. To make Eqn. 1 useable, we need to put it into a format that Excel can understand. Recall that multiplication is '*', division is '/', powers are '^' and the square root function is SQRT. So, you will enter in cell G4 "=F4*SQRT((0.01/B4)^2+(0.01/C4)^2+(0.01/D4)^2)" (without the quotes).
6. Since density equals mass divided by volume it would seem logical to graph mass as a function of volume so that the resulting straight line's slope would be the density. Although this looks good, it is wrong. This is because when one performs a linear regression on a data set the quantity with the smallest error must be on the x-axis. If this is not done then the results may be

statistically meaningless. In our case we know the mass much more precisely than we know the volume and so the mass must be on the x -axis.

7. Graphs, charts, spectra, or spreadsheet analyses should be affixed to the pages of the notebook with tape or glue (to both the original and duplicate pages of duplicating notebooks). Label the space where this material is to go with a description of the item and the results it contained. This way, if it is removed, there will be a record of it. Make no notes on the inserted material.

8. *CRC Handbook of Chemistry and Physics, 64th Edition*; Weast, R. C., Ed.; CRC Press: Boca Raton, FL, 1984, p. B-11.

9. You will actually enter "Volume H₂O" in cell C3, "Density (g/cm³)" in cell D4, "g/cm³" in cell E17, and "D =" in cell D21 (all without the quotes), and use *Format, Cells* to change the font. The procedure is similar for cells E4, F4 and D18.

10. An absolute reference is denoted using the '\$' character. So, in cell E5 enter (without the quotes) "=D5-\$D\$17" (where \$D\$17 is the cell containing the average) to calculate the difference between each data point and the average.

11. To raise a number to a power in Excel use the '^' character. So, in cell F5 enter (without the quotes) "=E5^2" to square the contents of cell E5.

Preparation and Analysis of Alum

D. L. McCurdy, V. M. Pultz and J. M. McCormick*

Introduction


One of chemistry's goals is to be able to transform any set of substances (the reactants) to another set of substances (the products) through a chemical reaction. As we have discussed in class, there are rules, such as the Law of Conservation of Mass, by which chemical reactions occur, and it took chemists a long time to understand these basic rules. Even though we know a great deal about chemical reactions, chemists are still finding new chemical reactions and new ways of assembling atoms into molecules and molecules into more elaborate structures. In this and the next laboratory exercise you will learn some of the basics of how chemists carry out chemical reactions and how they characterize the chemical substances involved in these reactions.

To fully describe a chemical reaction one needs to know the identities of both the products and the reactants, and the proportions in which the reactants combine and the products form. While it may seem a trivial exercise to identify the reactants, this is not always the case. Needless to say, identification of the reactants in a complex reaction mixture can be very difficult, and so we will only work with chemical reactions where the reactants are known.


The description of a chemical reaction consists of a series of steps: 1) carrying out the reaction, 2) isolating the product(s), 3) purifying the product(s), and 4) characterizing the product(s) and determining its (their) purity. The isolation and purification of the products are based on physical properties such as the ability to form crystals, boiling point, melting point, solubility, etc. Characterization of the products may be either quantitative or qualitative. In a quantitative characterization, the chemical formula and the structure (i. e., how the atoms are connected) are determined. The former is usually accomplished using elemental analysis, mass spectroscopy, X-ray crystallography or some spectroscopic method. Sometimes it is sufficient to show only that certain ions or elements are present in a sample, and in this case a chemist will perform a qualitative test. Qualitative tests often use chemical reactions that result in a visible change (formation of an insoluble solid, a color change, or evolution of a gas) as a way to quickly show whether a particular chemical species is present or not.

Once the chemical reaction's products are fully characterized, and the balanced chemical reaction is known, we can compute a theoretical and a percent yield. We do these final characterizations of the reaction because it is important to know how efficiently the reaction converts reactants to products. Chemists are always trying to strike a balance between the cost of the reactants, the value of the products, the time a reaction requires and the cost of any unwanted by-products that must be handled as hazardous waste. A reaction, even though it gives a valuable product, may be unusable because it has a low yield, takes too much time or generates too much waste.

In this experiment you will prepare and characterize alum (potassium aluminum sulfate dodecahydrate, $\text{KAl}(\text{SO}_4)_2 \cdot 12\text{H}_2\text{O}$). The first step in this synthesis, which you will perform during Week 1, is to react metallic aluminum with a concentrated solution of potassium hydroxide (KOH) to form the potassium salt of the tetrahydroxoaluminate complex ion, $[\text{Al}(\text{OH})_4]^-$. The balanced chemical equation for this oxidation-reduction reaction is



The second step of the procedure is to convert the KAl(OH)_4 to alum by addition of sulfuric acid (H_2SO_4) in an acid-base reaction. Under the experimental conditions, the alum has a limited solubility in water, and so it precipitates from the solution. The balanced chemical reaction that occurs in this step is



The overall balanced chemical reaction for the conversion of aluminum to alum, shown below, can be obtained by adding together the balanced chemical equation for each step (Help Me).



The second and third weeks of this exercise will be devoted to characterizing the alum. Alum is an ionic compound, which means its melting and boiling points are likely to be too high to be measured conveniently. Also, most spectroscopic methods would not yield useful information. Therefore, we will rely on chemical means to show that we did, in fact, form alum in our reaction. This procedure duplicates how chemists characterized chemical reactions until the late 20th century, and in some cases chemical means of characterization are still the only methods available.

In Week 2 you will perform qualitative tests to demonstrate the presence of K^+ , and sulfate ion (SO_4^{2-}) in the alum. You will also perform a quantitative determination to determine the percent water by mass in alum.

The qualitative test for sulfate uses the insolubility of barium sulfate (BaSO_4). When an aqueous solution of a barium salt (usually BaCl_2) is mixed with an aqueous solution containing sulfate, a white precipitate of insoluble BaSO_4 forms according to the net ionic equation:

$\text{Ba}^{2+} (\text{aq}) + \text{SO}_4^{2-} (\text{aq}) \rightarrow \text{BaSO}_4 (\text{s})$. A positive test for SO_4^{2-} , therefore, is the observation of a white precipitate when an aqueous BaCl_2 solution is mixed with the aqueous test solution.

When placed in a flame many elements give the flame a distinctive color; an effect that can be used to determine both which elements, and how much of each one, is present in a sample. Potassium produces a distinctive lavender flame that we can use as a qualitative test for the presence of potassium. Potassium's flame is often difficult to see because sodium, which is often present as an impurity, has an intense yellow flame that masks other colors. The potassium flame can be seen in the presence of sodium by viewing the flame through a dark-blue cobalt-glass filter, which absorbs the yellow light from Na, but allows the light from K to pass. When placed in a flame, aluminum does not change the flame's color, and so a visual flame test cannot be used to show the presence of Al.

Alum is a hydrate, which means that it is a compound that has water molecules trapped within the solid. Hydrates will release some, or all, of their "waters of hydration" upon heating. If the chemical reaction between Al and KOH does produce alum as a product, we would expect that heating the product should result in a decrease in the sample's weight corresponding to the loss of 12 water molecules per formula unit of alum. Thus, if one knows the starting mass of alum, and the amount (mass, and therefore number of moles) of anhydrous alum remaining after all of the water has been driven off, one can calculate the amount of water that was present in the alum (by the Law of Conservation of Mass). A comparison of the experimentally determined waters of

hydration and the number expected from the chemical formula can then be used as evidence for the formation of the desired product. The process by which the waters of hydration are driven off is described by the chemical equation shown below, where the “ Δ ” written below the arrow indicates that heat was applied to the reactant(s).



A quantitative analysis for Al^{3+} will be made in Week 3. Normally, Al^{3+} is colorless, which means that it does not absorb light in the visible portion of the spectrum. So, we will add a dye called aluminon that will react with the Al^{3+} in solution to give a colored complex ion. For a sufficiently dilute solution, the amount of light absorbed by a chromophore (a chemical species that absorbs light) present in the solution is given by Beer's Law, $A = \epsilon \cdot b \cdot C$, where A is the absorbance (how much light the sample absorbs compared to a solution that does not contain the chromophore), ϵ is the molar absorptivity (also known as an extinction coefficient; ϵ depends on the compound and the wavelength of light), b is the path length (how much sample the light must pass through) and C is the concentration of the chromophore.¹ According to Beer's Law the intensity of the color is linearly dependent on the amount of aluminon- Al^{3+} complex present. So if we knew ϵ for the complex ion formed between Al^{3+} and aluminon, we could make a single absorbance measurement and know the $[\text{Al}^{3+}]$ in a solution, and therefore, how much Al was in the original alum sample. Unfortunately, this is neither a precise nor accurate way to make this determination. It is imprecise because it is only a single measurement, and it is inaccurate because 1) we don't know the stoichiometry for the reaction between aluminon and Al^{3+} and 2) the commercially available dye is not pure (ϵ cannot be determined). So, we need a way to increase the method's precision and to overcome the problem with accuracy.

The problems with the colorimetric method are solved by using a calibration curve, which gives the relationship between absorbance and concentration. A calibration curve is constructed by preparing samples with known concentrations of the analyte (in our case, Al^{3+}) and then measuring the absorbance of these samples. If Beer's Law holds, a calibration curve is a straight line, for which we can obtain an equation from a regression analysis. Now if we measure the absorbance of a sample containing an unknown amount of analyte, it becomes a simple matter of substituting this value into the equation for our calibration curve and solving for the concentration. Because more than one measurement was used to construct the calibration curve, we improve our precision. A calibration curve also improves accuracy because only the analyte's concentration changes (everything else, such as the stoichiometry between aluminon and Al^{3+} and the dye's purity, is constant).

When you set up your laboratory notebook for this exercise, treat each week of the exercise as a separate experiment. So, each week will have its own title, statement of purpose, etc. Note that some of your results will actually be determined during a subsequent week. Be sure to carefully read the experimental procedure and be aware that there are a number of potential hazards. Also, there are several places in this exercise where you will be waiting for something to happen. You can substantially shorten your time in lab by working on another section of that week's exercise during these times. Be sure that you have completed all of the calculations for a given week's work before coming to laboratory; if you do not come to laboratory prepared, you will not be able to complete the Week 2 and Week 3 exercises in the allotted time.

Experimental

During the first week of the exercise you will be working with a corrosive substance (KOH) and generating a flammable gas (H_2). It is, thus, imperative that you work in the hood, avoid open flames in the lab while the reaction is in progress and avoid contact with the KOH, the reaction mixture or the mist given off.

Week 1²

Synthesis of Potassium Aluminum Sulfate Dodecahydrate

Obtain a piece of aluminum foil weighing about 0.5 g and weigh it precisely (to the nearest 0.001g). Cut the weighed foil into many small pieces. The smaller the pieces the faster the reaction will go because of the increased surface area exposed to the KOH solution.

Place the small pieces of aluminum in a 100-mL beaker. Add enough hot water to a Styrofoam cup so that, when the 100-mL beaker is placed inside, the beaker is completely surrounded by water, but the water does not spill out of the cup or into the beaker. If water from the hot water bath spills into the beaker the result will be a drastic decrease in the yield of alum.

Place the 100-mL beaker containing the aluminum into the hot water in the Styrofoam cup and transfer everything to the hood. Slowly and carefully add 25 mL of the 1.4 M KOH solution to the aluminum. **CAUTION!** No open flames can be present in lab while the reaction between KOH and Al is taking place. Stir the solution with your glass stirring rod and cover it with a watch glass. Repeat the stirring every few minutes until all of the aluminum dissolves. If the reaction slows down, replace the water in the bath with fresh hot water. If the reaction becomes too vigorous, remove the beaker from the water bath until the reaction subsides. **CAUTION!** Avoid inhaling the gas evolved during this reaction. The gas is not toxic at this concentration, but a fine mist of the corrosive KOH solution is formed by the gas evolution.

When the aluminum has completely dissolved (do not be concerned if the solution appears cloudy or contains black specks), gravity filter the reaction mixture into a 50-mL beaker through fluted filter paper (the instructor will demonstrate). Dispose of the used filter paper in the laboratory garbage. **CAUTION!** The filter paper will be wet with the corrosive KOH solution, so wash your hands after handling the wet filter paper.

Obtain approximately 5 mL of 9 M H_2SO_4 in your 10 mL graduated cylinder. Use a plastic pipet to slowly and carefully add the H_2SO_4 solution to the 50-mL beaker containing your filtered solution. Do not dip the pipet into the filtered solution! Continue adding the H_2SO_4 solution until no more precipitate forms. This should require no more than about 5 mL of the H_2SO_4 solution. Do not add too much H_2SO_4 , or your yield will suffer.³ After the H_2SO_4 addition, carefully stir the new mixture with your stirring rod and record your observations. **CAUTION!** The H_2SO_4 solution is very corrosive and the reaction between the H_2SO_4 and KOH is very exothermic (gives off heat).

Prepare an ice bath. Place the 50-mL beaker with the filtered reaction solution in the ice bath. Do not introduce any of the water from the ice bath into the beaker. Also place a test tube containing 15 mL of 95% ethanol in the ice bath. The ethanol solution will be used to wash residual H_2SO_4 and H_2O from the alum crystals.

After a crop of crystals has formed, set up a vacuum filtration apparatus as shown in Fig. 1. Do not under any circumstances push the rubber tubing more than 1/4" on to the side-arm of the filter flask and do not leave the tubing attached to the flask while the flask is unclamped.



Figure 1. A properly configured vacuum filtration apparatus. Note the clamp that holds the filtration flask to the ring stand which prevents the filter flask from tipping over.

While the vacuum is on, carefully remove some of the supernatant (the solution above a solid) from your crystals using a pipette and wet the filter paper. This will help the paper adhere to the filter and prevents leaks. **CAUTION!** The solution is corrosive. Remove the 50-mL beaker from the ice bath, swirl it gently to suspend the crystals and pour it into the Büchner funnel. Use your glass stir rod to remove any crystals that adhere to the side of the beaker. Once the aqueous solution has been filtered completely (leaving the crystals on the filter paper), place 2 - 3 mL (the plastic pipets hold about 3 mL) of the cold ethanol solution in the 50-mL beaker. Use your glass stirring rod to loosen any remaining solid that clings to the side. Swirl to suspend any crystals remaining in the beaker, and pour the suspension into the filter. Once the ethanol has been filtered away, repeat this washing several times. After the last ethanol wash, allow the vacuum to run for a minute or two to draw air through the crystals to help them dry.

After the ethanol solution has stopped draining from the funnel, inspect the product. If it looks dry, gently prod it with your metal spatula. If it is dry enough to remove from the filter, the solid will not be very sticky and will have the consistency of fine sand. Break the vacuum by removing the vacuum hose from the side-arm of the filter flask, and then turn off the aspirator. Transfer the solid and the filter paper from the funnel to a pre-weighed watch glass with the help of your metal spatula, as your instructor will demonstrate. Carefully scrape any alum that adheres to the side of the Büchner funnel onto the watch glass.

If the alum is dry, the filter paper will separate from the crystals and you can remove the filter paper. Gently scrape any crystals adhering to the filter paper onto the watch glass. If the alum is still too wet, leave the filter paper and remove it next week.

Obtain the mass of the wet alum. You will need to have about 2 g of wet alum (3 g if the mass includes the filter paper) so that you will have enough for the next two weeks. If you don't have enough, collect the second crop of crystals and/or redo the synthesis. Keep the crystals from different crops and syntheses separate. Cover the container holding the crystals with a piece of paper towel, and place it in your drawer to dry.

You may notice that more crystals formed in the filter flask during the washings. This second crop of crystals may also be collected, but if you choose to collect these crystals, they should be kept separate from the main crop. It is always good laboratory practice to keep different crops of crystals separate until the identity and purity of each crop is determined (second crops almost always contain more impurities than the first crop and the time needed to purify them sometimes far outweighs the additional yield). Collect the second crystal crop by vacuum filtration; wash with several small portions of the cold ethanol solution and dry, as described above.

Waste Disposal

All materials used in the synthesis of alum may be disposed in a sink if washed down with copious amounts of water.

Week 2

In this second week of the exercise you will be working with a toxic substance (BaCl_2), a corrosive and noxious substance (NH_3), hot materials and open flames. Care should be exercised in the laboratory at all times to prevent injury.

You are not required to bring your laptop/tablet this week, but you may find it useful to have it in lab to record the class data.

Before doing anything else in the laboratory obtain the mass of each crop of alum to the nearest milligram (three decimal places). Make observations on the crystalline product (color, texture, etc.), and record all of your observations in your laboratory notebook. Share your results with your classmates.

Qualitative Chemical Tests

Perform the following qualitative tests for SO_4^{2-} and potassium on your sample. If you collected a second crop of alum crystals, you should perform the sulfate and potassium qualitative tests on both the first and second crops to determine whether your two crops qualitatively the same.

Sulfate Test

Place a few crystals of your alum in a 6-inch test tube. Add distilled water dropwise while stirring until the alum dissolves. Add one drop of 0.5 M BaCl_2 (barium chloride). Record your results. Does alum contain sulfate?

Potassium Flame Test

The instructor will demonstrate the proper techniques for using the Bunsen burner and heating the needle. In the hood, heat the provided needle in the flame to remove impurities. Once the needle is clean, carefully scoop up a small amount of alum on the end of the hot needle. Place the alum in the flame and heat it until the crystals begin to melt and the solid glows. Note the color of the flame. If your flame is bright yellow (indicating the presence of sodium), try cleaning your needle again, or use the cobalt glass filter. Does this sample contain potassium?

Quantitative Determination of Waters of Hydration

Before beginning this section be sure that your alum sample is powdered and that you have weighed your alum sample and recorded its mass in your notebook!

Set up a ring stand, ring clamp and porcelain triangle, as your instructor will demonstrate.

Clean your crucible by placing a few drops of 1 M NH_3 solution in the empty crucible and scrubbing with a paper towel. **CAUTION!** This ammonia solution has a strong odor and is corrosive. Rinse the crucible with distilled water and place the empty crucible on the porcelain triangle supported by a ring and ring stand.

With the majority of the flame remaining below the bottom of the crucible, heat the crucible until its bottom glows a dull red. After heating for five minutes, remove the flame and let the crucible cool to room temperature on the triangle. **CAUTION!** Do not touch the crucible with your hand. It is extremely hot and will remain hot for several minutes. Remember that a hot crucible looks exactly the same as a cool crucible. When cooled, you can move the crucible to the bench top using the crucible tongs. Do not set a hot crucible on the bench top, because the temperature differential may cause the crucible to shatter. Once you have cleaned the crucible, it is important that you handle it only with the crucible tongs. This prevents burns and will eliminate a systematic error caused by the weight of your fingerprints.

Weigh the cooled crucible (and its cover) to the nearest milligram (three decimal places) and record this mass in your notebook. If the balance does not show three decimal places, notify the instructor. Place about 1.0 g of your alum sample in the crucible. Obtain the mass of the crucible, its cover, and the alum to the nearest milligram and record this in your notebook.

Return the crucible to the porcelain triangle and set the cover slightly ajar so that the water vapor can escape. For the first few minutes gently heat (only the light blue portion of the flame touches the crucible) the crucible by holding the Bunsen burner off to the side. **Take care!** The water can violently leave the alum at this point, if it is heated too strongly. Move the Bunsen burner such that the tip of the inner blue cone is approximately 3 cm below the crucible. Heat until the crucible glows red and continue heating for 10 minutes. If at any time you observe white smoke being given off, or smell an acrid odor, discontinue heating immediately (the sulfate is being decomposed to SO_2).

Remove the heat and completely cover the crucible with the lid. Cool the crucible to room temperature on the triangle (this takes about ten minutes). Weigh the cooled crucible (including its cover and the contents) to the nearest milligram (three decimal places). Using the tongs, move the crucible and contents back to the triangle and repeat the heating step for 10 minutes. When this heating step is over, cover the crucible and allow it to cool on the triangle to room temperature, and then reweigh the crucible, cover, and its contents. Record this second mass in your notebook. If the second mass is within a 50 mg of the mass after the first heating, then you have driven off all of the water. If the masses are not within 50 mg, then repeat the heating procedure until two subsequent masses agree.

Waste Disposal

Once you have made your final weighing, invert the crucible and the anhydrous alum should fall out. If it does not, add some water from a wash bottle and use your metal spatula gently to

dislodge it. The anhydrous alum may be disposed of in the trash or in the sink with plenty of water. Rinse the crucible with distilled water and dry it before returning it to your drawer.

The sulfate test solution contains BaSO_4 and a small amount of BaCl_2 . BaSO_4 is very insoluble in water and so poses little threat the environment (it is even administered to patients who undergo upper- and lower-GI imaging), and the amount of BaCl_2 present in the test solution is very small. Therefore, the sulfate test solution can be disposed of in a laboratory sink followed by a large amount of water.

Week 3

The materials used in the third week are not particularly hazardous, but care should be taken to minimize your exposure by working in the hood whenever possible.

Before coming to the laboratory you must have completed the following: 1) prepare a table, like Table 1, in your notebook's *Results* section in which to write your data for the calibration curve, 2) set up the calculations to calculate the $[\text{Al}^{3+}]$ in Table 1 (the number of significant figures in each volume is given in Table 1 and in Table 2), 3) prepare and save an Excel spreadsheet to graph the calibration curve, and 4) familiarize yourself with the Vernier spectrometer before laboratory (see the [*Operation of the Vernier LabQuest 2*](#) section for details); your instructor will review spectrometer operations before you begin work.

You will need to bring your laptop/tablet to lab this week, because there is a substantial amount of data analysis required that is cumbersome to perform on the LabQuest.

Your instructor will demonstrate how to prepare solutions using volumetric glassware and will review the protocols for using the balances.

Table 1. Example of a table that could be used to present the data for the calibration curve.

Solution Number	Volume of Al^{3+} Stock Solution Used (mL)	Final Volume of Solution (mL)	$[\text{Al}^{3+}]$ (M)	Absorbance at 525 nm
1	0.00	50.00		
2	1.00	50.00		
3	2.00	50.00		
4	3.00	50.00		
5	5.00	50.00		

Colorimetric Determination of Aluminum^{4, 5}

In this exercise you will be using electrically-powered laboratory equipment. Be sure to check the power cables for fraying and breakage before using any piece of equipment. Bring any such problems to the attention of your instructor before attempting to use the equipment.

Preparation of the Aluminum Stock Solution

Precisely weigh out (to the nearest milligram) about 0.1 g $\text{AlCl}_3 \cdot 6\text{H}_2\text{O}$ using an analytical balance. Quantitatively transfer this solid to a 100-mL volumetric flask (assume the flask's volume is 100.0 mL). Add about 10 mL of distilled water and swirl to dissolve the $\text{AlCl}_3 \cdot 6\text{H}_2\text{O}$. If the solid does not dissolve, carefully add small amounts of distilled water, swirling between each addition, until

it does. Add distilled water to bring the level of the solution in the flask to the mark on the neck (this procedure is called “diluting to the mark”). Mix thoroughly by stoppering the flask, and then inverting and shaking the flask. Repeat if necessary.

Pipet 3.00 mL of the aluminum solution that you just made into a 25-mL volumetric flask. Dilute to the mark and mix thoroughly. This is the aluminum stock solution that you will use to construct the calibration curve.

Construction of the Calibration Curve

Number five 50-mL volumetric flasks 1 to 5. Do not add any of the aluminum stock solution to flask 1. To flask 2 add 1.00 mL of the aluminum stock solution; to flask 3 add 2.00 mL; add 3.00 mL to flask 4 and 5.00 mL to flask 5. These measurements must be precise, and so you must use volumetric pipets.

To each flask then add 20 mL of the acetate buffer solution and 5 mL of the aluminon solution (in that order!) and swirl gently to mix. These volume measurements do not need to be highly precise. So, you can use your 50-mL and 10-mL graduated cylinders here. Dilute all the solutions to the mark by adding distilled water and mix thoroughly. Allow the solutions to sit for 20 min while monitoring the solutions’ colors. Note any changes in your notebook.

Follow the spectrometer’s operating instructions to ready the instrument for use, as set forth in the [Operation of the Vernier LabQuest 2](#) section. Fill the cuvette with the buffer solution to use as a blank (**IMPORTANT!** you must use the same cuvette for both the blank and for your samples). Remove any bubbles by gently tapping the cuvette with your finger. **Under absolutely no circumstance are you to tap a cuvette on a table top.** Do **not** handle the cuvette by the clear window (your fingerprints will cause an error in the measurement). Before placing the cuvette into the spectrometer, be sure to thoroughly wipe the clear sides of the cuvette with a Kim-Wipe (do **not** use a paper towel). When placing the cuvette in the spectrometer, be sure that clear sides are aligned with the light beam and that the cuvette is placed in the spectrometer the same way every time. The major sources of error when using these spectrometers come from poor technique, and you can avoid these by following these guidelines every time you make a measurement with the spectrometer.

A representative spectrum of a solution of the Al^{3+} -aluminon complex is shown in Fig. 2. The spectrum should exhibit a broad peak near 525 nm. If the shape of your spectrum looks dramatically different than that in Fig. 2, consult your instructor (it is likely that you forgot to perform one of the dilutions and therefore your sample is too concentrated). Transfer the absorbance spectrum to your laptop/tablet via Wi-Fi or USB device for later printing from Excel. Be sure to record the filename in your laboratory notebook.

Measure the absorbance at 525 nm for solutions 1 through 5.⁶ Graph the absorbance at 525 nm as a function of $[\text{Al}^{3+}]$ in Excel and perform a linear regression of the data by inserting a trendline on the graph (see the [Guide to Excel](#) section). Show your graph to your instructor; once he or she has approved it, you may proceed to the next section.




Figure 2. Representative absorbance spectrum for a dilute solution of the Al^{3+} -aluminon complex in acetate buffer.

Determination of Aluminum in Alum

Precisely (to three decimal places) weigh out about 0.2 g of your alum. Quantitatively transfer, as your instructor will demonstrate, to a small beaker and add distilled water to bring the volume to about 15 mL. Place the beaker on a hot plate in the hood, cover with a watch glass and heat to boiling, stirring occasionally with your glass stirring rod. After stirring rinse the glass rod into the beaker with a small amount of distilled water from your wash bottle. While the mixture is heating, clean and dry (exterior only) your volumetric flasks. Also prepare for a gravity filtration directly into the 100-mL volumetric flask using a long-stemmed glass funnel.

Remove the beaker from the hot plate just as the solution starts to boil. **CAUTION!** The beaker, the watch glass and the hot plate's top are all hot. Your instructor will demonstrate the safe way to remove the beaker from the hot plate. Immediately pour the hot solution into the funnel. As the solution is being filtered rinse the beaker, the bottom of the watch glass and the stirring rod each with several small washes of distilled water into the funnel. When the solution is completely filtered, remove the funnel from the volumetric and rinse the volumetric's neck with several small portions of distilled water. The volumetric should now be cool to the touch, but if it is not, wait until it is. Dilute the solution in the flask to the mark. Transfer 3.00 mL of this solution to the 25-mL volumetric flask and dilute as before.

Pipet 3.00 mL of the alum solution that you just prepared into a 50-mL volumetric. Add 20 mL of acetate buffer and 5 mL of aluminon solution using graduated cylinders and then dilute to the mark with distilled water. Wait 20 min and measure the absorbance at 525 nm, as you did for the other solutions. Your reading should be somewhere between the maximum and minimum values on your calibration curve. If it is not, consult your instructor (it is likely that you forgot to perform

one of the dilutions). When you are satisfied that the result is reasonable, record this absorbance value in your notebook.

Waste Disposal

All solutions that you use this week may be drain disposed with a copious amount of water.

Results and Analysis

Week 1

From the amounts of the reactants that you actually used, calculate the theoretical yield of alum.⁷

Week 2

Calculate the percent yield of alum from the theoretical yield you determined last week and the amount of alum that you actually obtained. Share your percent yield with your classmates.

Determine the percent water by weight in alum and the number of waters of hydration in the alum. Share these numbers with the rest of the class, as instructed. Perform a *Q*-test on the class data, and discard the discordant datum, if there is one. From the class data, calculate the average percent water by weight in alum, the standard deviation of the data and determine the confidence limits at 95% confidence. Consult the [Guide to Excel](#) and [Introduction to Statistics in Chemistry](#) sections, as needed. Based on the known formula for alum, determine the expected value of the percent water by weight in the sample. Calculate a percent error for the class average and for your result. Record all data and calculated results in your notebook. You may do the calculations in Excel, and if you do, you will need to tape copies of your output in your laboratory notebook.

Week 3

Calculate the concentration of the aluminum stock solution (the solution that you had after the second dilution) and the concentration of aluminum in each of the solutions that you prepared from the stock solution. Write these values in your table (Table 1, above) in the *Results* section of your notebook. In your calculations assume that the volumes of the flasks and pipets are as shown in Table 2.

Table 2. Nominal volumes of the volumetric glassware used in this exercise.

Volumetric	Volume (mL)
1-mL Pipet	1.00
2-mL Pipet	2.00
3-mL Pipet	3.00
5-mL Pipet	5.00
25-mL Flask	25.00
50-mL Flask	50.00
100-mL Flask	100.0

From the equation for the best-fit line for the absorbance at 525 nm as a function of $[\text{Al}^{3+}]$ determine the percent aluminum by weight in alum⁹ and share your results with the class. Perform a *Q*-test on the class data and then calculate the average percent Al by weight in alum, the standard deviation of the data and finally find the confidence limits at 95% confidence. Again, consult the [Guide to Excel](#) and [Introduction to Statistics in Chemistry](#) sections, as needed, if you are unsure

how to perform any of these tasks. Determine what the true percent Al by weight is for alum and then calculate a percent error for the class average and for your result. Record your calculations in your notebook, as you did for the Week 2 calculations, and include any spreadsheet output.

Conclusions

The first week of this exercise was a synthesis. Therefore, your *Discussion of Conclusions* section for this week should follow the [Outline for Synthesis Experiments](#). Note that you will not be able to discuss your results for Week 1 until after you have obtained the mass of your product and done the qualitative tests on it. It is advisable to reserve two or three pages in your notebook for the Week 1 *Discussion of Conclusions* and *Error Analysis* when you prepare for Week 2. One important question that you will need to address in your *Discussion of Conclusions* section for Week 1 is why is your percent yield of alum less 100% with specific references to what you did and observed.

Weeks 2 and 3 are both measurement exercises (see [Outline for Measurement Experiments](#)). In the Week 2 *Discussion of Conclusions* and *Error Analysis* you should include a brief discussion of the qualitative test results. In both Week 2 and Week 3 you gather evidence for the identity and purity of your alum. So, you must include a short discussion of whether your quantitative results support your purported synthesis of alum. Although a propagation of error analysis is possible, we won't perform one here. However, you should be able to identify where your major sources of uncertainty are and qualitatively discuss how they affected your results.

Summary of Results

Week 1

Use Table 3 to report your results for Week 1.

Table 3. *Summary Table* for the first week.

Mass of Al Used (g)	Theoretical Yield of Alum (g)	Mass of Alum Obtained (g)	% Yield of Alum

Week 2

Summarize your results for Week 2 using Tables 4 and 5. In the second column of Table 4, write either “positive”, or “negative”, as appropriate. Don't forget to report the confidence interval on the class data in Table 5.

Table 4. *Summary Table* for the qualitative tests.

Test for potassium:	
Test for sulfate:	

Table 5. *Summary Table* for quantitative determination of water in alum.

	Initial Mass of Alum (g)	Mass of Anhydrous Alum (g)	% Water by Mass in Alum	% Error in % Water by Mass	Number of Waters of Hydration
Our Results					
Class Results	-----	-----			

Week 3

Table 6 should be used to summarize the results for the third week's work. Remember to include the confidence interval on the class average % Al by mass in alum.

Table 6. *Summary Table* for Week 3.

	Mass of Alum Used (g)	Slope of Calibration Curve ($\text{AU}\cdot\text{M}^{-1}$)	Intercept of Calibration Curve (AU)	Absorbance of Alum Solution (AU)	% Al by Mass in Alum	% Error in % Al by Mass
Our Results						
Class Results	-----	-----	-----	-----		

References and Notes

1. The absorbance has no units (although sometimes "absorbance units" are used, abbreviated "AU"). The concentration's unit is molar, M, and the pathlength's unit is usually cm. Therefore, the unit of the molar absorptivity is $\text{M}^{-1}\cdot\text{cm}^{-1}$.
2. You will notice that in the synthesis portion of this exercise (Week 1) we measure volumes with graduated cylinders and beakers, but in Week 3 we use volumetric pipets and flasks. This is no accident! In the first week we are less concerned with precision than we are in the third week. This is because of the inherent uncertainty of most synthetic procedures which results from side-reactions and other uncontrollable factors. In the quantitative measurement in Week 3 we also perform a chemical reaction, but one which we know gives a specific answer. And so, we can be much more precise. Since our result, and the conclusions we draw from it, critically depends on how well we made our measurement, we use the more precise volumetric flasks and pipets.
3. In very acidic conditions sulfate exists predominantly as HSO_4^- , the hydrogensulfate ion, which does not combine with Al^{3+} to form an insoluble compound.
4. Smith, W. H.; Sager, E. E. and Siewers, I. *J. Anal. Chem.* **1949**, 21, 1334-1338.
5. Marczenko, Z. *Spectrophotometric Determination of Elements*; Ellis Horwood Ltd.: Chichester, England, 1976, p. 116-117.
6. The spectrometer may not read exactly 525.0 nm. Any wavelength within a few tenths of a nanometer will work just as well. Be sure that you use the same wavelength for all the solutions and that you note the exact wavelength in your notebook.

7. To determine the percent yield of a product in a chemical reaction we need to know the amount of all reactants used, the amount of the product formed and the balanced chemical reaction. From the balanced chemical reaction and the amount of reactants, we determine first the limiting reagent and then theoretical yield of the product. The percent yield is then simply the actual amount of product obtained divided by the theoretical yield times 100.

For the reaction of Al with KOH to form alum the balanced chemical reaction is as follows:



To simplify things we have told you that the Al is the limiting reagent (if you wish, you can check this). Since we do not need to determine the limiting reagent, our first step is to determine the amount of alum that can theoretically be formed from the amount of Al that we have. Let's assume that we used 0.475 g of Al and that we obtained 2.930 g of alum. Note that throughout this discussion extra insignificant figures will be carried along in the calculation to prevent rounding errors (indicated as subscripts at the end of each number).

First, we must determine the moles of Al in 0.475 g of Al. You should get 0.0176₀₄ moles of Al.

$$0.475 \text{ g Al} \left(\frac{1 \text{ mole Al}}{26.982 \text{ g Al}} \right) = 0.0176_{04} \text{ mole Al}$$

Now convert moles of Al to moles of alum using the stoichiometric factor⁸ from the balanced chemical equation. You should have found that the reaction could form 0.0175₀₄ moles of alum.

$$0.0176_{04} \text{ mole Al} \left(\frac{2 \text{ mole alum}}{2 \text{ mole Al}} \right) = 0.0176_{04} \text{ mole alum}$$

Calculate the mass of alum (in grams) from moles of alum. This is the theoretical yield.

CAUTION! The molar mass of alum includes 1 K, 1 Al, 2 S and 8 O and the twelve H₂O! You should get 8.35₁₃ g of alum.

$$0.0176_{04} \text{ mole alum} \left(\frac{474.39 \text{ g alum}}{1 \text{ mole alum}} \right) = 8.35_{13} \text{ g alum}$$

Determine the percent yield. Your result should be 35.1% to the correct number of significant figures, although this would often be reported as 35%.

$$\% \text{ yield} = \left(\frac{2.930 \text{ g alum}}{8.35_{13} \text{ g alum}} \right) \times 100 = 35.1\%$$


The calculation from the mass of Al to the mass of alum could be done in a single step.

$$0.475 \text{ g Al} \left(\frac{1 \text{ mole Al}}{26.982 \text{ g Al}} \right) \left(\frac{2 \text{ mole alum}}{2 \text{ mole Al}} \right) \left(\frac{474.39 \text{ g alum}}{1 \text{ mole alum}} \right) = 8.35_{13} \text{ g alum}$$

8. The stoichiometric factor is the ratio of the stoichiometric coefficients of two substances in a balanced chemical equation. The stoichiometric coefficients are the numbers in front of the substances and represent the number of moles of each that are required for, or are formed in, the

reaction. In this case, 2 moles of Al are required contents and 2 moles of alum are formed. So, the stoichiometric factor that relates moles of alum to moles of Al is 2 moles alum/2 moles Al.

9. We are trying to find the percent Al in our alum sample by mass. To calculate this we need the mass of Al in the sample and the sample's total mass. If our alum sample was pure then this would be trivial because we know the chemical formula of alum. But, our sample is not pure and so we must find the mass of Al in our alum sample. We do this by diluting our original alum sample, adding a dye that will bind with Al^{3+} and measuring the absorbance of the dye solution. Using the absorbance of the diluted alum solution and the equation of the calibration curve, we can find the $[\text{Al}^{3+}]$ in the diluted alum solution and then by working backwards find the mass of Al in the alum sample. This is shown in Scheme 1.



Scheme 1. Steps for determining the mass of Al in the alum sample. Items in boxes are quantities that are either measured or calculated and items to the right of an arrow are things that get us from one box to another.

None of the steps in this scheme are individually difficult, but when taken together they can be intimidating. To help you better understand the method, we will work an example where we start from a calibration curve and guide you step by step through the calculations to the final % Al by mass in the sample. However, it is in your best interest to try to set up the equation and perform the calculation before looking at the answer for a particular step.

Concentration of Al^{3+} in the Dilute Solution from the Absorbance Reading. The first step in determining the % Al by mass in alum is to extract the $[\text{Al}^{3+}]$ in the diluted alum solution from the absorbance reading of the diluted solution and your calibration curve. Your calibration curve should look like the one shown in Fig. 3. Note that your absorbance and $[\text{Al}^{3+}]$ values may differ slightly from those shown below, but your absorbance values should be less than about 1.1 and the $[\text{Al}^{3+}]$ values should be on the order of 10^{-5} M. To find the $[\text{Al}^{3+}]$ from a measured absorbance, one simply substitutes the absorbance of the solution (y) into the equation for the best fit line and solves for the concentration (x).




Figure 3. Typical calibration curve obtained for the colorimetric determination of Al^{3+} .

If we assume that the final diluted alum solution had an absorbance of 0.712, what is the $[\text{Al}^{3+}]$? The answer is that the $[\text{Al}^{3+}] = 2.93_{37} \times 10^{-5} \text{ M}$ (note that we can only keep three significant figures and that we will retain two extra digits, shown as subscripts, to prevent rounding error).

Concentration of Al^{3+} in the Original Alum Solution. We have just found that the $[\text{Al}^{3+}]$ in the final diluted solution is $2.93_{37} \times 10^{-5} \text{ M}$. The next step is to determine the number of moles of Al^{3+} in the final diluted solution using the measured concentration and volume (50.00 mL). You should have calculated that there were $1.46_{68} \times 10^{-6} \text{ moles } \text{Al}^{3+}$ in that solution.

$$50.00 \text{ mL} \left(\frac{1 \text{ L}}{1000 \text{ mL}} \right) \left(\frac{2.93_{37} \times 10^{-5} \text{ mole } \text{Al}^{3+}}{1 \text{ L}} \right) = 1.46_{68} \times 10^{-6} \text{ mole } \text{Al}^{3+}$$

We know that all of this Al^{3+} came from the 3.00 mL of the second solution that we made. From the moles of Al^{3+} above and the volume, you should find that the $[\text{Al}^{3+}]$ in the second solution was $4.88_{95} \times 10^{-4} \text{ M}$.

$$\frac{1.46_{68} \times 10^{-6} \text{ mole } \text{Al}^{3+}}{3.00 \times 10^{-3} \text{ L}} = 4.88_{95} \times 10^{-4} \text{ M } \text{Al}^{3+}$$

Since all of the Al^{3+} in the second solution came from the 3.00 mL that we took from the original alum solution, all we need do is repeat what we just did to find the $[\text{Al}^{3+}]$ in the original alum solution, except using 3.00 mL and 25.00 mL. This gives a $[\text{Al}^{3+}]$ of $4.07_{46} \times 10^{-3} \text{ M}$.

$$25.00 \text{ mL} \left(\frac{1 \text{ L}}{1000 \text{ mL}} \right) \left(\frac{4.88_{95} \times 10^{-4} \text{ mole } \text{Al}^{3+}}{1 \text{ L}} \right) = 1.22_{24} \times 10^{-5} \text{ mole } \text{Al}^{3+}$$

$$\frac{1.22_{24} \times 10^{-5} \text{ mole } Al^{3+}}{3.00 \times 10^{-3} \text{ L}} = 4.07_{46} \times 10^{-3} \text{ M } Al^{3+}$$

There is a shortcut for this process. Remember that we performed a serial dilution of the initial solution, first by taking 3.00 mL of it and diluting it to 25.00 mL and then taking 3.00 mL of this new solution and diluting it to 50.00 mL. Instead of knowing the concentration of the original solution and finding the dilute solution, here we know the concentration of the dilute solution and need to find the original concentration.

$$2.93_{37} \times 10^{-5} \text{ M } Al^{3+} = \left(\frac{3.00 \text{ mL}}{25.00 \text{ mL}} \right) \left(\frac{3.00 \text{ mL}}{50.00 \text{ mL}} \right) C_{initial}$$

Mass of Aluminum in the Alum Sample. With the $[Al^{3+}]$ in the original alum solution we can now calculate the mass of Al in original alum sample. The procedure is simple. Since we know the $[Al^{3+}]$ in the original alum solution and we know its volume (100.00 mL), we can calculate moles of Al^{3+} present. There is a one-to-one relationship between moles of Al^{3+} and moles of Al and we know the molar mass of Al. So, we can calculate the mass of Al present in our original alum sample, which is 0.0109_{94} g .

$$100.00 \times 10^{-3} \text{ L} \left(\frac{4.07_{46} \times 10^{-3} \text{ mole } Al^{3+}}{1 \text{ L}} \right) \left(\frac{1 \text{ mole } Al}{1 \text{ mole } Al^{3+}} \right) \left(\frac{26.982 \text{ g } Al}{1 \text{ mole } Al} \right) = 0.0109_{94} \text{ g } Al$$

% Aluminum by Mass of Aluminum in Alum. Now we have the mass of Al in the alum sample. Let's assume that we used 0.200 g of alum to prepare the original solution. This means that there is 0.0109_{94} g Al in 0.200 g alum, and that alum is 5.50% Al by mass.

$$\% Al = \frac{\text{mass } Al}{\text{mass alum}} \times 100 = \frac{0.0109_{94} \text{ g } Al}{0.200 \text{ g alum}} \times 100 = 5.50\%$$

This page intentionally left blank.

Hydrolysis of Oil of Wintergreen

D. Afzal, A. E. Moody and J. M. McCormick*

Introduction

Oil of wintergreen is an essential oil obtained from wintergreen leaves or sweet birch bark. The primary constituent of oil of wintergreen is methyl salicylate (its structure is shown in Fig. 1), which has a fragrant smell and is responsible for the wintergreen flavor in foods and beverages. Because of the high demand in the food industry for methyl salicylate most of it is made synthetically, which is both cheaper and easier than extracting it from the natural sources.






Figure 1. Structure of methyl salicylate.

Methyl salicylate undergoes hydrolysis in the presence of H^+ or OH^- . A hydrolysis reaction is where something is broken apart by water (*hydro-* = water, *-lysis* = splitting). In the experimental procedure that you will follow, the methyl salicylate will be first reacted with a concentrated NaOH solution (the source of OH^-) to give compound Y, which we will not isolate (compound Y is an example of a synthetic intermediate). A sulfuric acid (H_2SO_4) solution will be added as a H^+ source to convert compound Y into compound X, which we will collect and characterize. The reaction's two steps are shown in Scheme 1.



Scheme 1. Unbalanced chemical equations for the hydrolysis of methyl salicylate to give compound X.

When these two steps are added together, we get the overall chemical equation for the hydrolysis of methyl salicylate shown in Scheme 2.¹ It is important for you realize that none of the reactions in either scheme are complete! We are missing products, and we do not know the identity of either compound Y or compound X.



Scheme 2. Overall unbalanced chemical equation for the hydrolysis of methyl salicylate.

You will hydrolyze oil of wintergreen in Week 1 of this exercise and then in Week 2 you will use chemical, physical and spectroscopic means to identify compound X. First you will demonstrate that compound X is an acid, and then you will titrate a known amount of X with a standardized base to determine its molar mass. You will also compare the melting point of impure and recrystallized X. Pure substances have unique and distinct melting points, while impure substances (mixtures) usually do not have a unique melting point; rather they melt over a range of temperatures. Thus, a melting point determination is a quick and easy way to determine the purity of a substance, if its melting point is neither too high (as is the case with many ionic compounds) nor too low (as for most gases).

In this exercise and in the previous one, we have seen how compounds are characterized by their chemical and physical properties, and until the 1960's these were the primary ways to characterize new compounds. Starting in the 1960's, new methods based on the interaction of matter with electromagnetic radiation revolutionized chemistry. These spectroscopic techniques were faster than the older methods and gave much more information on the substances being analyzed, and so they have almost entirely supplanted the older methods. In this exercise you will be introduced to one of the most widely used and powerful spectroscopic techniques, nuclear magnetic resonance (NMR) spectroscopy.

NMR spectroscopy uses the fact that certain nuclei behave like very small magnets, which in a magnetic field can either line up with the field or against it. The alignment of the nuclei can be flipped when they absorb electromagnetic radiation of the correct frequency (in NMR spectroscopy the frequency is expressed as a chemical shift with units of parts per million, ppm). The frequency of radiation that is needed to perform this flip depends on the nucleus and the environment around the nucleus. And so, information on how atoms in molecules and polyatomic ions are arranged and the bonding interactions between them can be determined. Typical chemical shift values for various arrangements of hydrogen and carbon atoms are shown in Table 1. Note that these are typical values and that an atom's chemical shift might be outside the stated range depending on its environment and what is atoms are bonded to it.

The ^1H nucleus (protium) and the ^{13}C nucleus are the most common nuclei studied with NMR spectroscopy. Protium is the most abundant isotope of hydrogen, and so it is likely that in any particular molecule two or more ^1H nuclei are near enough to each other to interact. This interaction between ^1H nuclei leads to a splitting of the NMR line for each ^1H nucleus (see Fig. 2). The splitting is predictable, and can be used to determine the molecule's structure, but it will not be discussed further in this class.

Because the natural abundance of ^{13}C is low, it is unlikely that two ^{13}C nuclei will be next to each other, so we do not need to worry about the interaction between ^{13}C nuclei. A ^{13}C nucleus can, however, interact with nearby ^1H nuclei and this interaction leads to a splitting of the ^{13}C line. However, the ^{13}C NMR experiment is usually set up such that this interaction between the ^1H and ^{13}C nuclei is suppressed. The result is that a ^{13}C NMR spectrum consists of a set of lines, each corresponding to a unique carbon in the molecule.

A final piece of information that we can extract from a ^1H NMR spectrum is the number of hydrogens present. We do this by find the area under a peak (in the parlance of NMR, we integrate the peak), which is proportional to the number of hydrogen atoms that flip at a particular energy. This integration gives the relative number of hydrogens, but not the absolute number, so it is similar to an empirical formula. A ^{13}C NMR spectrum cannot be integrated, because of the way the experiment is performed, but we can determine the number of different carbon atoms present in the molecule by simply counting the number of peaks in the spectrum.


There is one serious problem with NMR spectroscopy. Consider this: if we are trying to find a ^1H NMR signal from a solute dissolved in a solvent that contains ^1H , how are we ever going to find it in the midst of all the solvent's ^1H ? The answer is "not easily." To avoid this problem NMR spectra are obtained in solvents where ^1H have been replaced by ^2H (deuterium, often given the atomic symbol "D"). Although ^2H has an NMR signal, it is at a different frequency than that of ^1H , and doesn't interfere. Because 100% deuterated solvents are seldom used, there are usually small peaks in the spectrum corresponding to residual un-deuterated solvent that can be ignored in the analysis of the spectrum.

The ^1H NMR spectrum of methyl salicylate is shown in Fig. 2 (note that you should print two copies of both the ^1H and ^{13}C NMR spectra of methyl salicylate and compound X for you laboratory notebook, one of each pair for an original page and the other for a duplicate page). We see six signals; two are single, un-split peaks (called singlets in NMR jargon) and four are multiple peaks (called multiplets). The singlet at 3.927 ppm integrates as 3 hydrogens, and so it must arise from the hydrogens in the $(\text{C}=\text{O})\text{-O-CH}_3$ group of atoms. The peaks in the 6 to 8 ppm range each integrate as 1 hydrogen and their chemical shift is correct for the hydrogens on the ring of carbons (called an aromatic ring). That leaves only the singlet at 10.766 ppm, which integrates as 1 hydrogen, that must be the hydrogen in the -OH group attached to the ring.

Table 1. Typical chemical shifts for selected hydrogen atoms in ^1H NMR spectra (left) and carbon atoms in ^{13}C NMR spectra (right). In each entry the chemical shift is for the atom shown in red.²

Type of Hydrogen (in red)	Typical Chemical Shift (ppm)	Type of Carbon (in red)	Typical Chemical Shift (ppm)
	0.9		13 – 16
	1.3 – 1.5		16 – 25
	1.5 – 2.0		25 – 38
	2.0 – 2.3		20 – 30
	3.8		50 – 90
	6 – 8		125 – 150
	3.8		170 – 175
	9 – 11		177 – 185
	6 – 12		205 – 210

A)



B)




Figure 2. A) Full ^1H NMR spectrum of methyl salicylate, and B) enlarged view of the peaks in the 6 - 9 ppm region.

Methyl salicylate's ^{13}C NMR spectrum is shown in Fig. 3. Counting signals we find 8 singlets and 1 multiplet. The multiplet (near 79 ppm) arises from residual non-deuterated chloroform, while the other signals arise from the eight unique carbons in the molecule. The peak at 52.167 ppm is the carbon in the $-\text{CH}_3$ group, while the five peaks between 110 and 140 ppm are five of the six ring carbons based on their chemical shifts. The peak at 170.533 ppm is in the range expected for a carbon double bonded to one O and singly bonded to both a C and another O. This leaves the

peak at 161.622 ppm as the remaining ring C (note that it is slightly outside the expected range, but we really can't assign it any other way).




Figure 3. ^{13}C NMR spectrum of methyl salicylate in CDCl_3 . The multiplet near 77 ppm arises from the solvent.

In your laboratory notebook, treat both weeks of this exercise as one experiment. That means you will have one *Statement of Purpose*, one *Procedural Outline* (this may be a bit long), one *Results* section, one *Calculations* section and one *Discussion of Results and Conclusions*. Don't worry about preparing a discussion of the first week's results. Wait until you have all of the results at the end of the second week before you start trying to make sense of everything (it will be a lot easier!).

Experimental³⁻⁵

The primary hazards in this week's laboratory exercise are from handling hot flasks (during the hydrolysis reaction itself and the product's recrystallization) and caustic/corrosive solutions. The NaOH and the 3 M H_2SO_4 solution used in the hydrolysis reaction is moderately caustic and corrosive, respectively, and skin contact should be avoided.

Week 1

Place 50 mL of water and a magnetic stir bar (or a few boiling stones) in a 250 mL Erlenmeyer flask. Carefully add 10 g NaOH and swirl to dissolve. **CAUTION!** Sodium hydroxide is very caustic (causes burns) and the solution will become quite hot (dissolution of NaOH is exothermic). Once the solution has cooled add 5 mL (measured to the nearest 0.1 mL) of oil of wintergreen. Record this volume in your notebook and calculate the mass of oil of wintergreen added from its density (1.18 g/mL).

A white solid may appear when the oil of wintergreen is added to the NaOH solution, but will disappear as the solution is heated. Place a watch glass on top of the Erlenmeyer flask. This will provide a place for water vapors to condense and fall back into the flask, maintaining the level of liquid in the reaction (this is called refluxing). Covering the flask also prevents splashing of the caustic reaction mixture out of the flask as it boils. Gently heat the solution until it just boils (hotplate set at a maximum of 3 or 4 on the dial). If the watch glass rattles violently on top of the flask, then you have the hotplate set too high; turn it down. Continue gently heating for 20 minutes.

Remove the flask from the hotplate and allow it to cool to room temperature. Carefully add 3 M H_2SO_4 to the solution with stirring until the solution is acidic (test by touching the solution with your glass stirring rod, and then touching the glass rod to pH paper). As much as 40 mL of the sulfuric acid solution may be required to make the solution acidic. Once the pH paper indicates a slightly acid solution, add an extra 5 mL of the 3 M H_2SO_4 . The white product should precipitate at this point. Cool the mixture in an ice bath to assure complete precipitation. Recover the product by vacuum filtration, using a Büchner funnel. Rinse the solid product well with ice cold water, in the Büchner funnel, with the suction flowing.

Save approximately 10 mg of the product in a small container (such as a vial or crucible) and cover with a Kimwipe to dry. Recrystallize the remainder of the product from water.³⁻⁶ Add about 100 mL water (it is better to use less than 100 mL than more than 100 mL), magnetic stir bar or a boiling stone and your crude product to a 250 mL beaker. Heat the mixture until the product dissolves. If the mixture boils and the solid has not dissolved, carefully add small portions of water, returning the solution to a boil between additions, until it does dissolve. Remove the solution from the heat using a paper towel as a hot pad (as your instructor will demonstrate). Place the beaker on a paper towel and cover it with your watch glass. Allow the solution to cool slowly without disturbance. Needle-like crystals should form fairly quickly (if they don't try scratching the bottom of the beaker with your glass stirring rod). When the solution has returned to room temperature (cool to the touch), vacuum filter the solution to recover the crystals of compound X. Transfer the product from the Büchner funnel to a watch glass, cover with a paper towel and place it in your desk to dry.⁷

Waste Disposal

All chemicals used in the hydrolysis of oil of wintergreen can be flushed down the drain with water.

Week 2

The acid and base solutions used this week are not particularly hazardous, but should be washed off the skin immediately. The melting point apparatus used will get quite hot during the course of the melting point determination. Also, check the power cord for any breaks before using the melting point apparatus.

Before coming to laboratory prepare tables in your notebook to record your titration data. These tables should include places to enter the buret readings, the volume and molarity of the titrant used, and the calculated values of the quantities sought. Remember that you will be doing two sets of titrations (one to standardize the NaOH and one to determine the molar mass of compound X) and that each titration should have its own table. See Tables 2 and 3 for suggestions on how these tables should look.

Table 2. Notebook table to be used for recording the titration data for the standardization of the NaOH solution.

	Trial 1	Trial 2	Trial 3	Trial 4
Final buret reading (mL)				
Initial buret reading (mL)				
Volume NaOH solution used (mL)				
Molarity of NaOH solution (mole/L)				
Average molarity of NaOH solution (mole/L)				

Table 3. Notebook table to be used for recording the titration data for the determination of compound X's molar mass.

	Trial 1	Trial 2	Trial 3	Trial 4
Mass of compound X (g)				
Final buret reading (mL)				
Initial buret reading (mL)				
Volume NaOH solution used (mL)				
Molar mass of compound X (g/mole)				
Average molar mass of compound X (g/mole)				

During this laboratory period you will be doing three things: measure the melting point of the pure and impure compound X, titrate compound X to determine its molar mass and interpret compound X's NMR spectrum. Each of these tasks takes some time, so it is imperative that you plan what you are going to do in lab so that you use your time wisely.

The first thing to do once you are in the laboratory is determine the mass of your dry product and record it in the Results section of your notebook.

Determining the Melting Point of Compound X

Pack a melting point capillary with enough compound X so that the bottom 2-3 mm of the tube are filled. Also prepare a separate sample of impure compound X. Be sure you know which is which. Your instructor will demonstrate the correct method for filling a melting point capillary.

Check that the Mel-Temp is at room temperature before proceeding. Never place a sample into a hot Mel-Temp. Carefully, place the melting point capillary in the Mel-Temp's stage. If it doesn't go in easily, do not force it. Look through the magnifying lens on the front of the Mel-Temp to make sure that the capillary is seated properly. There are three sample slots on the stage, which means you can run up to three samples simultaneously (this requires extra vigilance). Please refer to the [Melting Point Determination Apparatus](#) section for more information on Mel-Temp operation.

Switch on the apparatus and start heating the sample by turning the knob on the front of the Mel-Temp to "4". **CAUTION!** The capillary and the sample stage will become very hot during the course of this experiment. The pure product's melting point range is reported to be 158-160 °C, and at the recommended setting it will take about 10 min for the Mel-Temp to reach this

temperature. During the first few minutes you do not need to monitor the sample closely, but you must not forget about it. When the temperature reaches about 130 °C, decrease the power level so that the temperature change is about 1 °C per minute. This heating rate generally gives the most accurate melting point ranges.

As you pass 150 °C observe the sample more closely and record any changes in its appearance. The temperature at which you observe the first drop of liquid forming is the start of your melting point range. Continue closely monitoring the sample until the last bit of solid melts. This is the end of the melting point range. Record the melting point range in your notebook (e. g., 152.0-155.3 °C). If your sample is pure, your melting point range will be less than 0.1 °C (the melting point is said to be “sharp”). If you miss the start or the end of the melting point range, you will need to start over.

When you are finished, set the knob to “0” and shut off the Mel-Temp. Remove your sample and discard it in the broken glass box. Leave the thermometer in the Mel-Temp.

Testing the Acidity of Compound X


In a small test tube dissolve a small amount (~10 mg) of compound X in a 15% by volume aqueous ethanol solution. Test the acidity of the resulting solution by transferring a small amount of the solution to litmus paper using a stirring rod. An acidic solution will cause blue litmus paper to turn red, while a basic solution will turn red litmus paper blue. Verify that compound X is an acid.

Titration of Compound X

Preparation and Standardization of the Base Solution

Place approximately 30 mL of 1 M NaOH in a clean 600-mL beaker and add 270 mL of distilled water. Swirl, and then cover the beaker with a watch glass. This will become your standard base solution, but it is not critical at this point if we know the exact concentration, because we will standardize it, but calculate the approximate concentration of the NaOH solution.

The exact concentration of the NaOH will be found by titrating a solution that contains a known amount (in this case about 0.350 g) of the weak acid potassium hydrogen phthalate (chemical formula $\text{KC}_8\text{H}_5\text{O}_4$, abbreviated KHP) via the reaction shown in Scheme 4. The titration's end point will be determined using phenolphthalein as the indicator. An indicator is a substance that does not participate in a reaction, but changes in some way in response to the reaction conditions. In this case, phenolphthalein is colorless in acidic solution (where only the monoprotonated form of KHP is present) and bright pink in basic solution (where only the fully deprotonated form of KHP is present). As long as there is any monoprotonated KHP present the solution is acidic and the phenolphthalein is colorless, but as soon as all the monoprotonated KHP has reacted there will be an excess of NaOH present. The solution will then be basic and the phenolphthalein will be pink. The trick, therefore, to getting the most precise and accurate results from a titration using phenolphthalein is to titrate until the solution is the slightest pink.



Scheme 4. The neutralization of potassium hydrogen phthalate (KHP) with NaOH.

Clean and rinse a buret with distilled water. Then rinse with three small portions (about 5 mL) of the base solution. Drain the rinse solution through the buret tip into a waste beaker between rinses. Fill the buret with the base solution and check for air bubbles in the tip. If there are any bubbles present, drain some solution through the tip by opening the buret's stopcock. Refill the buret, if necessary and drain solution from the buret until the solution's meniscus is below the 0.00 mL line on the buret. Record the initial buret reading to the nearest 0.01 mL.

Clean and rinse an Erlenmeyer flask (drying is not necessary). Mass and quantitatively transfer between 0.300 g and 0.400 g of KHP into the Erlenmeyer flask, and then add approximately 50 mL distilled water and 2 drops of phenolphthalein. Be sure to record the mass of the KHP to nearest milligram (by this we mean record the balance reading to the third decimal place).

Run approximately 50% of the volume that you calculated above to be the equivalence point from the buret into the flask as you swirl the flask. This can be done relatively rapidly; however if the solution turns pink and stays pink as you stir, stop the addition. Make a note of the approximate volume at which the color changed and start again, but this time use this volume as your target equivalence point.

As you approach your calculated equivalence point, add the NaOH solution from the buret at a slower rate as you continue to swirl the Erlenmeyer flask. As you near the end point of the titration the pink color will persist for a longer time. Slow the addition of the NaOH further as you near the calculated equivalence point. Add the NaOH solution dropwise with swirling until the first hint of a persistent (lasting more than 30 sec) pink color. If a drop is left hanging on the buret tip, gently touch the drop to the flask's side to dislodge it and rinse it down into the solution with a little distilled water. Record the final volume of the NaOH solution in the buret to the nearest 0.01 mL.

Discard the titrated solution in the sink and then clean and rinse your flask (you don't have to dry it!). Repeat this titration until you have three concentrations that agree within 2%. If you expect that enough NaOH solution remains in the buret to perform another titration, do not refill the buret. It saves time and it may increase the accuracy of your result. Note that this means that the final buret reading of one trial would be the initial buret reading of the next trial.

Titration of Compound X

Weigh accurately between 0.100 and 0.150 g of compound X, and place it in a clean 250-mL Erlenmeyer flask. Add approximately 50 mL of the 15% aqueous ethanol solution with swirling. Note that the solid will not dissolve completely, but that this will not matter. Add 2 drops of phenolphthalein indicator to the solution in the flask.

Refill the buret with the now standardized NaOH solution (the solution in the 600 mL beaker). Record the initial level of the NaOH in the buret. Titrate the solution of compound X.

IMPORTANT! Compound X will dissolve as it reacts with the NaOH. You must be careful during the titration not to add too much NaOH at any one time and to constantly swirl the flask. From time to time rinse the sides of the flask with a little distilled water to dislodge any compound X crystals from the side. Ensure that all solid has dissolved prior to declaring the end point. From the volume of base used, calculate the molar mass of compound X (helpful hint: set these calculations up before lab). Repeat the titration (discarding the titrated solution and rinsing the Erlenmeyer between analyses) until you obtain three molar masses for compound X that agree with each other within 2%.

Waste Disposal

All chemicals used in this week's exercise can be flushed down the drain with water.

Results and Analysis

Melting Point Determination

Pure compounds have very sharp melting points (small melting point ranges). Based on this information, is your compound X pure?

Titration


Assume that compound X is a monoprotic acid (one H^+ per molecule of X will react with NaOH), and calculate the molar mass of compound X from your titration results. Do this for each titration separately, and then calculate an average molar mass and a standard deviation for the molar mass. Calculate the average molar mass for the class as a whole, perform a *Q*-test (if warranted), and calculate the standard deviation of the class data and the confidence interval at the 95% confidence limit. Compare your result to that of the whole class.

NMR Spectral Analysis

The NMR peaks for compound X are given in Table 4 (along with the integration of the ^1H spectrum) and the NMR spectra themselves are given in Figs 4 and 5. You should print out two copies of each spectrum and the table to put in your laboratory notebook (one on the original page and one on the duplicate page). By comparing the NMR spectra of compound X and methyl salicylate assign the spectra of compound X to the best of your ability. Here are some questions to ask as you examine the spectra: 1) how many peaks are there? 2) have you lost any relative to those you started with in methyl salicylate? 3) have you gained any peaks relative to methyl salicylate? 4) what are the similarities between compound X's NMR spectra and those of methyl salicylate? Look carefully at the integration in the ^1H spectrum, it will help you.

Once you have assigned the NMR, you should be able to draw a structure (like that of methyl salicylate shown in Fig. 1) of compound X. You should also have a direct measure of the number of C and H atoms (and by default O atoms, too) in compound X. Calculate the molar mass of compound X based on this information and compare it to your titration results.

A)



B)




Figure 4. A) Full ^1H NMR spectrum of compound X and B) enlarged view of the 6–8 ppm region.




Figure 5. ^{13}C NMR spectrum of compound X.

Table 4. Peaks observed in the ^1H and ^{13}C NMR spectra of compound X.

^1H NMR Chemical Shift (ppm)	Integration	^{13}C NMR Chemical Shift (ppm)
10.41	1	174.899
7.94	1	162.426
7.54	1	137.170
7.28	Solvent	131.164
7.02	1	119.812
6.95	1	118.058
~7	1 (very broad and weak)	111.573
		~77 (solvent)

Since you know the molar mass and chemical formula of compound X, you can balance the chemical equation shown in Scheme 2 (hint: there is only one other product in addition to compound X, use the Law of Conservation of Mass to determine its formula). And with the balanced chemical equation you can calculate a percent yield for the reaction. Share your percent yield with others in the class.

Conclusions

This exercise is primarily a synthesis experiment; use the basic [Outline for Synthesis Experiments](#) as you write your conclusions. When you discuss the properties of the prepared material (section A of the [Outline for Synthesis Experiments](#)), give a brief outline of how you assigned the NMR

spectra and used them to identify compound X. You are also to discuss your titration results (follow section A of the [Outline for Measurement Experiments](#)) and the melting point range. What must be made clear to the reader is how each of these means of characterization allowed you to identify compound X; your logic is important!

References and Notes

1. Note that we simplified the reaction by changing H^+ (aq) + OH^- (aq) on the reactants side to H_2O (l).
2. Ege, S. *Organic Chemistry*; Heath: Lexington, MA; 1984, p. 678 and p. 703.
3. Mohrig, J. R.; Hammond, C. N.; Morrill, T. C. and Neckers, D. C. *Experimental Organic Chemistry*; Freeman: New York; 1998, p. 32-36.
4. Pavia, D. L.; Lampman, G. M. and Kriz, G. S., Jr. *Introduction to Organic Laboratory Techniques: A Contemporary Approach*, 2nd Ed.; Saunders: New York; 1982, 93-94.
5. Byrd, H. and O'Donnell, S. E. *J. Chem. Educ.* **2003**, *80*, 174-176. This article is available as a PDF file at <http://pubs.acs.org/doi/pdf/10.1021/ed080p174> for Truman addresses and *J. Chem. Educ.* subscribers.
6. In a recrystallization we generally want to add only enough solvent to just dissolve the substance when the solution is hot. However, this requires some trial and error to find the best solvent and the optimal volume, which can be a time-consuming and frustrating experience. To save time in this exercise, we will give you that approximately 100 mL of water is typically sufficient for recrystallizing the product.
7. Note that normally in a recrystallization we would gravity filter the hot solution. This would require us to preheat the funnel and the filter paper and perform the filtration very slowly; adding only small amounts of solution to the filter at any one time. Because of time constraints and because we know that it isn't necessary in this case, we will dispense with the filtration step.

Vitamin C Analysis

B. K. Kramer, V. M. Pultz and J. M. McCormick*

Introduction

Vitamins are a group of small molecular compounds that are essential nutrients in many multi-cellular organisms, and humans in particular. The name “vitamin” is a contraction of “vital amine”, and came about because many of the first vitamins to be discovered were members of this class of organic compounds. And although many of the subsequently discovered vitamins were not amines, the name was retained. In this exercise you will be studying vitamin C, also known as ascorbic acid.




Figure 1. Structure of vitamin C (ascorbic acid).

Ascorbic acid ($C_6H_8O_6$) is a water-soluble vitamin, whose structure is shown in Fig. 1. Vitamin C is easily oxidized, and the majority of its functions in vivo rely on this property. It plays a key role in the body's synthesis of collagen and norepinephrine by keeping the enzymes responsible for these processes in their active reduced form.¹ Vitamin C may also play a role in detoxifying by-products of respiration. Occasionally during respiration O_2 is incompletely reduced to superoxide ion (O_2^-) instead of being reduced completely to its -2 oxidation state (as in H_2O). Normally an enzyme called superoxide dismutase converts O_2^- to H_2O_2 and O_2 , but in the presence of Fe^{2+} the hydrogen peroxide may be converted into the highly-reactive hydroxyl radical ($\cdot OH$). The hydroxyl radical can initiate unwanted and deleterious chemistry within a cell when it removes a hydrogen atom ($H\cdot$) from an organic compound to form H_2O and a new, potential more reactive free radical. Ascorbic acid can donate a hydrogen atom to a free radical, and thus stop these reactions from occurring.²

The human body cannot produce ascorbic acid, and so it must be obtained entirely through one's diet. A vitamin C deficiency in humans results in the disease called scurvy, whose symptoms include hemorrhaging (especially in the gums), joint pain and exhaustion.¹ In its final stages scurvy is characterized by a profound exhaustion, diarrhea, and then pulmonary and kidney failure, which result in death.³ A very small daily intake of vitamin C (10-15 mg/day for an adult) is required to avoid deficiency and stave off scurvy.⁴ However, there has been, and continues to be, vigorous debate on what the optimum daily intake of vitamin C is. Some have argued that 200 mg/day is an optimal daily intake for adult humans. Others have suggested 1-2 g/day is best,² this despite studies that show that the blood is saturated with vitamin C at 100 mg/day, and any excess is excreted in the urine. In an attempt to balance the competing claims, and ensure the general

population's good health, the Federal Food and Drug Administration has adopted a the recommended dietary allowance (RDA) of 60 mg/day for adults (aged 15 or older), less for children, and more for pregnant and lactating women.²

Fruits, vegetables, and organ meats (e.g., liver and kidney) are generally the best sources of ascorbic acid; muscle meats and most seeds do not contain significant amounts of ascorbic acid.⁴ The amount of ascorbic acid in plants varies greatly, depending on such factors as the variety, weather, and maturity.⁵ But the most significant determinant of vitamin C content in foods is how the food is stored and prepared. Since vitamin C is easily oxidized, storage and the cooking in air leads to the eventual oxidation of vitamin C by oxygen in the atmosphere. In addition, ascorbic acid's water-solubility means that a significant amount of vitamin C present in a food can be lost by boiling it and then discarding the cooking water.


In this exercise you are to develop a testable hypothesis concerning the amount of vitamin C in a sample, design a procedure to test your hypothesis and then execute your plan in the laboratory. You and your laboratory partner will need to develop this hypothesis and submit a one-page typed proposal, which must be approved by your instructor before you can begin work. Include in your proposal a list of the samples that you will analyze (you will supply the samples) and outline any sample preparation that is not described in the experimental section of this exercise. If you need to analyze many samples, you may collaborate with another group in the laboratory (submit one proposal and one report for everyone involved).

Please remember that certain activities are not easily performed in the laboratory (i. e., cooking food). Therefore, you may need to perform some sample preparations, such as cooking, before coming to the laboratory. Also remember that if you wish to look at the change in vitamin C over time then you must plan ahead to allow for this.

Some possible topics are given below, but you are not limited to this list (creativity will be rewarded). In general, projects that seek to only measure the amount of vitamin C in two or more samples (for example, comparing apple and orange juice) will not be graded as highly because they have a weak hypothesis. The only restriction on your choice of topic is that your work may not involve vertebrate animals or any samples of human origin.

- What factors are important in decomposition of ascorbic acid in a particular system?
- What are the effects of different food preparation methods on ascorbic acid content?
- Are there differences in ascorbic acid content in various parts of a fruit or vegetable?

The amount of vitamin C in a sample will be determined by redox titration using the reaction (shown in Scheme 1) between ascorbic acid and 2, 6-dichloroindophenol (DCIP).^{7, 8} DCIP is used as the titrant because it should 1) only oxidize ascorbic acid and no other substances that might be present, and 2) because it will act as a self-indicator in the titration. To be a self-indicator a substance must be one color in the presence of excess analyte (i. e., ascorbic acid) and another color when the analyte has all reacted.



Scheme 1. Redox reaction between ascorbic acid (vitamin C) and 2, 6-dichloroindophenol (DCIP).

In acidic solutions DCIP is red, but if ascorbic acid is present, it will be reduced to a colorless substance. The solution will remain colorless as more DCIP is added until all of the ascorbic acid has reacted. As soon as the next drop of DCIP solution is added at the solution will be light red, due to the excess DCIP and the end point of the titration has been reached.

There are several limitations to this method that you must consider when you are designing and performing your experiment. First, the presence of particles (as in fruit juice) can interfere with your ability to see the end point, and so cloudy juices will need to be filtered first. Second, when analyzing solid substances you will need to dissolve the vitamin C before you can perform the titration. This may require blending or crushing the material and it may also require the addition of known amounts of water to the material. And finally, red materials cannot be analyzed by this method because it is impossible to see the end point.

There are several other things to consider as you design your experiment. First, you will need to have an estimate of how much titrant (DCIP) will be used in your titration. Look up values for the amount of vitamin C in your sample and from the known concentration of DCIP (approximately 250 mg/L), determine what volume of DCIP will be required to titrate your sample. To achieve the most precise and accurate results, it is preferable that you do not use more than one full buret (40 mL) of your titrant (DCIP) for a single titration. Conversely, a sample that requires only a few milliliters of titrant to reach an endpoint can also have very imprecise results (your instructor will not accept endpoints which occur with less than 1 mL of DCIP). If you find that you need to completely empty and refill the buret during a single run of you sample or your endpoint is reached too soon, you may want to change your titration procedures for that sample so that a more reasonable volume (between 10 and 20 mL) of titrant is required. Finally, remember that all solutions must be prepared using volumetric glassware!

Experimental

Preparation of the 2, 6-Dichloroindophenol Solution

The DCIP solution will be prepared for you by dissolving 0.250 g of 2, 6-dichloroindophenol in about 500 mL of water. Sodium bicarbonate (0.21 g) is then added and dissolved. The resulting solution is finally diluted to 1 L with distilled water. The approximate concentration of this solution is 250. mg DCIP/L. Because the exact concentration can vary from day-to-day, you will need to determine the actual [DCIP] by the standardization procedure given below.

Standardization of the 2, 6-Dichloroindophenol Solution

Before using DCIP to quantitatively measure vitamin C, you must know the concentration of the DCIP solution. You can use the reaction of the DCIP solution with a solution of ascorbic acid with a known concentration to find the concentration of the DCIP solution. This is known as “standardizing” the solution, and it must be done once a day for the DCIP solution.

Each pair will separately standardize the DCIP. Carefully pipet 5.00 mL of a standard ascorbic acid solution into a 250-mL Erlenmeyer flask. Be sure to record the concentration of the standard ascorbic acid (it should be about 0.500 g of vitamin C per 1.000 L). Add 2 mL of the sulfuric acid mixture and about 25 mL of distilled water to the flask.⁹ Swirl the flask to mix the solution. Fill a 50-mL buret with the DCIP solution. Use the DCIP to titrate the ascorbic acid until a permanent (lasting more than 30 sec) light red or pink color appears. Record the volume of DCIP needed to oxidize all of the ascorbic acid. Repeat the procedure on two additional samples of standard ascorbic acid. If you are working with another group, or groups, check that all [DCIP] are in agreement.

Using the balanced equation for the oxidation-reduction reaction between ascorbic acid and DCIP, determine the concentration (in mg/L) of the DCIP solution found with each of your titrations. If the results from your three runs are not within 5% of each other, you should repeat the standardization until you have three runs that are.

Determination of Vitamin C in a Sample

The procedure for determination vitamin C in a sample is identical to that for standardizing the DCIP solution. Simply substitute your sample for the standard ascorbic acid solution. You may need to adjust the volume that you pipet in, based on your calculations of the approximate amount of vitamin C present (pipets with the following volumes will be available for you to use: 1.00, 2.00, 3.00, 4.00, 5.00, 10.00, 25.00 mL). If you use a different volume of sample, you will also need to change the amount of the sulfuric acid solution used such that the proportion between the solutions remains the same. Perform at least three titrations, and repeat as necessary until at least three runs are within 5% of each other.

Waste Disposal

Solid waste (fruit rinds, for example) should be placed in the laboratory trash cans. All other solutions used in this exercise can be flushed down a sink drain with plenty of water.

IMPORTANT! DCIP is an expensive reagent so you must take care not to take more than you need (thus making the excess waste). If you find that you have taken too much material, please see if others in the laboratory can use it before discarding it.

Results and Analysis

Determine the average and the standard deviation for the DCIP concentration and for the amount of vitamin C in your sample. Calculate the uncertainty in the amount of vitamin C in your sample at the 95% confidence limit. Collect your values, and those of other groups that worked with you, in a single table in your laboratory notebook.

Conclusions

The *Discussion and Conclusions* section of your notebook for this exercise should follow the [Outline for Measurement Experiments](#).

In addition to your laboratory notebook, you will also prepare a short written report (one report per group) on your experimental results. Detailed instructions for writing a formal laboratory report can be obtained from [The Laboratory Report](#) section this laboratory manual and an example of a completed formal laboratory report is available in the [Example Formal Report](#) section.

A Laboratory Report template is available on the ChemLab.Truman site (<http://chemlab.truman.edu/formal-lab-report/>). Save the template to your disk. The template is a Microsoft Word file that you should be able to open from any university computer.

When you write your lab report using the template, replace underlined text with your specific information in the appropriate face ((i. e., not underlined) and then delete the underlined text. For example, where it says Title, you would enter your title (not underlined, but in bold) and then delete the word Title.

Any information in parentheses is there for information purposes only; delete it in your final document.

You should delete sections that your instructor has told you that you don't need, or that you have nothing for.

The following are some additional helpful hints for you to consider as you prepare your report.

In your report's *Introduction* section you should explain the problem under investigation and how you tried to solve it. Include enough background that a person who knows something about chemistry, but not this particular problem, can understand why your problem is important/interesting and how you planned on solving it. Your proposal is a good starting point for your report's Introduction.

The *Experimental* section is a relatively detailed description of the procedures that you used. A helpful hint, when you write the *Experimental* section be concise; it should not be a step-by-step recounting of everything that you did.

The report's *Results* section presents what you found. This should succinctly summarize your data and any results calculated from your data; tables and figures may be the best way to present this information.

In the report's *Discussion* section you will try to make sense of your results. This section is almost identical to the *Discussion and Conclusions* section of your notebook, although it is often

longer. For this exercise, use the [Outline for Measurement Experiments](#) as a guide to writing your formal discussion.

The *Conclusions* section of a formal laboratory report may be very short and simply summarizes what you found.

The final section that you should write, although it will be the first thing your reader sees, is the *Abstract*, which is a short summary of the problem to be addressed, how it was addressed and what was learned.

Part of your grade will be based on the originality and/or difficulty of your proposal. Include any modifications, especially ones prompted by your experience in the laboratory. Include all data in a table. Part of your grade will be based on the accuracy and precision of your data, or at least whether you attempted to obtain accurate and precise data.

References and Notes

1. Brody, T. *Nutritional Biochemistry*; Academic Press: San Diego, CA, 1994; pp. x and 450-9.
2. Pauling, L. *Vitamin C, the Common Cold, and the Flu*; W. H. Freeman: San Francisco, 1976, pp. x, 4-5, 21-2, 33, 60-1, 145.
3. Kallner, A. *Annals of the New York Academy of Sciences* **1986**, 498, 418-423.
4. Combs, Jr., G. F. *The Vitamins: Fundamental Aspects in Nutrition and Health*; Academic Press, San Diego, CA, 1992; pp. 4-6 and 24-5 and 223-249.
5. Chaney, M. S.; Ross, M. L. and Witschi, J. C. *Nutrition*, 9th Ed.; Houghton Mifflin: Boston, MA, 1979; pp. 283-295.
7. Thompson, S. *Chemtrek: Small-Scale Experiments for General Chemistry*; Allyn and Bacon: Boston, MA, 1990; pp. 194-212.
8. Boyer, R. F. *Modern Experimental Biochemistry*; Addison-Wesley: New York, 1986; pp. 515-521.
9. The sulfuric acid mixture contains sulfuric acid, acetic acid and EDTA (ethylenediamine tetraacetic acid). The EDTA binds strongly to ions, such as Cu^{2+} and Fe^{3+} , preventing them from catalyzing the oxidation of ascorbic acid by atmospheric oxygen.

Enthalpies of Solution

B. D. Lamp, T. Humphry, V. M. Pultz and J. M. McCormick*

Introduction

Thermochemistry Background

Thermochemistry investigates the relationship between chemical reactions and energy changes involving heat. It was born out of the practical problem of cannon making and today continues to play an important role in almost every facet of chemistry. Practical applications of thermochemistry include the development of alternative fuel sources, such as fuel cells, hybrid gas-electric cars or gasoline supplemented with ethanol. On a fundamental level, thermochemistry is also important because the forces holding molecules or ionic compounds together are related to the heat evolved or absorbed in a chemical reaction. Therefore, chemists are interested in the thermochemistry of every chemical reaction, whether it be the solubility of lead salts in drinking water or the metabolism of glucose.

The amount of heat generated or absorbed in a chemical reaction can be studied using a calorimeter. A simplified schematic of a calorimeter is shown in Fig. 1. The “system” (our chemical reaction) is placed in a well-insulated vessel surrounded by water (surroundings). A thermometer is used to measure the heat transferred to or from the system to the surroundings. Ideally, only the water would be the “surroundings” in the thermodynamic sense, and the vessel would not allow heat to pass. In reality, the vessel does allow heat to pass from the water to the rest of the universe, and we will need to account for that (*vide infra*).




Figure 1. Schematic representation of a calorimeter.

There are two types of calorimeters: constant-pressure and constant-volume.^{1,2} In the constant-volume calorimeter, the chemical reaction under study is allowed to take place in a heavy-walled vessel called a bomb. Because the reaction in the bomb takes place at constant volume, the heat that is generated by the reaction (mostly exothermic reactions are studied in a constant volume calorimeter) is actually the change in the internal energy (ΔU) for the reaction. Although ΔU is a useful quantity, for chemists the change enthalpy (ΔH) is more relevant. However, we can convert ΔU to ΔH using Eqn. 1, if we know the change in the number of moles of gas (Δn) in the reaction and the temperature (T).

$$\Delta H = \Delta U + \Delta n(RT) \quad (1)$$

In a constant-pressure calorimetry experiment, like the one that you will be performing, the energy released or absorbed is measured under constant atmospheric pressure. A constant-pressure calorimeter is simpler to assemble than a constant-volume calorimeter and a wider range of chemical reactions can be studied with it. Also, because the reaction is run at constant pressure, ΔH is equal to the amount of heat a reaction generates or absorbs,³ and one need only measure the temperature change when the reactants are mixed to obtain ΔH for the reaction.

Constant-pressure calorimetry is normally conducted with liquids or solutions that have the same temperature.⁴ When a solid is used, it is usually assumed that the solid's T_{initial} is the same as the ambient temperature. After the T_{initial} measurement is made, the reactants are quickly placed into the constant-pressure calorimeter. If the reactants mix and react instantaneously, and the thermometer responds perfectly to the change in temperature, the change in the temperature (ΔT) would simply be $T_{\text{final}} - T_{\text{initial}}$, as shown in Fig. 2. Note that if the calorimeter is perfect (no heat leaks) the temperature inside the calorimeter will not change, and the graph of temperature as a function of time will be flat, also as shown in Fig. 2.




Figure 2. Graph of temperature as a function of time for an exothermic reaction in a perfect calorimeter.

Unfortunately, no calorimeter is perfect, and instantaneous mixing and reaction are not always achieved (even with efficient mixing). In this case, the graph of temperature as a function of time looks more like Fig. 3. We can still find ΔT , but now we must extrapolate back to when the solutions were mixed (time, t , equals zero). This is most easily done by performing a linear regression on the sloped portion of the graph (where, for exothermic reactions, heat is leaking out of the calorimeter) and obtaining T_{final} from the y-intercept.




Figure 3. Graph of temperature as a function of time for an exothermic reaction in a real calorimeter showing extrapolation back to the ideal T_{final} at the time of mixing ($t = 0$).

Some other experimental problems with real calorimeters that we need to account for are: 1) real calorimeters can absorb heat, and 2) although the species that undergo the chemical change result in a release/absorption of thermal energy, it is the entire solution that changes its temperature. Luckily both of these problems can be accounted for by measuring a calorimeter constant, C , which is essentially a specific heat capacity for the calorimeter and everything in it (with units of $\text{J}\cdot\text{g}^{-1}\cdot^{\circ}\text{C}^{-1}$). As long as we work with dilute aqueous solutions and the nature of the solutions does not change significantly from one experiment to another (e.g., the solutions are all dilute and aqueous), the calorimeter constant may be used for many different experiments in the same calorimeter. The calorimeter constant is most easily determined by performing a reaction with a known enthalpy change (ΔH_{rxn}). For this exercise we will use the neutralization reaction $\text{HCl (aq)} + \text{NaOH (aq)} \rightarrow \text{H}_2\text{O (l)} + \text{NaCl (aq)}$ to determine the calorimeter constant (see the [Using Enthalpies of Formation to Calculate Enthalpies of Reaction](#) section).

To relate ΔH_{rxn} to the calorimeter's temperature change, we need to use the First Law of Thermodynamics (see [The First Law and the Sign Convention used in Thermodynamics](#) section for additional information). The heat that the chemical reaction puts out, or takes up, (q_{rxn}) is simply the moles of the limiting reagent, $n_{limiting\ reagent}$ times ΔH_{rxn} (recall that this is how an enthalpy change was defined), as given by Eqn. 2.

$$q_{rxn} = n_{limiting\ reagent} \Delta H_{rxn} \quad (2)$$

The solution (including the reactants and the products) and the calorimeter itself do not undergo a physical or chemical change, so we need to use the expression for specific heat capacity to relate their change in temperature to the amount of heat (q_{cal}) that they have exchanged (Eqn. 3). In Eqn. 3, m is the mass (mass of the reactants + mass of water + mass of calorimeter), C is the calorimeter

constant (specific heat capacity) and ΔT is the change in the temperature of the solution (and calorimeter).

$$q_{cal} = m \cdot C \cdot \Delta T \quad (3)$$

By the First Law of Thermodynamics, q_{rxn} must be equal in magnitude to q_{cal} , but opposite in sign (if the reaction gives off heat, the calorimeter must take it in). This leads to Eqn. 4, which is the starting point for all of the calculations in this exercise. It is then simply a matter of algebraic manipulation to put it in the form that we need (either solve Eqn. 4 for C , when we are determining the calorimeter constant, or for ΔH_{rxn} when we are trying to find the enthalpy change for a salt dissolving in water, ΔH_{sol}).


$$n_{limiting\ reagent} \Delta H_{rxn} = -m \cdot C \cdot \Delta T \quad (4)$$

Periodic Trends

Trends related to the positions of the elements on the periodic table are a well-established fact. There are vertical trends, horizontal trends, and some properties can trend both ways. Atomic size is an example of a trend. Generally, the lower in a group an element is, the larger it is, so a potassium atom is larger than a sodium atom. Also, the farther right in a row an atom is, the smaller it tends to be, so a carbon atom is smaller than a sodium atom.

In this exercise, calorimetry will be used to investigate whether there is a periodic trend in the enthalpies of formation for the common cations of some metallic elements in aqueous solution. The formation of an aqueous cation from an element in its standard state is a fairly abstract multi-step process, but it relates directly to the oxidation-reduction reactivity of the element and to the solubility of ionic compounds. So, this is an important chemical process!

Each salt investigated in the lab has a metallic cation. Using the known ΔH_f^0 of the solid inorganic salt (including any waters of hydration) and the known ΔH_f^0 of the aqueous anion (from the [Selected Enthalpies of Formation](#) table), the ΔH_f^0 of the aqueous cation can be calculated using the dissolution equation of the salt and the enthalpy of dissolution measured in the experiment. For example, the dissolution equation for aluminum chloride hexahydrate is



The ΔH_f^0 of the Al^{3+} ion can, therefore, be found from the ΔH_f^0 of the other species present in the reaction and the ΔH_{soln} found experimentally. Note that in this example the water changes from being bound in the aluminum chloride crystal lattice to being free liquid water. There is an enthalpy change associated with this process! So, it is important to know whether the solid salt is a hydrate or not, and if so, how many waters are present.

If a periodic trend in the enthalpy of formation of the aqueous cation is present down a column or across a row, it should become apparent from the results. In order to visualize any periodic trend in ΔH_f^0 (if there is one), it is helpful to write the enthalpies of solution out on a periodic table. If no trend is present, that should also be readily apparent.

Experimental

The acid and base solutions used to standardize the calorimeter are moderately corrosive and caustic, respectively. The salts used should be treated as toxic, although many of them present minimal hazards.

You are required to bring your laptop/tablet to lab this week. If you wish to download the LoggerPro software on to your computer and use it to acquire data through the LabQuest, you are welcome to do so, however you assume all responsibility for any damage to your computer incurred in the laboratory. Instructions for installing LoggerPro on your personal computer are given in the [Basic Data Analysis Using Logger Pro](#) section.

Before coming to the laboratory be sure that you have determined ΔH_{rxn} for the reaction of aqueous HCl with aqueous NaOH using the tabulated ΔH_f^0 (see [Selected Enthalpies of Formation](#) section) Helpful hint: you might find this reaction's net ionic equation, $H^+(aq) + OH^-(aq) \rightarrow H_2O(l)$, an easier way to calculate ΔH_{rxn} . It is highly advised that you have set up all of the equations that you will need during the laboratory in your notebook before lab. The main cause of people not finishing this exercise on time is being ill-prepared!



Figure 4. Experimental setup of the constant-pressure calorimeter (shown without the cover in place).

Calorimeter Setup

In this experiment, you will use a computer-based data collection system to record solution temperature as a function of time, and a magnetic stirrer to ensure mixing of reagents (see experimental setup in Fig. 4). Before beginning, read the [Operation of the Vernier LabQuest 2](#) section to review basic operation of the data acquisition system. You will be using two stainless steel temperature probes; one in channel 1 of the LabQuest interface and the other in channel 2. Set the software to collect data every second for 4 minutes (240 sec) and adjust the displayed precision to two decimal places (see the Operation of the Vernier LabQuest 2 section for instructions on how to do this).⁵

Determination of the Calorimeter Constant

Measure and record the mass of a clean, dry Styrofoam cup. Place a dry magnetic stir bar in the cup and record the new mass. This cup will be your calorimeter for the day. Do **not** change cups! Otherwise, you will need to re-determine the calorimeter constant for the new cup.

Measure 50.0 mL of ~2 M NaOH with your graduated cylinder and place it into the cup. Assuming the solution has a density of 1.00 g/mL, determine the mass of the solution. Record the mass of the cup and the solution it contains in your notebook. Do **NOT** place a wet cup or a cup filled with liquid on the balance! Doing so can cause severe damage to the balance. Helpful hint: use the density and the volume to calculate the mass. Be sure that you also record the molarity of the NaOH used in your notebook. Calculate and record the number of moles of NaOH used.

Place 51.0 mL of ~2 M HCl in another clean, dry coffee cup. Again, assuming the density of the HCl solution is 1.00 g/mL, determine the mass of the solution that was used. Do **NOT** place a wet cup or a cup filled with liquid on the balance! Rather, use the solution's density and volume to find mass. Record the molarity of the HCl used. Calculate and record the number of moles of HCl used. Determine whether NaOH or HCl is the limiting reagent.

Assemble the calorimeter apparatus, as shown in Fig. 5, by positioning the cup containing the NaOH solution and stir bar on the magnetic stirrer. Your instructor will assist you in positioning the cover, if needed. Begin gently stirring the solution (a setting of 1 or 2 on the magnetic stirrer is a good starting point).

Rinse the channel 1 temperature probe with distilled water into a beaker and pat dry with a KimWipe. Place the probe in the cup being careful that the stir bar does not strike the probe. Gently clamp the temperature probe in place. Rinse the channel 2 temperature probe and then place it in the HCl solution. **CAUTION!** The temperature probe should not sit in the HCl solution for longer than one minute. If the probe stays in an acidic solution any longer than this, the steel will be irrevocably corroded. The LabQuest will display the temperature of both solutions in real-time in the upper left-hand corner of the window. Monitor the temperatures over the next several minutes. While the temperatures are equilibrating, make sure that the LabQuest software is ready to start recording data. When the temperatures of the NaOH and HCl solutions no longer change, record the temperature of each. Calculate the average of the two temperatures, which will be $T_{initial}$ of the mixture. Remove the probe from the HCl solution and rinse it well with distilled water into a waste beaker.

Move the cover to the side and then rapidly, but carefully, pour as much of the HCl solution as possible into the calorimeter and simultaneously initiate data collection on the LabQuest. Slide the cover back into place. While continuing to stir, record the solution's temperature every second over the next 4 minutes.

It is possible to perform the analysis described below, which was done in LoggerPro, on the LabQuest, however you will want to transfer your data via Wi-Fi or using a USB drive to your own laptop or tablet for the analysis. This will allow you to both save the data for each run and to work up one run's data while the next one is being acquired.

Since we used two temperature probes, the LoggerPro file will contain temperature versus time data for both. However, only the channel 1 probe will change, and it will be the only one that we

will analyze. Be sure that you select channel 1 to analyze. Note that for clarity the signal from channel 2 has been omitted in all figures shown below.





Figure 6. Typical trace of temperature as a function of time for an exothermic as shown in the LoggerPro software. Note that data from only one channel is shown.

The trace shown in Fig. 6 is fairly typical for an exothermic process, where the temperature of the solution rises rapidly before slowly diminishing as the system returns to room temperature. Since the temperature probe cannot respond instantaneously to a rapid change in temperature and the reaction may not take place instantaneously, the first portion of the data may exhibit some curvature before reaching a maximum. However, the data to the right of the curve's maximum should be fairly linear. Use the linear fit icon () to draw the best-fit line extending it back to the time of mixing, i.e., time = 0 min (see Fig. 7). The ideal final temperature of the mixture, T_{final} , is the temperature given by the best-fit line at the time of mixing. If you determine a best-fit linear line based on data to the right of the curve's maximum, the intercept is T_{final} . Calculate the ideal temperature change $\Delta T = T_{final} - T_{initial}$.

Determine the total mass of the calorimeter, m (includes the mass of the cup and everything in it), by adding the mass of the dry cup and stir bar, the mass of HCl and the mass of NaOH. Using the total mass, ΔH_{rxn} , moles of the limiting reagent and ΔT , calculate the specific heat capacity of the calorimeter, C .

Repeat the above procedure twice more. Calculate an average specific heat capacity of the calorimeter and its associated 95% confidence interval.

Determination of a Heat of Solution

In this portion of the exercise, you will use the calorimeter from the previous portion to determine the heat of solution (ΔH_{soln}) for an inorganic salt. Your specific salt will be assigned by your

instructor in the laboratory; all measurements are to be conducted using your assigned salt. Since you will not be using the second temperature probe, you can disconnect it.




Figure 7. Same data as shown in Fig. 6, but now with the results of the linear regression shown. As in Fig. 6, data from only one channel is shown.

Clean and dry the coffee cup that you used for the calorimeter in the first part. Place 50.0 mL of distilled H₂O in the cup. Assuming the water has a density of 1.00 g/mL, determine the mass of distilled water used. **Do NOT place a wet cup or a cup filled with liquid on the balance!**

Assemble the calorimeter apparatus, insert the magnetic stir bar and begin gentle stirring. Rinse the channel 1 temperature probe with distilled water and pat dry. Place the probe in the water, as you did before, and note the temperature of the water over the next several minutes. When the temperature no longer changes, record it as $T_{initial}$.

Grind your assigned salt to a fine powder with a clean, dry mortar and pestle. Place approximately 2.0 g of the powdered salt into a clean, dry weigh boat and record the mass. The salt should be at room temperature, which we will assume is the same as the temperature of the water.

Begin stirring the water in the calorimeter. This should be fairly vigorous, but not so vigorous that water splashes out of the calorimeter or there is excessive cavitation in the water. Slide the cover out of the way, initiate data collection and then rapidly, but carefully, add the salt to the stirring water in the calorimeter. Slide the cover back over the cup's mouth. While continuing to stir, record the solution's temperature every few seconds over the next 15 minutes. The time required to obtain the maximum/minimum temperature may be as short as 5 minutes and as long as 40 minutes (if the sample was not ground finely enough); adjust your acquisition parameters as required. The appearance of your data will depend on how exothermic or endothermic the dissolution of your salt is. As with the HCl/NaOH data, draw the best-fit line through the data

points which are approaching room temperature. The ideal final temperature of the mixture, T_{final} , is the temperature where the best-fit line crosses the time of mixing. If your data looks really strange, you might approximate T_{final} by the lowest temperature, for an endothermic reaction, or the highest temperature, for an exothermic reaction, that is achieved. Calculate ΔT .

Using the total mass of the solution (mass of cup and stir bar from first part, mass of water added and mass of salt) m , the number of moles of solute, ΔT , and the previously established specific heat capacity of the calorimeter, calculate the heat of solution, ΔH_{soln} for your salt.

Store the latest run and repeat the analysis of your salt two additional times (don't forget to save your data!). Calculate the average ΔH_{soln} , with its associated 95% confidence interval for your salt. Before you leave the laboratory, report your results to the rest of the class. Copy one run each for the HCl/NaOH and ΔH_{soln} portions of the experiment into Excel and include a printout of a plot of each dataset in your notebook.

Waste Disposal

Everything used in this exercise **EXCEPT** lithium salts, barium salts and their solutions can be flushed down a drain. All lithium salts, barium salts and their solutions **MUST** be collected for proper disposal in the *Secondary Hazardous Waste Accumulation Center* located in the laboratory.

Results and Analysis

Determination of the Calorimeter Constant

Determine the average C for your calorimeter from your three runs. Determine the estimated standard deviation and the 95% confidence interval for C . You will use the average C in your calculation of ΔH_{soln} , but we will not do a propagation of error analysis.

Determination of a Heat of Solution

From your three runs determine an average ΔH_{soln} for your salt. Also calculate the estimated standard deviation and the 95% confidence interval for ΔH_{soln} . Report your ΔH_{soln} and its 95% confidence interval to the class. From your ΔH_{soln} and the tabulated ΔH_f^0 , determine ΔH_f^0 for the cation in the salts. Be careful how you write the reaction that describes the salt dissolving (hydrates are different than anhydrous salts!). There is no need to propagate the uncertainty here (so there will be no confidence interval on ΔH_f^0 for the cations).

Conclusions

For your conclusions use the [Outline for Measurement Experiments](#). Examine class data as a whole. Do you see any trends (for example, how does ΔH_f^0 for the cations change as a function of an element's place on the periodic table)? Predict ΔH_f^0 for the cations that were not studied (e. g., Rb^+ and Sr^{2+}). Evaluate whether your results are likely to be accurate, based on the uncertainty in the measurements that you and your classmates made. Model your *Summary Table* after Table 1, below.

Table 1. Example of the *Summary Table* for this exercise. Fill in your values, and remember to include the 95% confidence interval for each ΔH_{soln} .

Salt	ΔH_{soln} (kJ/mole)	ΔH_f^0 (kJ/mole)for Cation

References and Notes

1. Zumdahl, S. S. *Chemical Principles*, 4th Ed.; Houghton-Mifflin: New York, 2002; chapter 9.
2. Atkins, P. *Physical Chemistry*, 6th Ed.; W. H. Freeman: New York, 1998; chapters 2 and 3.
3. Because constant-pressure calorimeters are often open to the atmosphere (the source of the constant pressure), there might be expansion work. If a gas is evolved by the reaction, the gas will do work as it pushes the atmosphere out of the way. This does not affect ΔH for the reaction, so we usually ignore it. It does, however, affect ΔU .
4. If the reactants are not at the same temperature, then there will be an additional heating/cooling process occurring that is not part of the chemical reaction for which we are attempting to determine ΔH . If we do not correct our data for this process, then our value of ΔH will be in error. Since it is usually not possible to account for this additional heating/cooling, we need to make sure that the temperatures are the same.

Selected Enthalpies of Formation

Aqueous Ions

Chemical Species	ΔH_f° (kJ/mole)
H ⁺ (aq)	0.000
OH ⁻ (aq)	-229.99
C ₂ H ₃ O ₂ ⁻ (aq)	-486.01
Br ⁻ (aq)	-121.55
Cl ⁻ (aq)	-167.16
I ⁻ (aq)	-55.19
NO ₃ ⁻ (aq)	-205.0
CO ₃ ²⁻ (aq)	-677.14
HCO ₃ ⁻ (aq)	-691.99
SO ₄ ²⁻ (aq)	-909.27

Bromide Salts

Chemical Species	ΔH_f° (kJ/mole)
KBr (s)	-393.80
NaBr (s)	-361.06

Iodide Salts

Chemical Species	ΔH_f° (kJ/mole)
KI (s)	-327.90
NaI (s)	-287.78

Nitrate Salts

Chemical Species	ΔH_f° (kJ/mole)
Al(NO ₃) ₃ ·9H ₂ O (s)	-3757.06
Ba(NO ₃) ₂ (s)	-992.07
KNO ₃ (s)	-488.96
LiNO ₃ (s)	-483.13
Mg(NO ₃) ₂ ·6H ₂ O (s)	-790.65
NaNO ₃ (s)	-461.278
NH ₄ NO ₃ (s)	-365.56

Molecular Compounds

Chemical Species	ΔH_f° (kJ/mole)
H ₂ O (l)	-285.830
HC ₂ H ₃ O ₂ (aq)	-484.5

Hydroxide Salts

Chemical Species	ΔH_f° (kJ/mole)
Ba(OH) ₂	-946.30
KOH	-424.76
LiOH	-484.93
NaOH (s)	-425.61
NaOH (aq)	-470.1

Chloride Salts

Chemical Species	ΔH_f° (kJ/mole)
AlCl ₃ ·6H ₂ O (s)	-2691.6
BaCl ₂ (s)	-858.6
BaCl ₂ ·2H ₂ O (s)	-1460.13
CaCl ₂ (s)	-795.8
CaCl ₂ ·2H ₂ O (s)	-1402.9
CuCl (s)	-137.2
CsCl	-442.83
KCl (s)	-436.75
LiCl (s)	-408.61
MgCl ₂ (s)	-641.32
NaCl (s)	-411.15
NH ₄ Cl (s)	-314.43
SrCl ₂ ·6H ₂ O (s)	-2623.8

Sulfate Salts

Chemical Species	ΔH_f° (kJ/mole)
K ₂ SO ₄ (s)	-1437.71
MgSO ₄ (s)	-1284.9
MgSO ₄ ·7H ₂ O (s)	-3388.71
Na ₂ SO ₄ (s)	-1375.653
Na ₂ SO ₄ ·10H ₂ O (s)	-4327.26

Carbonate Salts

Chemical Species	ΔH_f° (kJ/mole)
Li ₂ CO ₃ (s)	-1215.9
Na ₂ CO ₃	-1130.77
Na ₂ CO ₃ ·H ₂ O (s)	-1431.26

References

Afeefy, H. Y.; Liebman; J. F. and Stein, S. E. "Neutral Thermochemical Data" in *NIST Chemistry WebBook, NIST Standard Reference Database Number 69*, Eds. Linstrom, P. J. and Mallard, W. G., National Institute of Standards and Technology, Gaithersburg MD, <http://webbook.nist.gov>, last accessed October 24, 2013.

Atkins, P. and de Paula, J. *Physical Chemistry*, 9th Ed. W. H. Freeman: New York, 2010.

Kotz, J. C. and Treichel, Jr., P. *Chemistry and Chemical Reactivity*, 4th Ed. Saunders College Publishing: New York, 1999.

Wagman, D. D.; Evans, W. H.; Parker, V. B.; Schemm, R. H.; Halow, I.; Bailey, S. M.; Churney, K. L. and Nuttall, R. L. *J. Phys. Chem. Ref. Data* **1982**, 11 (Suppl. 2).

Using Enthalpies of Formation to Calculate Enthalpies of Reaction

Enthalpies for chemical species cannot be measured directly. However, we can use standard enthalpies of formation, ΔH_f^0 , as functional equivalents of a substance's enthalpy. This is because enthalpy is a state function and we define ΔH_f^0 of an element to be zero when it is in its standard state.

To calculate ΔH for a reaction (ΔH_{rxn}^0), we simply need to use the equation

$$\Delta H_{rxn}^0 = \sum \Delta H_f^0(\text{products}) - \sum \Delta H_f^0(\text{reactants})$$

So, for example ΔH_{rxn}^0 for the reaction $\text{HCl (aq)} + \text{NaOH (aq)} \rightarrow \text{H}_2\text{O (l)} + \text{NaCl (aq)}$ would be

$$\begin{aligned} \Delta H_{rxn}^0 = & (1 \text{ mole}) \Delta H_f^0(\text{H}_2\text{O, l}) + (1 \text{ mole}) \Delta H_f^0(\text{NaCl, aq}) \\ & - (1 \text{ mole}) \Delta H_f^0(\text{HCl, aq}) - (1 \text{ mole}) \Delta H_f^0(\text{NaOH, aq}) \end{aligned}$$

We could also use the net ionic equation to calculate ΔH_{rxn}^0 . The ΔH_f^0 for NaCl (aq), HCl (aq) and NaOH (aq) include ΔH_f^0 for the ions into which they dissociated in aqueous solution. But since the Cl^- and Na^+ are spectator ions, their ΔH_f^0 do not need to be included in the overall calculation of ΔH_{rxn}^0 . So, ΔH_{rxn}^0 for the reaction $\text{HCl (aq)} + \text{NaOH (aq)} \rightarrow \text{H}_2\text{O (l)} + \text{NaCl (aq)}$ is the same as for the reaction $\text{H}^+ (\text{aq}) + \text{OH}^- (\text{aq}) \rightarrow \text{H}_2\text{O (l)}$, which is

$$\Delta H_{rxn}^0 = (1 \text{ mole}) \Delta H_f^0(\text{H}_2\text{O, l}) - (1 \text{ mole}) \Delta H_f^0(\text{H}^+, \text{aq}) - (1 \text{ mole}) \Delta H_f^0(\text{OH}^-, \text{aq})$$

The First Law and the Sign Convention used in Thermodynamics

The version of the First Law of Thermodynamics with which you are probably most familiar is “energy can neither be created nor destroyed”. One mathematical statement of the First Law is given by Eqn. 1, where ΔU represents the change in a system’s internal energy (the combined kinetic and potential energies of the particles within the system), q is the heat that the system exchanges with its surroundings and w is the work done on, or by, the system.

$$\Delta U = q + w \quad (1)$$

The system’s internal energy is a *state function*. This means that U is completely defined by the state variables of pressure (p), volume (V), temperature (T) and the amount of material in the system (n). It also means that upon going from one state (one set of state variables) to another state (a different set of state variables) the change in U (ΔU) depends only on the state variables in each state and not on how the system was taken from one state to another. For example, imagine a system that begins with state variables p_1, V_1, T_1, n_1 . Let us then say that the system undergoes a number of processes that gives it a different set of state variables, but it is ultimately brought back to the original conditions (i. e., p_1, V_1, T_1 and n_1). No matter what was done to the system from beginning to end, the final internal energy of the system will be equivalent to its initial value ($\Delta U = 0$).

In general it is not possible to determine a system’s total energy. However, it is possible to determine *differences* in a system’s internal energy as it undergoes a process going from one state to another. The mathematical version of the First Law makes clear that there will be no change in a system’s internal energy unless heat either leaves/enters the system and/or work is done by/on the system. By measuring how much heat leaves/enters a system and/or how much work is done by/on the system, one will be able to determine how much the internal energy of the system is changed.

The sign convention for designating whether thermal energy is leaving or entering a system is simple. Assuming no work is done, if heat leaves a system, the total energy content of the system decreases. Thus ΔU is negative and q is negative, and the process is said to be *exothermic*. In an *endothermic* process, heat enters a system causing the system’s total energy content to increase and so ΔU and q are positive. The sign conventions for work are similar, but require some explanation. The problem is that we need to clearly define what “work done by/on a system” means.

A constant force (F) applied through a displacement (Δd) gives work as given by Eqn. 2.

$$w = -F \times \Delta d \quad (2)$$

There are a number of different types of work, including electrical work. However, the easiest to understand, and the relevant to our everyday experience is expansion work (also known as pressure-volume work), which is work associated with a system expanding or contracting against its surroundings. Under even moderate changes in pressure, the volumes of solids and liquids do not change appreciably and they are said to be virtually incompressible. And so, we will limit our discussion of expansion work to gases.

Consider a system consisting of a gas confined in a cylindrical container of cross-sectional area A , the volume of which is determined by the position of a weightless piston, as shown in Fig. 1. The particles of gas within the cylinder will exert a pressure, p , on the piston,




Figure 1. Ideal thermodynamic system consisting of an ideal gas trapped in a cylindrical container separated from the constant pressure of the atmosphere (p_{ext}).

The particles of gas within the cylinder will exert a pressure, p_{gas} , on the piston, which will be positioned such that the pressure of the gas is equal to a constant external pressure, p_{ext} . Pressure is defined to be a force applied to an area, so the force associated with the work done when the gas expands/contracts is given by Eqn. 3.

$$F = -p_{\text{ext}} \times A \quad (3)$$

When the system of gas undergoes a change that would alter its pressure, the piston will move so that, once again, the pressures inside and outside the system are equivalent. Alternatively, the position of the piston may be forcibly changed, which would cause the gas to respond accordingly. A change in the position of the piston, Δd , causes a change in the volume of the gaseous system, and because volume is the product of height and cross-sectional area, the displacement associated with the work done when the gas expands/contracts can be written as Eqn. 4.

$$\Delta d = \frac{V}{A} \quad (4)$$

The expansion work done when the system undergoes a change can then be related to the change in volume, as shown in Eqn. 5.

$$w = -(p \times A) \left(\frac{\Delta V}{A} \right) = -p\Delta V \quad (5)$$

If work is done by the system, it pushes against the surroundings and expands. This means that ΔV is positive and that w and ΔU are negative. In this case the system's internal energy falls as some of its energy is transferred to the surroundings as work. If work is done on the system, the surroundings push against the system and it contracts. Now ΔV is negative, w is positive, and ΔU is positive, which means that energy is entering the system in the form of work done on it by its surroundings.

The difference in a system's internal energy upon undergoing a chemical change (e. g., combustion of a hydrocarbon) can be directly measured by using a constant-volume calorimeter, or "bomb" calorimeter. The chemical reaction takes place in a heavy-walled container so that ΔV must equal zero. If no other work is done, then expansion work ($w = -p \cdot \Delta V$) must be zero, and by the First Law of Thermodynamics, the heat that the system exchanges with its surroundings is equal to ΔU directly (see Eqn. 6), where the "v" subscript on q reminds us that the heat was measured at constant volume.

$$\Delta U = q_v \quad (6)$$

For many reactions it is not convenient to use a bomb calorimeter to determine the changes in internal energy. Many chemical reactions take place at constant pressure. However, since pressure-volume work may be difficult to measure, it is convenient to measure enthalpies of reaction instead of changes in internal energy. Enthalpy (or "heat content"), represented by H , is defined to by Eqn. 7. Like internal energy, enthalpy is a state function.

$$H = U + pV \quad (7)$$

Also like internal energy, it may not be possible to determine the enthalpy of a system at a particular state, but fortunately we are interested in the changes in enthalpy, ΔH , going from one state to another. If the pressure is constant ($\Delta p = 0$), the relationship between ΔH and ΔU can be written as Eqn. 8.

$$\Delta H = \Delta U + p\Delta V + V\Delta p = \Delta U + p\Delta V + 0 = \Delta U + p\Delta V \quad (8)$$

When only expansion work is considered, the change in internal energy is given by Eqn. 9, which when substituted into Eqn. 8 gives Eqn. 10, where the subscript "p" indicates that this heat is measured at constant pressure.

$$\Delta U = q - p\Delta V \quad (9)$$

$$\Delta H = (q - p\Delta V) + p\Delta V = q_p \quad (10)$$


Kinetics of Crystal Violet Bleaching

V. C. Dew and J. M. McCormick*

Introduction

Chemists are always interested in whether a chemical reaction can occur and exactly how it occurs. The first question is answered through thermodynamics, as you saw in a previous laboratory exercise, while the second is the domain of kinetics. In a kinetics experiment, a chemist attempts to understand the step-by-step transformation of reactants to products. Taken together these elementary steps give us the mechanism by which the reaction proceeds. Note that a reaction's kinetics are very much tied to the pathway the reactants take to the products (i. e., the mechanism), which is very different from the reaction's thermodynamic properties (i. e., ΔH , ΔS and ΔG) that do not depend on the path. While the thermodynamics and kinetics of a reaction may at times seem complementary, and at other times seem contradictory, it is always important to have a detailed understanding of both.

In this experiment you will determine the rate law for a chemical reaction. The rate law is a mathematical expression that relates the amount of time it takes a reaction to happen to the concentrations of the starting materials. The disappearance of reactant over time depends on the rate constant, and the concentration of each reactant raised to some power. This power is known as the order with respect to that reactant. The sum of the individual orders is the overall order of the reaction. The order of reaction with respect to each reactant, as well as the rate law itself, cannot be determined from the balanced chemical equation; each must be found experimentally.¹ The rate law is the basic equation of kinetics and it will be the standard against which we judge possible mechanisms.



Scheme 1. Reaction of crystal violet with OH^- .

In this experiment you will determine the rate law for the reaction of the dye crystal violet (CV) with OH^- in aqueous solution according to the balanced net ionic equation given in Scheme 1. We will define the rate of reaction as the disappearance of the colored CV over time, which can be expressed in differential form as $d[\text{CV}]/dt$. So, the rate law for this reaction can be written as shown in Eqn. 1 in terms of the concentrations of CV and OH^- and the rate constant for the

$$\text{rate} = -\frac{d[\text{CV}]}{dt} = k[\text{CV}]^x[\text{OH}^-]^y \quad (1)$$

reaction, k . In writing this equation we assume that both CV and OH^- are involved in the reaction (that is x and y are both not zero and are likely integers), but only the experiment will tell us whether these assumptions are valid.

The point of any kinetics experiment is to determine the order with respect to each reactant (i. e., find x and y) and to find the value of k . This is a problem if we have more than one reactant (as we do here), in which case the isolation method is often used. The isolation method entails making the concentration of all but one of the reactants very high (so their concentrations do not change appreciably over the course of the reaction). The order with respect to the isolated reactant is then determined. The process is then repeated, isolating each of the other reactants in turn, until all of the orders have been determined.

In this experiment we will make the $[\text{OH}^-]$ very large and, therefore, essentially constant. We can then simplify Eqn. 1 to Eqn. 2, where we have defined a new rate constant, k_{obs} , which is the observed rate constant at some specific $[\text{OH}^-]$. The relationship between k_{obs} and the intrinsic rate constant, k , for this reaction is given by Eqn. 3.

$$\text{rate} = k_{obs}[\text{CV}]^x \quad (2)$$

$$k_{obs} = k[\text{OH}^-]^y \quad (3)$$

Under conditions of high, constant $[\text{OH}^-]$, the order with respect to CV can be determined by graphically applying the integrated rate laws.² Since the absorbance of a CV solution is directly proportional to the concentration of CV, according to Beers' Law,³ the actual $[\text{CV}]$ can be replaced by A_{max} , the solution's maximum absorbance (somewhere around 600 nm). A graph of A_{max} as a function of time will give a straight line if the reaction is zero-order in CV ($x = 0$). If the reaction is first-order in CV ($x = 1$), then a graph of $\ln(A_{max})$ as a function of time is linear. And finally, if a graph of $1/A_{max}$ as a function of time is linear, it indicates that the reaction is second-order with respect to CV ($x = 2$). In each case, if a particular relationship is linear, then the slope of that graph can be used to determine k_{obs} . Note that only one of these three graphs will be linear!

In some instances it is not possible to isolate one of the reactants, because the concentration of that reactant must remain high for the system to behave predictably, as is the case here. However, the order of the reaction with respect to OH^- , and k , can still be found. First we need to change Eqn. 3 into an easily-graphed form by taking the logarithm of both sides to give Eqn. 4 (note that the natural logarithm would also work). To determine the order with respect to OH^- and k , we first perform the kinetics experiment at different, albeit still high, OH^- concentrations and then graph $\log(k_{obs})$ for these reactions as a function of $\log[\text{OH}^-]$. The slope of this graph is y , the order with respect to OH^- , and the intercept is $\log(k)$.

$$\log(k_{obs}) = y\log[\text{OH}^-] + \log(k) \quad (4)$$

In this experiment each you will work in a group with the other students at your laboratory table. Each group will measure the absorbance of CV as a function of time at a different assigned hydroxide concentration using a spectrometer. Each group will determine the order of reaction for

CV at their $[\text{OH}^-]$. The class data (k_{obs} at different $[\text{OH}^-]$) will then be pooled and used to determine the order with respect to OH^- and the intrinsic rate constant.

Experimental⁴

The majority of the substances used in this exercise have minimal toxicity. Please note, however, that the OH^- solutions are caustic and the crystal violet solutions will stain skin and clothing.

You must bring your laptop/tablet to the laboratory this week so that you can perform the required data analysis during the lab period.

Upon arrival in class you will be assigned hydroxide concentrations (either 0.050 M, 0.040 M, 0.030 M, 0.025 M, 0.010 M or 0.0050 M). Before you come to the laboratory, work out the dilution you will need to do to prepare 50.0 mL of each of these solutions by dilution of a 0.050 M NaOH stock solution. Also before coming to lab, review the Vernier spectrometer's operating instructions before (see the [Operation of the Vernier LabQuest 2](#) section).

Prepare your assigned concentration in a beaker using graduated cylinders for measuring the stock solution and the required volume of water. **CAUTION!** The NaOH solution is caustic. Stir the solution and cover the beaker with a watch glass.

Obtain approximately 40 mL of the stock CV solution in a clean and dry beaker. **CAUTION!** Crystal violet will stain skin and clothing! Record the concentration of CV (it should be around 2×10^{-5} M).

Prepare the LabQuest to collect absorbance data as described in the [Operation of the Vernier LabQuest 2](#) section (the spectrometer and the LabQuest should already be connected and ready to go when you arrive, if not, please consult your instructor). Be sure that the LabQuest is in the *Full Spectrum* acquisition mode before proceeding. Use distilled water as the reference in the same cuvette that you will use for the rest of your measurements. Dilute 10 mL of the CV solution with 10 mL of distilled water in a small beaker and mix well. Rinse a cuvette three times with small portions of this solution and then fill the cuvette approximately three-quarters full. Remove any bubbles by gently tapping the cuvette with your finger. Under **ABSOLUTELY NO** circumstances are you to tap a cuvette on a table top. Before placing the cuvette in the spectrometer, be sure to thoroughly wipe the cuvette's clear sides with a Kim-Wipe (do **not** use a paper towel). Obtain the solution's absorption spectrum. Use the cursor to determine the point of maximum absorbance, which should be near 600 nm. Record this value in your notebook. Send the absorption spectrum to your laptop/tablet (as described in the [Operation of the Vernier LabQuest 2](#) section to print later (from Excel) and include in your notebook.

Now set up the spectrometer to record the change in absorption at the wavelength of maximum absorption by changing the acquisition mode to *Time Based*. Set the *Interval* to 1 sample per second and the experiment's *Duration* to 200 seconds. It is important that you set the units on the time to seconds to facilitate the sharing of the results among the different groups.

Transfer 10 mL of your NaOH solution to a clean, dry beaker. Add 10 mL of the crystal violet solution and stir the resulting reaction mixture well. Rinse the cuvette three times with small portions of the reaction mixture and fill the cuvette three-quarters full with the reaction mixture. Wipe the outside of the cuvette with a KimWipe and place the cuvette in the spectrometer's cell holder. Tap the *Start* icon on the LabQuest to start the kinetics run. You will need to work

efficiently so as to not lose too much data while you are preparing the cuvette. **CAUTION!** Do not attempt to adjust any of the LabQuest settings during a kinetics run; it will ruin your data.

When the kinetics run is over, discard the reaction mixture as described below, and then clean and dry the beaker in which you ran the reaction. Show your graph of absorbance as a function of time to your instructor for approval. Once your results have been approved, obtain two more kinetics data sets with your assigned $[\text{OH}^-]$ (for a total of three). While you are performing these runs, you can be working up the previous data as described below. Note that one or two groups may need to change the length of their kinetics run, if practical, depending on their $[\text{OH}^-]$. What groups are they and why must they change the length of their runs?

Although you can perform all of the data manipulations on the LabQuest, it is far easier to do the analysis in Excel or LoggerPro. Therefore, transfer your data to your laptop/tablet via Wi-Fi or by using a USB memory device. This will also maximize your efficiency in lab because you can be working up one run while the next one is under way.

Waste Disposal

IMPORTANT! Do not discard your hydroxide or CV solutions until you have prepared all three graphs for all three runs. Once you are convinced that your results are valid, you may discard the solutions containing NaOH or crystal violet by flushing them down the drain with a copious amount of water. If you find that you have taken too much of the crystal violet or NaOH solutions, please see if you can share it with another lab group before discarding it.

Results and Analysis

Prepare three graphs for the first run: the first graph is A_{max} as a function of time, the second is $\ln(A_{\text{max}})$ as a function of time, and the third is $1/A_{\text{max}}$ as a function of time. One of these graphs will give a straight line; from this graph determine order with respect to crystal violet at your $[\text{OH}^-]$ and the rate constant, k_{obs} . Compare your results to other groups in the laboratory. Do you get results that are consistent (that is do all groups have the same graph being linear, does the rate of reaction go up with increasing $[\text{OH}^-]$, etc.)?

Prepare graphs for the other two runs, but you only need to prepare the graph that gave you a linear result for the first run (i. e., if the second-order integrated rate law graph was linear for the first run, then you only need to do a second-order integrated rate law graph for each of the other two runs). Average your three k_{obs} , determine the estimated standard deviation and the uncertainty at the 95% confidence level. Share your results with the rest of the class and then prepare Table 1 in your notebook using the class data.

From the class data at different $[\text{OH}^-]$, prepare a graph of $\log(k_{\text{obs}})$ as a function of $\log[\text{OH}^-]$ to determine the order of reaction with respect to OH^- . You do not need to determine the uncertainty in k .

Print out the absorbance spectrum of CV using Excel with absorbance (no units necessary) as a function of wavelength (units are nm). Set the x -axis so that only the absorbance in the range 450 to 900 nm is displayed. Set the y -axis so that the absorbance near 600 nm is clearly seen (choose the scale so that A_{max} is about three-quarters of the full y -axis).

Table 1. Model table to summarize the class data; your table may have up six rows of data. Note that k_{obs} 's units have been omitted; determine them and place them in your table (like shown for the column containing the $[\text{OH}^-]$).

$[\text{OH}^-]$ (M)	k_{obs}

Conclusions

The conclusion should use the [Outline for Measurement Experiments](#).

Use Table 2 as a model for your *Summary Table*.

Table 2. Model *Summary Table* for this exercise. Note that you will need to determine the units on k and put them in the column heading.

Order with respect to CV	Order with respect to OH^-	k

References and Notes

1. The thermodynamic state functions (ΔH , ΔG and ΔS) can be computed directly from the balanced equation for a reaction because they depend only upon the reactants and the products. In contrast, the rate law depends intimately on the path a chemical reaction takes, and therefore the kinetics of a reaction can only be determined experimentally. So, even if we know the reactants and products of a given reaction, we must do some experiments to determine exactly how the reactants become products.

2. For the reaction where some reactant, R, goes to products (i. e., $\text{R} \rightarrow \text{products}$), the rate law is given by the expression: $\text{rate} = k [\text{R}]^n$. The table below shows how to graph the data of $[\text{R}]$ as a function of time, $[\text{R}]_t$, for each of the possible overall orders of the reaction (n). The table also gives the meaning of the slope and intercept when these graphs are linear and the units on the rate constant, k .

n	Order of Reaction	Rate Law	Dependent Variable	Independent Variable	Slope	Intercept	Units on k
0	0 th Order	$\text{rate} = k$	$[\text{R}]_t$	t	$-k$	$[\text{R}]_0$	$\text{M}\cdot\text{s}^{-1}$
1	1 st Order	$\text{rate} = k[\text{R}]$	$\ln[\text{R}]_t$	t	$-k$	$\ln[\text{R}]_0$	s^{-1}
2	2 nd Order	$\text{rate} = k[\text{R}]^2$	$1/[\text{R}]_t$	t	$+k$	$1/[\text{R}]_0$	$\text{M}^{-1}\cdot\text{s}^{-1}$

3. Beers' Law states that for a sufficiently dilute solution, the amount of light absorbed by a chromophore (a chemical species that absorbs light) present in the solution is given by $A = \varepsilon \cdot b \cdot C$, where A is the absorbance (how much light the sample absorbs compared to a solution that does not contain the chromophore), ε is the molar absorptivity (also known as an extinction coefficient; ε depends on the compound and the wavelength of light), b is the pathlength (how much sample the light must pass through) and C is the concentration of the chromophore. According to Beer's Law the amount of color absorbed (and, therefore, the intensity of the color) is linearly dependent on the amount of material absorbing the light. Note that the absorbance has no units (although sometimes 'absorbance units' are used, abbreviated 'AU'). The concentration's

unit is molar, M, and the path length's unit is usually cm. Therefore, the unit of the molar absorptivity is $\text{M}^{-1} \cdot \text{cm}^{-1}$.

4. Holmquist, D. D. and Volz, D. L. *Chemistry with Computers: Using Logger ProTM*; Vernier Software: Portland, OR, 1997, p. 30-1 ff.

Spectrophotometric Determination of an Equilibrium Constant

B. K. Kramer, B. D. Lamp, D. L. McCurdy* and J. M. McCormick

Last Update: April 21, 2011

Introduction

Typically acid-base indicators are themselves weak acids or bases whose acid and base forms have different colors in solution. As the result of the reaction with excess titrant, we convert one form to the other causing a color change that indicates the endpoint of a titration. If we represent the indicator's acid form as HIn and its basic form as In^- , then the following equilibrium describes the chemical reaction that occurs as the $[\text{H}^+]$ is changed. If HIn and In^- have different colors, then the solution's color will change as a function of $[\text{H}^+]$ depending on which of the compounds is present in the greater amount.



The acid dissociation equilibrium constant (K_a) for the indicator that describes this reaction is given by Eqn. 1, in terms of the concentrations of the hydrogen ion, In^- and HIn . Because we are working in aqueous solution, it is convenient to rearrange Eqn. 1 to Eqn. 2 by taking the negative base-ten logarithm of both sides.¹ By convention, we use the prefix "p" to denote the negative base-ten logarithm, and so $-\log K_a$ becomes $\text{p}K_a$ and the $-\log[\text{H}^+]$ becomes pH . This leads to a very convenient way of writing the very small $[\text{H}^+]$ that occur in aqueous solution as numbers that generally fall between 0 and 14. We can then rewrite Eqn. 2 as Eqn. 3 (which some of you may recognize as simply another version of the Henderson-Hasselbach equation). Note that Eqn. 3 predicts that the indicator's $\text{p}K_a$ corresponds to the pH of an indicator solution when the logarithmic term equals zero (i. e., when $[\text{In}^-]$ equals $[\text{HIn}]$).

$$K_a = \frac{[\text{H}^+][\text{In}^-]}{[\text{HIn}]} \quad (1)$$

$$-\log(K_a) = -\log[\text{H}^+] - \log\left(\frac{[\text{In}^-]}{[\text{HIn}]}\right) \quad (2)$$

$$\text{p}K_a = \text{pH} - \log\left(\frac{[\text{In}^-]}{[\text{HIn}]}\right) \quad (3)$$

A convenient way to determine the equilibrium constant of a reaction involving colored species and H^+ is to use absorbance spectroscopy. If we monitor a wavelength at which either one of the two species strongly absorbs we will see the absorbance as a function of pH change as that species' concentration in solution changes. From the equilibrium between HIn and In^- , given above, and considering Le Chatelier's principle, we can see that when the $[\text{H}^+]$ is large (low pH),² the equilibrium will shift completely to the left and the indicator will be completely in the HIn form.

Consider the experiment whose results are shown in Fig. 1. In this experiment, the absorbance of a solution of the indicator bromthymol blue is measured as the solution's pH is varied. As the solution pH changes from acidic to basic, we observe an evolution from spectra where the peak

centered around 430 nm (denoted by λ_1) is largest, to spectra where the peak at about 630 nm (denoted by λ_2) dominates. Further examination of the spectra reveals one wavelength, noted by a star, where the absorbance is essentially independent of pH. This is called an isosbestic point and results from the both forms of the indicator (HIn and In^-) having the same molar absorptivity at this wavelength. It is not uncommon to see one or more isosbestic points in a set of spectra.




Figure 1. Spectrum of the indicator bromothymol blue as a function of pH. The isosbestic point is marked with a *.

A closer examination of the trends in the spectra is shown in Fig. 2. When the solution is very acidic, as on the left side of Fig. 2, all of the indicator is in the HIn form, resulting in a large absorbance at λ_1 (labeled A_{max, λ_1}) but a small absorbance at λ_2 (since $[\text{In}^-]$ is small). At high pH, all of the indicator is in the In^- form giving a strong absorbance at λ_2 (labeled A_{max, λ_2}) and minimal absorbance at λ_1 . As the pH changes from acidic to basic, the position of the equilibrium changes such that HIn is converted to In^- in accordance with Eqn. 3. This conversion results in a decrease in [HIn] and a corresponding increase in $[\text{In}^-]$. Since the absorbance at each wavelength is directly proportional to concentration, we observe a decrease in the absorbance at λ_1 (because [HIn] decreases), and an increase in the absorbance at λ_2 (because $[\text{In}^-]$ increases). It is important to realize that even though we may only be collecting absorbance values at two wavelengths in our experiment, the entire spectrum is undergoing a transformation as we cause the relative amount of HIn to decrease and In^- to increase by increasing the pH of the system. In principle, we could choose any wavelength where HIn absorbs for λ_1 , and any wavelength where In^- absorbs for λ_2 in our analysis. In practice, however, it is beneficial to choose the wavelength of maximum absorbance for each species. Doing so provides the greatest sensitivity (ability to distinguish small changes in absorbance) for our measurement.

From Eqn. 3 it should be obvious that the pH where $[\text{HIn}] = [\text{In}^-]$ corresponds to the indicator's $\text{p}K_a$. This occurs when exactly half of the indicator is in the HIn form and half is present as In^- . In

terms of the experiment, this corresponds to the pH where the absorbance for each form is one half of its maximum, as shown by the dotted line in Fig. 2. Consequently, the pK_a of an indicator corresponds to the pH of the solution at the inflection point in a plot of absorbance as a function of pH.³




Figure 2. Dependence of absorbance HIn and In^- on pH. The red line shows how the absorbance changes at a wavelength (λ_1) where the acidic form of the indicator, HIn, absorbs strongly. The blue line indicates the behavior at a wavelength (λ_2) where In^- , the basic form of the indicator, absorbs strongly. The pH at which the inflection point in both lines occurs is the indicator's pK_a .

$$pH = \log \left(\frac{[In^-]}{[HIn]} \right) + pK_a \quad (4)$$

In practice it is difficult to precisely and accurately determine the inflection point in curves of this type. To get a more precise measure of the pK_a , Eqn. 3 is rearranged to give Eqn. 4. This equation gives a straight line when the solution's pH is graphed as a function of $\log([In^-]/[HIn])$. The slope of this line should be +1 and the y-intercept, where $\log([In^-]/[HIn])$ is zero (i. e. $[In^-] = [HIn]$), is the pK_a , as shown in Fig. 3.

To use Eqn. 4 to determine the pK_a of the indicator, it is necessary to know the pH of solutions that have different ratios of the two indicator species HIn and In^- . Since the pH of the solution will determine the amount of the total indicator that will be in each form, it can be difficult to control exactly how much of the indicator exists as In^- and HIn. It is possible, however, to use absorbance to obtain the ratio $[In^-]/[HIn]$. We will simply monitor the absorbance at two different wavelengths. The first wavelength (λ_1) is chosen where the acidic (HIn) but not the basic (In^-) form of the indicator strongly absorbs radiation. The second wavelength (λ_2) is chosen where the basic but not the acidic form strongly absorbs radiation.




Figure 3. Relationship between pH and $\log([In^-]/[HIn])$ for an indicator. The pK_a of the indicator corresponds to the intersection of the line with the pH axis.

$$A_{\lambda 1} = \varepsilon_{(HIn, \lambda 1)} \cdot b \cdot [HIn] \quad (5)$$

If Beer's law is obeyed, the absorbance at λ_1 is given by Eqn. 5 where $A_{\lambda 1}$ is the absorbance at λ_1 , $\varepsilon_{(HIn, \lambda 1)}$ is the molar absorptivity of HIn at λ_1 and b is the cell path length. Since the amount of the indicator that is in the form of HIn depends on the pH, $[HIn]$ can be difficult, or impossible, to determine. What is known, however, is the total concentration (C_{total}) of the indicator in both forms since a known amount of the indicator was added to the solution at the beginning. At any pH the indicator's C_{total} is given by Eqn. 6. In solutions where the pH is sufficiently low, all of the indicator is in the acidic form, and consequently $C_{total} = [HIn]$. Substituting $C_{total} = [HIn]$ into Eqn. 6 gives Eqn. 7.

$$C_{total} = [HIn] + [In^-] \quad (6)$$

$$A_{\lambda 1}(\text{lowest } pH) = \varepsilon_{(HIn, \lambda 1)} \cdot b \cdot C_{total} \quad (7)$$

Likewise, at λ_2 we can write an equation analogous to Eqn. 6 (Eqn. 8) where we have simply replaced the molar absorptivity of HIn with that of In^- and the concentration of $[HIn]$ with $[In^-]$. In basic solution, $C_{total} = [In^-]$, which leads to Eqn. 9.

The ratio $[In^-]/[HIn]$ at any pH can be obtained by combining Eqn. 5 and Eqn. 8 to give Eqn. 10. Substituting equations Eqn. 7 and Eqn. 9 into equation Eqn. 10 gives Eqn. 11, which simplifies to Eqn. 12.

$$A_{\lambda_2} = \varepsilon_{(In^-, \lambda_2)} \cdot b \cdot [In^-] \quad (8)$$

$$A_{\lambda_2}(\text{highest pH}) = \varepsilon_{(In^-, \lambda_2)} \cdot b \cdot C_{total} \quad (9)$$

So the first step of our experiment will be to determine the value of A_{λ_1} (*lowest pH*) and A_{λ_2} (*highest pH*). A_{λ_1} (*lowest pH*) is the maximum absorbance at λ_1 (A_{\max, λ_1}) in the solution with the lowest pH which should be only due to HIn. A_{λ_2} (*highest pH*) is the maximum absorbance at λ_2 (A_{\max, λ_2}) in the solution with the highest pH and it should be due to only In^- .

$$\frac{[In^-]}{[HIn]} = \frac{\left(\frac{A_{\lambda_2}}{\varepsilon_{(In^-, \lambda_2)} \cdot b} \right)}{\left(\frac{A_{\lambda_1}}{\varepsilon_{(HIn, \lambda_1)} \cdot b} \right)} \quad (10)$$

$$\frac{[In^-]}{[HIn]} = \frac{A_{\lambda_2} \left(\frac{C_{total}}{A_{\lambda_2}(\text{highest pH})} \right)}{A_{\lambda_1} \left(\frac{C_{total}}{A_{\lambda_1}(\text{lowest pH})} \right)} \quad (11)$$

$$\frac{[In^-]}{[HIn]} = \frac{A_{\lambda_2} \cdot A_{\lambda_1}(\text{lowest pH})}{A_{\lambda_1} \cdot A_{\lambda_2}(\text{highest pH})} \quad (12)$$

Until this point we have assumed that at λ_1 the measured absorbance is due only to any HIn present and at λ_2 all absorbance is due to In^- . In fact, the basic form may absorb somewhat at λ_1 and the acidic form may absorb at λ_2 . Because of this, the values of absorbance used in these equations must be corrected to take into account the amount of the absorbance that is due to the other species. To make this correction, we simply subtract the minimum absorbance measured at λ_2 from each of the other measurements made at that wavelength. For example, $A_{\lambda_2} = A_{\lambda_2}(\text{measured}) - A_{\lambda_2}(\text{minimum})$, where the value of $A_{\lambda_2}(\text{minimum})$ comes from the measurement made on the most acidic solution at λ_2 , the wavelength at which it should have the least absorbance. The same correction is to be made for λ_1 by measuring the absorbance of the most basic solution at λ_1 .

Experimental⁴

In this experiment, you will determine the pK_a of bromothymol blue (3',3''-dibromo-thymolsulfonephthalein) by the two methods which have been discussed. At a pH less than 6, the indicator is yellow and at a pH which is greater than 7.6, the indicator is blue. At an intermediate pH, the blue and yellow combine to yield a green solution. You will use the Spectronic 20 Genesys spectrometer (see the [Operating the Spectronic 20 Genesys Spectrometer](#) section for operating instructions) to measure the absorbance of a bromothymol blue solution at the two specified wavelengths as a function of pH. You will also be using a pH electrode interfaced to a LabQuest to measure the pH of each solution that you make (see the [Operation of the Vernier LabQuest 2](#) section for more details). Be sure that you have read and understand the operation of

both instruments before coming to the lab. Note that there will be only two or three spectrometers in the lab for this exercise and you may be sharing with up to three other tables.

You will need to bring your laptop/tablet to the laboratory for this exercise.

Before coming to the laboratory prepare a table in the *Results* section of your laboratory notebook like that shown as Table 1. You will also find it useful to prepare an Excel file for the analysis of the data.

Table 1. Example of a table to record the results of this experiment.

Flask Number	Absorbance at λ_1	Absorbance at λ_2	pH	Solution Color
1				
2				
3				
4				
5				
6				
7				
8				
9				

Pipet 1.00 mL of the bromothymol blue stock solution into each of two 25-mL volumetric flasks. To one of the flasks add 5 mL of distilled or deionized water and 4 drops of concentrated hydrochloric acid; label this flask “Flask 1”. Dilute the solution to the mark with water. The resulting solution should have a pH of approximately 1. To the second flask add 12 drops of 4 M sodium hydroxide solution (**CAUTION!** The sodium hydroxide solution is very caustic.) and fill the flask to the mark with water; label this flask “Flask 9”. The solution should have a pH of about 13.

Prepare the Spectronic 20 Genesys to obtain at 430 nm or at 630 nm according to the instrument’s operating instructions (see the [Operating the Spectronic 20 Genesys Spectrometer](#) section). You should discuss with the other groups sharing a particular spectrometer on how best make your measurements so as to use your time effectively. Obtain and record the absorbance of the bromothymol blue solution at pH 1 (Flask 1) and at pH 13 (Flask 9) at both wavelengths. Be sure that you use the same cuvette throughout this experiment! Remove any bubbles from the cuvette by gently tapping with your finger. Under absolutely no circumstances are you to tap a cuvette on a table top.

Label seven 25-mL volumetric flasks as “Flask 2” through “Flask 8”. Use a pipet to deliver 1.00 mL of the bromothymol blue solution to each of the flasks. Add the volumes of the 0.10 M Na_2HPO_4 solution and the 0.10 M KH_2PO_4 solution to each flask that are indicated in Table 2, below, using graduated cylinders. Dilute each solution to the mark with water.

Set up and calibrate the pH electrode as given in the [Operation of the Vernier LabQuest 2](#) section. Measure the pH of the solutions in all nine flasks, and record the pH of each solution in the table in your notebook. Be sure to copiously rinse the pH electrode with distilled water between each

measurement and pat (do not rub) the electrode dry. When you have finished your measurements, place the electrode in the standard pH 7 buffer solution that you used to calibrate the electrode.

Measure the absorbance of each of the solutions in flasks 2 through 8 at each wavelength. Record the absorbance values in your table.

Table 2. Volumes of 0.10 M KH_2PO_4 solution and 0.10 M Na_2HPO_4 solution to be added to each volumetric flask.

Flask Number	Volume of KH_2PO_4 Solution (mL)	Volume of Na_2HPO_4 Solution (mL)
2	5.0	0.0
3	5.0	1.0
4	10.0	5.0
5	5.0	10.0
6	1.0	5.0
7	1.0	10.0
8	0.0	5.0

Waste Disposal

All solutions generated in this exercise may be flushed down the drain with plenty of water.

Results and Analysis

Input your data from Table 1 into Excel. Subtract the minimum absorbance at each wavelength (the minimum absorbance corresponds to the absorbance of the pH 1 solution at λ_2 and to that of the pH 13 solution at λ_1) from the absorbance of each of the nine solutions at that wavelength. The resulting absorbance values are now corrected for background absorbance. Prepare a graph of absorbance as a function of pH for the nine solutions at each of the two wavelengths. This should look like Fig. 2. Connect the points with a smooth line and determine a value of the $\text{p}K_a$ of the indicator from the inflection point in each plot. It is up to you to determine the best way to find the inflection points.

Using Eqn. 12 and the corrected absorbances calculate the ratio $[\text{In}^-]/[\text{HIn}]$ at each of the nine pHs (also do this in Excel). Graph pH as a function of $\log ([\text{In}^-]/[\text{HIn}])$, as in Fig. 3. From the best-fit line through the data determine the y-intercept and thus the $\text{p}K_a$ of bromothymol blue. Share your value with your laboratory section. Perform a Q -test on the class data, and discard an errant datum, if warranted. Calculate the estimated standard deviation on the class average $\text{p}K_a$ and the confidence limits on the average $\text{p}K_a$ at the 95% confidence level.

Conclusions

Use the [Outline for Measurement Experiments](#) as a guide as you write your conclusions.

Your *Summary Table* should look like Table 3. Be sure to include the confidence interval for the class $\text{p}K_a$ value.

Table 3. Example of a *Summary Table* for this exercise.

Our pK_a	Number of Class Values	Class pK_a

References and Notes

1. The reason that we switch to pH (the negative logarithm of the hydrogen ion concentration) is historical. When the theories that explain equilibria in solution were first developed, there were no computers and so all complex multiplication and division required by the theories was done using logarithms (first with log tables, and then with slide rules). It was simply more convenient to express the small $[H^+]$ involved in aqueous equilibria directly in log form. The convention has been retained, even though it very easy to do the calculations with $[H^+]$ on a pocket calculator.
2. The somewhat counterintuitive relationship between pH and $[H^+]$ is a result of defining the pH as the negative logarithm of $[H^+]$. As the $[H^+]$ decreases, the exponent portion of the number (when written in scientific notation) becomes a bigger negative number. Thus, the pH becomes a larger positive number.
3. You might notice that these plots look very similar to weak acid and base titration curves that you may be familiar with from high school. In the case of a weak acid titrated with a strong base, the dependent variable is the pH and it is measured as a function of the amount of base added. As the acid is titrated, it is converted to its conjugate base. At any point along the curve, the pH of the solution can be said to depend on the ratio of the concentrations of the original acid and its conjugate base. In this sort of titration, the pK_a occurs when half of the acid has been converted to the conjugate base. In Fig. 1 the dependent variable is the absorbance, which is plotted as a function of pH. Again, the pH is a function of the ratio of the concentrations of the acid and its conjugate base. The pK_a also occurs at the point when the amounts of acid and conjugate base are equal. It may appear to be in a different location on the titration curve, but note that it is still halfway to the point where all of the solution present is in the basic form, this time as determined by the absorbance due to one of the species.
4. Braun, R. D. *Introduction to Chemical Analysis*; McGraw-Hill; New York, 1982, pp. 197-199.




Equipment and Techniques

Operation of the Vernier LabQuest 2

B. D. Lamp* and J. M. McCormick

Before using the LabQuest 2, familiarize yourself with the interface by viewing these links: <http://www.vernier.com/products/interfaces/labq2/hardware/> (hardware overview) and <http://www.vernier.com/products/interfaces/labq2/software/> (software overview).

As you use the LabQuest app, you will spend most of your time in one of three screens (or “modes”): *Meter*, *Graph*, or *Data Table*. You can switch between *Meter*, *Graph*, or *Data Table* modes by tapping on the appropriate icon on the top of the screen with the stylus attached to the LabQuest.

Mode	Icon	Available Menus in this Mode
Meter		File, Sensors
Graph		File, Graph, Analyze
Data Table		File, Table

Common Operations of the LabQuest

Powering Up and Configuring the LabQuest

Power the LabQuest on by pressing the power button located on the top of the device.

Plug a sensor, or multiple sensors, into the appropriate port (this can be done before or after powering on). Most sensors used in our chemistry labs are analog sensors that plug into the ports on the left side of the device. If the sensor does not seem to fit in the port; you are using the wrong port. The LabQuest should recognize the sensor automatically and display the sensor's current reading on the screen in *Meter* mode. The sampling mode, sampling rate and duration will also be displayed. Note that switching probes on the fly by unplugging one and plugging in another (especially if one is a spectrometer) may result in neither being recognized by the LabQuest, or one of them living on as a zombie. In this case, you will need to power cycle the LabQuest (turn it off, then on) or reboot the LabQuest by tapping on the *Reboot* button in the *System* folder of the LabQuest *Home Page* (push the house button on the right side of the LabQuest or tap on the house button on the LabQuest).

To configure the sensor, tap **Sensors** to open a drop-down menu that allows the user to adjust a variety of parameters. The most common adjustments will occur in the **Data Collection** section, where the collection mode and acquisition parameters (sampling rate, sample duration, etc.) can be adjusted as needed for the experiment. This can also be done by tapping on **Mode** on the upper right of the *Meter* screen or on the rectangle that shows what sensor is attached, to which port the sensor is attached and its current reading. Tapping on the rectangle will cause a pop-up to appear with the various options, depending on the sensor.


A sensor will display and store a factory-set number of significant figures, even though it records many more. The number of digits to the right of the decimal place displayed and stored can be set by switching to the *Data Table* mode and then tapping on the column corresponding to the measured value whose significant figures or number of digits to the right of the decimal that you want to change. Set the *Displayed Precision* to whatever number that you want and select either *Decimal Places* or *Significant Figures*. Note that these often give similar results, but they are not

the same! This should set the number of displayed (and stored) digits until the sensor is unplugged.

Data Acquisition and Manipulation

To initiate data collection, tap the green arrow (think “play button”) on the bottom left of the screen or press the collect button. During data collection, the green arrow becomes a red square (think “stop button”) and data is plotted graphically in *Graph* mode. Data collection will continue until the experiment duration has elapsed or until the red square icon is pressed.

Once a run is completed, it is wise to save the data. Students may save runs to the device itself and later transfer the data files via e-mail or USB drive. See the section entitled *Data File Transfer* for more details.

Multiple runs may be saved in a single data file. After collection, runs can be stored when in the *Graph* mode by tapping the **Graph** menu and selecting **Store Run**. As a shortcut, a run can be stored by tapping the file cabinet icon () on the right side of the screen. Individual runs (or all runs) can be displayed by tapping the *Run #* button to the left of the file cabinet icon. (Storing a run only keeps the run in memory; you must eventually **Save** the file that contains the run to avoid losing the data).

Curve fitting can be accomplished through the **Analyze** menu in the *Graph* mode. Entire runs or selected portions of runs can be fit to a variety of relationships. Note that it may be more convenient to export the data file and use Logger Pro for such an analysis (see [Basic Data Analysis Using Logger Pro](#)).

To fit only a portion of a run, select the portion (**Graph, Graph Options** from the menu bar at the top of the screen in *Graph* mode) of the run before tapping on the *Analyze* menu. A portion of the data may be selected by first tapping on the graph. Then touch the small circle that appears where you tapped with the stylus, and leaving the stylus on the screen, drag it across to the other end of the data that you wish to select. Tap **Graph, Zoom In** to enlarge it. **Graph, Zoom Out** will return you to the previous view (**Graph, Autoscale once** will return the original view).

Individual runs may be renamed or deleted using the *Data Table* mode. In this mode, tapping on the name of a run will allow you to edit the name. To delete a run, tap on **Table**, select **Delete Run** and choose the run to be deleted. Column calculations can also be done in *Data Table* mode.

Shut Down

For the introductory laboratories at Truman, you will usually not need to shut down the LabQuest at the end of your laboratory period. However, you should check that you have logged off the wireless network before you leave.

You only need to press the power button to shut down the LabQuest (shut-down can also be accomplished by tapping on the Shut Down icon in the *System* folder of the LabQuest’s home screen). You will receive a warning if you have any unsaved data, and you should choose the appropriate option. After power-down is complete, the sensors may be unplugged and they and the LabQuest can be returned to their storage location, as directed by your instructor.

Data File Transfer

Students may use the LabQuest File menu to save your data on the LabQuest. For further data analysis and permanent storage, you will need to either transfer the file to a USB thumb drive or e-mail the file to yourself (or both if you are a little paranoid!). Files will be periodically removed from the LabQuest devices so be sure to transfer your data!

USB Drive Instructions

You can transfer the data file to a USB flash drive by inserting the drive into the USB port on the LabQuest and tapping “File” and “Save”. The screen that appears should have a USB drive icon on the top left of the screen, if it does not, be sure the USB drive has been inserted completely. Tap on the USB icon to select the drive and tap “Save” to save the file to the drive.

E-Mail Instructions

To e-mail the file to yourself, you must first log in to the wireless network. Tap on the Wi-Fi icon next to the battery icon on the bottom right of the screen to enter the *Connections* screen. Tap on the “gear” icon on the top right side of the screen that appears to open the Network Settings screen. On the *Network Settings* screen, select the TrumanSecureWireless network and use your Truman credentials to log in. Once you are logged in, tap “OK” and tap the “X” in the connections window to return to the data screen. From the data screen, tap “File” and select “Email” and “Data File”. In the screen that appears, enter your e-mail address in the “To” field, tap “Done” and “Send”. In your inbox you should receive a message from trumalabquest@gmail.com that has the LoggerPro data file attached. To log out of the wireless network, tap on the Wi-Fi icon again to enter the *Connections* screen and tap on the gear icon in the top right of the screen to open the *Network Settings* screen. Click on the gear icon on at the top right corner of the screen to open the *Network Details* and tap “Forget Network”. Tap *Ok* to exit the screen and then tap on the “X” in the upper right corner to return to the data screen.

Using a Temperature Probe with the LabQuest

Plug the temperature probe into a port on the left side of the LabQuest, as described above.

The LabQuest2 will display a large red rectangle showing what channel the probe is in and the temperature it is reading (note if more than one sensor is attached, there will be other boxes of different colors, one for each of the other probes). The acquisition mode will be set to *Time Based* by default when the LabQuest is powered up, but it will be set for whatever the last acquisition mode was if the LabQuest is already on.

The acquisition parameters can be set as described above. The temperature probes do not require calibration.

See the online user’s manual at <http://www.vernier.com/files/manuals/tmp-bta/tmp-bta.pdf> for more information on the Vernier temperature probe.

Using a Spectrometer with the LabQuest

Set-Up

Plug a Vernier spectrometer, shown in Fig. 1 into the USB port on the left side of the LabQuest (note that Ocean Optics and other mini-spectrometers will not work with the LabQuest). The LabQuest will automatically recognize the spectrometer and display a box labeled USB: Abs.

Calibration

First set the acquisition mode to *Full Spectrum* (although there are many modes available, but the only ones of relevance for a spectrometer are *Full Spectrum* and *Time Based*).



Figure 1. The Vernier SpectroVis Plus spectrometer. The wavelength range of this spectrometer is from ~380 nm – ~900 nm.

From the **Sensors** menu in the *Meter* display on the LabQuest, select **Calibrate** and choose the USB: Spectrometer. This can also be done by clicking on the box labeled USB: Abs and selecting **Calibrate** in the pop-up. With either method, it may take a few seconds for the calibration routine to start; don't get impatient.

Once the calibration routine has started, the LabQuest will wait 90 seconds for the spectrometer to warm up. If the spectrometer has been in use, this step may be skipped by tapping *Skip Warmup*. You will then be asked to place a blank cuvette in the device; tap *Finish Calibration*. When the calibration is finished the LabQuest will display "Calibration completed," at which point you can click *Ok* to return to the data screen.

Obtaining Data

If you are interested in obtaining a spectrum on the full wavelength range you ready to acquire data once the spectrometer is calibrated. Note in the *Full Spectrum* mode you will need to stop acquisition manually by tapping the *Stop* button.

For a time-based measurement, change the acquisition mode to *Time Based*, as described above and enter the rate of acquisition or the interval between data points, the duration of the experiment and the time units, as needed. Click *Ok*. If you are switching between *Full Spectrum* and *Time Based* modes you may receive a warning message about saving your data, if you acquired in the first mode. Choose whichever option you want. The LabQuest will now display the absorbance reading at some wavelength (the format will be USB: Abs @ xxx nm and then below that the absorbance at that wavelength will be given). This wavelength can be changed by tapping on the absorbance reading and selecting **Change Wavelength** in the pop-up. Enter the wavelength and

tap *Ok*. Data acquisition (absorbance at the selected wavelength as a function of time) can be initiated by tapping the green *Go* icon.

For more information on the Vernier SpectroVis spectrometer see the online user's manual at <http://www.vernier.com/files/manuals/svis-pl/svis-pl.pdf>.

Using a pH Sensor on the LabQuest

A pH sensor attaches to the LabQuest in the same way as a temperature probe and the information given in that section applies equally to the pH sensor. Normally, we are only interested in the pH value at the time of measurement, and so it is not necessary to change the acquisition mode or to take data over a set time period (although you could do so, if you wished). However, like the spectrometer, the pH sensor must be calibrated before use.

Calibration a pH Sensor

Once the pH sensor is attached to the LabQuest, it will need to be calibrated before it can be used. Calibrate the pH sensor using two pH buffers (7 and 4 if working in acidic solutions or 7 and 10 if working with basic solutions)

From the **Sensors** menu in the *Meter* display on the LabQuest, select **Calibrate** and choose the pH Sensor.

Immerse the pH Sensor in the pH 7 buffer. Tap *Calibrate Now*. Observe the value under *Live voltage*, once it has stabilized, enter 7.00 for "Value 1" and tap *Keep*.

Remove the pH sensor from the pH 7 buffer, rinse it with distilled water, and immerse it into the second buffer. Again observe the value under *Live voltage*, once it has stabilized, enter the corresponding pH for "Value 2" and tap *Keep*.

Tap *Ok* to exit the calibration routine. The pH sensor is now ready for use. It should be removed from the buffer solution and placed into your unknown solution.

This calibration process may be repeated as needed to assure accurate results.

See the online user's manual at <http://www.vernier.com/files/manuals/ph-bta/ph-bta.pdf> for more information on the Vernier pH probe.

Other Features of the LabQuest

By pressing the home button on the side of the LabQuest or tapping the home icon you can access other LabQuest functions. Most of the settings and apps here should not be modified, but you might find the **Periodic Table** app, the **Calculator** app (in the *Accessories* folder), and the **Quick-Start Guide** app in the *System* folder to be useful. To return to the *Home* screen, tap on the house icon and to return to the data acquisition and analysis modes described above, tap on the **LabQuest App** icon.

Basic Data Collection Using Logger Pro

B. D. Lamp* and J. M. McCormick

Vernier's Logger Pro is a software package that can be used to control a simple data acquisition system attached to a USB port of the computer. This LabPro interface is an analog to digital converter capable of being applied to several measurement scenarios, from temperature and pH measurement, to instrument interfacing. Once data is collected, the Logger Pro software provides several options for basic data analysis. This document provides the basics that you will need to set up the interface, connect a sensor (or sensors) and acquire data. More detailed instructions are available (<http://www2.vernier.com/manuals/LP3QuickRefManual.pdf>).

Hardware Setup

Verify that the LabPro interface is plugged into a USB port on the computer, and that the power supply for the LabPro is plugged in (green light is lit).

Connect the appropriate sensor(s) to one of the four channel inputs on the left side of the LabPro. Depending on the specific measurement being done, additional connections or preparation may be necessary.

Software Setup:

Start the Logger Pro software. The opening screen should appear as shown in Fig. 1. The left portion of the screen is a spreadsheet-like table which will hold data that is collected. Multiple




Figure 1. Opening screen for the LoggerPro software.

runs can be stored in a single table. Explore the *Data* menu to see how. Data can also be copied and pasted to and from Excel for plotting or manipulation. The right portion of the screen will present a graphical display of the data. The button and menu bars across the top of the window contain many useful data handling and analysis functions. It will likely be worth your time to explore these (don't forget the *Help* menu!)

Open a new file by selecting *File, New* from the menu bar.

It is very likely that the LabPro interface and Logger Pro software will have automatically recognized the sensor you have installed on the device. However, it is a wise idea to verify this by selecting *Experiment, Setup Sensors* and choosing the appropriate LabPro interface. Once selected, the window shown in Fig. 2 should appear:





Figure 2. LoggerPro *Setup Sensors* window.

In the example above, a Stainless Steel Temperature probe has been connected to channel 1 (CH1). If your sensor is not automatically detected, simply find it in one of the scroll-down lists, click on the sensor icon and “drag and drop” it on the appropriate channel. Once all sensors are being recognized, select *Close*.

Select *Experiment, Data Collection* or click on the *Data Collection* icon () to define your data collection parameters. The window shown in Fig. 3 will appear:

Select the appropriate sampling mode, duration, and sampling rate for your application. Once setup is complete, click *Done*.

You are now ready to collect data. Near the top of the screen (above the data table), you should see a live display of the signal(s) at the probe(s).






Figure 3. LoggerPro *Data Collection* window.

Data collection can be started and stopped by clicking on the  **Collect** and  **Stop** buttons, or by pressing *Enter* on the keyboard. As data is collected, it is recorded in the table on the left of the screen and plotted on the graph on the right. The table can be used to store several individual runs, the most recently collected data set is called “Latest”. If you intend to save your data for later use, it is important to “store” this run before collecting a new dataset, because it will be overwritten by the new data. Selecting *Experiment, Store Latest Run* places the most recently acquired data into a new column in this spreadsheet (usually named “Run #”). Overwritten data is lost and cannot be recovered. Individual datasets can be renamed by selecting *Data, Data Set Options* and selecting the appropriate data set from the menu.

Basic Data Analysis Using Logger Pro


B. D. Lamp* and J. M. McCormick

Vernier's Logger Pro is a software package that can be used to control the LabPro data acquisition system via the USB port of the computer. The LabPro is what is called an analog to digital converter and is capable of being applied to several measurement scenarios, from temperature and pH measurement, to instrument interfacing. Once data is collected, the Logger Pro software provides several options for basic data analysis. Logger Pro can also be used as a stand-alone data analysis package. Because our introductory chemistry laboratories use LabQuests for data acquisition, only the basics of data analysis with Logger Pro will be described. More information is available at <http://www2.vernier.com/manuals/LP3QuickRefManual.pdf>.

Student Access to LoggerPro

Aside from being installed on the university image, students may install a copy of LoggerPro on their personal computers through our site license with Vernier. Instructions are located at: <http://wp-internal.truman.edu/chemistry/installing-loggerpro/>.

Data Management and Analysis

Multiple data sets can be open in Logger Pro and they will be displayed simultaneously on the screen, but this can be somewhat messy. Data files can be renamed, hidden, or displayed using the *Data Browser*, activated by clicking . Once activated, a window such as the one shown in Fig. 1 will appear. Double-clicking on the name for a dataset will bring up a dialog box with a checkbox for hiding the data set, a place to change the data set's name, as well as a place to enter notes regarding the data. Data that has been hidden is not erased; it is simply not displayed on the graph or in the data table.





Figure 1. The LoggerPro *Data Browser* window.

Many of the data analysis functions, such as integration and linear fitting, require that a portion of the data set be selected prior to executing the function. To select a portion of the data, simply left-click on one end of the desired range and drag the mouse to the other end of the range. The selected area will appear in grey. Only this selected area will be subject to a function such as integration or linear fitting. If no area is selected, the entire data set will be used. If multiple data

sets are being displayed, you will be given the option of applying the function to any or all of the data.

The data files generated by a LabQuest or LoggerPro can only be read by a LabQuest or LoggerPro. To save your file in a format that Excel can read, select *File, Export As* from the menu bar In LoggerPro and then select the *Text* option. It is also possible to cut and paste from LoggerPro directly into Excel and vice versa. Record the file name in your notebook.

Adjusting the Graphical Display

Many of the features of the displayed graph can be adjusted by double-clicking on the appropriate part of the graph. For example, the x and y axes can be adjusted manually by clicking on one of the minima or maxima on the axes and typing in a new value. Double-clicking anywhere else on the graph brings up a *Graph Options* window which presents a wide variety of options for modifying the graph. The arrows on either side of an axis title can be used to “slide” the plot along the chosen axis. To quickly auto-scale the X and Y axes, click the  icon.











Helpful Hints

Right-clicking on an item will probably bring up a menu allowing you to change the properties of that item. If all else fails, try the *Help* menu!

You should save your file frequently to prevent accidental loss.

The data in the spreadsheet can be copied and pasted into Excel.

Below are some useful menu buttons. These buttons, as well as the others below the menu bar provide several useful features; you will find it worth your while to explore these buttons.

Button	Operation (left-click these buttons to)
	Initiate data collection.
	Stop data collection.
	Open the Data Collection setup window.
	Open the Data Browser.
	Turns on the “Examine” function which allows you to use a cursor to examine individual data points.
	Integrates the selected area.
	Calculates a linear fit over the selected data range. Outputs slope, intercept, and correlation coefficient for the linear fit.
	Auto scales the graph axes.
	Zooms graph to the selected area.
	Zooms graph out

Operating Ocean Optics and Vernier Spectrometers with LoggerPro

B. D. Lamp* and J. M. McCormick

LabQuest User's Manual: http://www2.vernier.com/manuals/labquest2_user_manual.pdf

Basic Set Up and Calibration

First, set up the spectrometer and allow it to warm up. Start LoggerPro from the desktop shortcut, or from the Start menu on the lower left-hand side of the screen. If you can't find LoggerPro by either of these methods, please consult your instructor.

LoggerPro will start and the screen shown in Fig. 1 will be displayed. If status line displays the message "No device connected," as it does in Fig. 1, you will need to connect the spectrometer. Do this by selecting *Experiment* from the menu bar and then select *Connect Interface, Spectrometer* and *Scan for Spectrometer*.




Figure 1. Starting screen for LoggerPro.

If the program finds spectrometer, the "No device connected" message will be replaced by the spectrometer icon and the message "Absorbance = x.xxx" where "x.xxx" is some number, as shown in Fig. 2. You will also hear a click, if you are using one of the one-piece spectrometers.





Figure 2. LoggerPro screen that is displayed once the Ocean Optics spectrometer has been found by the software.

Once the computer has located the spectrometer, the software must be configured to acquire the data. Click on the spectrometer icon  Absorbance = 0.102 (or select from menu bar *Experiment*, *Set Up Sensors* and select *Spectrometer:1*). The *Set Up Sensors* pop-up shown in Fig. 3 will





Figure 3. The *Set Up Sensors* pop-up box for the Ocean Optics spectrometers.

appear. You probably won't have to change anything there. However, if you do change something and want to return the spectrometer to its factory-default settings, simply click on the *Restore Defaults* button. Note that if you change something in this box after you have taken a background (*vide infra*), you will need to reacquire the background. The meaning of each parameter in the Set-Up Sensors pop is summarized in Table 1.

Table 1. Summary of the parameters and their purpose in the *Set-Up Sensors* window of LoggerPro when it is in the spectrometer mode.

Parameter	Purpose
Integration Time	This is how long the computer reads the signal coming in from the spectrometer for at each data point. Longer integration times reduce the noise on the data, but reduce how fast one can acquire data.
Wavelength Smoothing	Averages absorbance (or transmission) values at nearby wavelengths. This can help reduce noise in the measurement, but can also lead to distortions of the signal.
Samples to Average	How many individual measurements that the computer will make at each data point and will then average together to get the final value that is displayed.
Wavelength Range	Wavelength (in nm) range over which the computer will record data. The default setting is the range of the spectrometer, but you can select a narrower range to exclude extraneous data.

Before continuing, prepare a cuvette containing your blank.¹ **IMPORTANT!** If you are using a spectrometer that has the detector and light source separated, you **MUST** place a thick piece of paper in the sample holder to block the light beam before proceeding.

In the *Set Up Sensors* pop-up click on the spectrometer icon  and a new pop-up will appear. Your options here are *Calibrate*, *Configure Collection* and *Current Units* (under *Current Units*, *Absorbance* should be selected; do not change). Select *Calibrate*. If you have a one-piece spectrometer, it will click and the pop-up shown in Fig. 4 will be displayed. If you have a spectrometer where the detector and source are separate, you will not hear a click (make sure the light beam is still blocked).




Figure 4. Calibrate Spectrometer pop-up immediately after it has been activated.

After the 1 minute the warm-up is complete, the *Calibrate Spectrometer* pop-up will look like that shown in Fig. 5. If you are using a one-piece unit, simply insert the blank into the spectrometer's sample holder and click *Finish Calibration*. The light travels across the spectrometer in the long direction, so make sure the cuvette is arranged so that the light travels through the unfrosted sides of the cuvette. If you have a spectrometer with the detector and light source separated, first remove the piece of paper that you used to block the beam, then insert the blank into the sample holder and click *Finish Calibration*. **IMPORTANT!** If you don't have a one-piece spectrometer and did not block the beam, click *Cancel* and start the calibration again. Close the pop-up by clicking *OK* when finished.










Figure 5. Calibrate Spectrometer pop-up after warm up is complete.

Simple Absorbance Measurements

If you have successfully complete all of the steps above, the *Collect* button  should now be active and you are ready to acquire absorbance data as a function of wavelength (the default for the spectrometer). Just click on the *Collect* button  to acquire the spectrum and *Stop* button  when you are done. If wavelength is not on the x-axis, please check that the spectrometer is not in the *Absorbance vs. Wavelength* mode, as described below. The function of the other buttons on the LoggerPro tool bar are describe in the [Basic Data Analysis Using Logger Pro](#) section.

Kinetics Measurements

Before configuring the spectrometer for a kinetics experiment, you must first calibrate the spectrometer, as described above.

Open the *Set Up Sensors* pop-up box again, click on spectrometer icon  and select *Configure Collection* or click on the *Configure Spectrometer Data Collection*  icon on the tool bar. The window shown in Fig. 6 will open. Set the data collection mode to *Abs(orbance) vs Time* (for kinetics experiments) or *Abs(orbance) vs Wavelength* (used to obtain a spectrum). When *Absorbance vs Time* is selected, the check boxes under *Full Spectrum* are activated and you can choose which wavelength, or wavelengths, to follow as a function of time. When you check the *Absorbance vs Time* box, the program will automatically select a wavelength to monitor. If you used the *Examine Data* button , the software will default to the wavelength that the cursor was last set to, otherwise the software automatically selects the first point in the list. With *Absorbance vs Time* selected, you will be able to use the check boxes under the *Full Spectrum* heading to

select the wavelength or wavelengths that you wish to monitor. Click *OK* when finished to close this window.





Figure 6. The *Configure Spectrometer Data Collection* window.

Upon closing the *Configure Spectrometer Data Collection* window, the main LoggerPro page should now have absorbance on the y axis and time on the x axis, as shown in Fig. 7. The current absorbance at the wavelength(s) that you specified will be displayed in the lower left-hand portion of the screen.

Now set the data acquisition parameters by clicking the *Data Collection* button  on the tool bar (or from the menu bar by selecting *Experiment, Data Collection*). If you have successfully switched the spectrometer to time-based measurement mode, then the window shown in Fig. 8 will open. If not, then the data collection window for an absorbance as a function of wavelength measurement will appear (Fig. 9) and you will need to reset the collection mode. Note that if you try to change between a time-based measurement and a wavelength-based measurement without changing the data collection mode as described above, then you will see an error message and things will not work properly.

In the *Data Collection* window for a time-based measurement, you need to set four parameters. First, select the *Length* of the experiment (enter the number in the first box and the units in the second). Check the *Sample at Time Zero* check box if you want the computer to measure the absorbance immediately at the start of the experiment (a good idea in most cases). Select *Sampling Rate* either using the slide bar or by typing a number into one of the boxes below the slide bar. Some important things to know about the sampling rate: 1) the units on the sampling rate are the same as those you selected in the *Length* of the experiment, 2) when you type a value in one the sampling rate boxes, the computer automatically calculates the other and 3) if you select

a sampling rate that is faster than the spectrometer can actually do, a message will be displayed and you will not be allowed to proceed.




Figure 7. LoggerPro screen after the data collection mode has been set to absorbance as a function of time.




Figure 8. The *Data Collection* window for time-based measurements.




Figure 9. The *Data Collection* window for wavelength-based measurements.

To have the spectrometer wait before it begins to record data (i. e., set a delay), select the *Triggering* tab and the window should now look like that shown in Fig. 10. Click on the




Figure 10. The *Triggering* control window.

Triggering check box to activate triggering and select either the *On Keyboard* or *On Sensor Value* options. *On Keyboard* means that the computer will wait until you press a key on the keyboard (the *Enter* key is a good choice) before it starts taking data, while *On Sensor Value* means that the computer will wait until the reading from the sensor is a particular value before beginning. Of the two, the *On Keyboard* option is usually more useful.


If you enter “0” in the *Collect ___ Samples before Trigger* box, the software will start acquiring data only after you have started the run (see below) and hit a key on the keyboard. If you enter any other number in the *Collect ___ Samples before Trigger* box, the software behaves somewhat

differently. Once you have started the run, you cannot start acquisition until the number of points that you specified have been acquired. Then and only then can you hit a key to initiate acquiring the rest of the data.

When you finish setting up the collection and triggering parameters for your experiment, click *Done*. You should now be ready to make a run. To initiate the run, push the *Collect* button



. If you selected the *On Keyboard* trigger then the message “Waiting for keyboard trigger . . .” will be displayed in the center of the graph. If you selected to acquire no data before the trigger, simply press a key on the keyboard to start acquisition. If you elected to acquire samples before the trigger, you will have to wait the amount of time required to obtain these data before pressing a key on the keyboard. Once acquisition has been triggered, the “Waiting for keyboard trigger . . .” message will go away and the first data point will be displayed.

To stop data acquisition simply click the *Stop* button .



Notes

1. A blank is usually a solution that is identical in composition to the sample except that it does not contain the chemical species of interest. By running a blank first, we can effectively remove the contribution of everything in the sample that absorbs light, except for what we are interested in.

pH Measurements with Logger Pro

B. D. Lamp* and J. M. McCormick

Set up the Logger Pro A/D converter as described in the [Basic Data Analysis Using Logger Pro](#) section with the pH electrode on Channel 1. When the Logger Pro software is started on the computer, the electrode should be automatically detected and the pH will be displayed on the status line in the upper left-hand corner of the window.

If the pH electrode has been stored for more than 24 hours it will be in a storage solution (prepared by dissolving 10.0 g KCl in 100 mL of standard pH 4 buffer). Unscrew the cap on the storage bottle before removing the electrode from the cap and the rubber gasket. Place the storage bottle where it will not be tipped over.

Rinse the electrode with distilled water and pat dry with a Kim Wipe. Be careful not to rub the glass bulb at the electrode's end. Place the electrode in pH 7 buffer (or the most precise measurements, it is better if this buffer has been freshly-prepared). When the electrode is in the solution, it will be top-heavy. So, be sure that the electrode is secured against tipping.

Before the electrode can be used it must be calibrated. Click on the “pH =” in the status window in the upper left-hand corner of the window. This brings up the sensor configuration window shown in Fig. 1. Click on the pH electrode on channel 1 and select *Calibrate* from the pop-up menu.




Figure 1. Logger Pro sensor configuration window.

The sensor setting window, shown in Fig. 2, will then be displayed. *Live Calibration* is the default setting, and should be displayed as in Fig. 2. If you want to do the calibration with only one buffer, select the *One Point Calibration* check box. If you want to perform the standard two-point calibration, click *Calibrate Now*.




Figure 2. Logger Pro sensor settings window.

Once you have clicked on the *Calibrate Now* button, the window will update to that shown in Fig. 3. If the electrode is immersed in standard pH 7 buffer, and the *Channel Input: (volts)* reading is stable (changes by less than ~0.01 V in 30 sec), you can type in 7.00 in the *Enter Value:* box and click *Keep*. If you want to use another buffer for the first reading, remove the electrode from its current buffer, rinse copiously with distilled water, pat dry and place it in the new buffer. Once the voltage reading has stabilized, enter the pH of the buffer and click *Keep*.




Figure 3. Logger Pro sensor settings window after the *Calibrate Now* button has been pushed.

Once the first buffer's pH has been entered, remove the electrode from the buffer and rinse it copiously with distilled water. Pat the electrode dry and place it in the new buffer. Wait for the voltage to stabilize and then enter the buffer's pH in the *Enter Value:* box under *Reading 2*. Click *Keep* and then *Done*.

The electrode is now calibrated and ready to measure the pH of other solutions. Before transferring the electrode to a new solution, rinse the electrode copiously with distilled water and remove excess water by patting dry.

When you are finished making measurements, the electrode may be stored in the standard pH 7 buffer, if the electrode will be stored for less than 24 hours. If the electrode will not be used again within 24 hours, it should be placed in the pH 4 KCl storage solution.

Operating the Spectronic 20 Genesys Spectrometer

J. M. McCormick

Basic Set Up

The instrument is shown in Fig. 1. Before attempting any measurements, the instrument must be turned on (switch is located on the back near the power cord) and allowed to warm up for at least 10 minutes. When switched on, the instrument will display various messages and perform several self-tests. If the power-up is successful, the default wavelength and the uncalibrated absorbance at that wavelength will be displayed.



Figure 1. Spectronic 20 Genesys spectrometer.








Operation

The keypad for the Genesys is shown in Fig. 2. For simple absorbance measurements, only three keys are used (two to adjust the wavelength and one to set the blank). The control keys and their function are summarized in Table 1.



Figure 2. The Spectronic 20 Genesys control pad.

Table 1. Spectronic 20 Genesys keys and their functions.

Key	Function
	Decreases the measurement wavelength.
	Increases the measurement wavelength.
	Sets the absorbance to 0 or the %transmittance to 100% at the current wavelength.
	Toggles between <u>A</u> bsorbance, % <u>T</u> ransmittance and <u>C</u> oncentration measurement modes.
	Allows the user to set various default parameters. There is no need to use this key for routine operations.
	Print the results. This key is not functional on our instruments because we have no printer.
	Used to set the absorbance at a specific wavelength. Should not be used for routine measurements. If you accidentally press one of the arrow keys, which will activate this mode, simply press the left key (the word Esc will be shown in the display directly above it) to exit.

Once the spectrometer has warmed up, set the wavelength using the keypad. Open the hinged sample compartment (on the lower right side of the instrument) and place a cuvette containing your blank in the sample holder (shown in Fig. 3). Be sure that you have wiped the cuvette carefully with a Kimwipe to remove any drops of solution or other contaminants from the outside

of the cuvette before placing it in the sample holder. Close the sample compartment door and press the 0 ABS/100% T button to blank the spectrometer at your chosen wavelength. The display will read 0.000 A (or 100.0 T). If not, press to blank button again.



Figure 3. Enlarged view of the Genesys' sample compartment showing the sample holder. The light path is from the back of the instrument to the front (top to bottom in this view) as indicated by the arrow on the right side of the sample holder.

You are now ready to make an absorbance measurement on your sample. Simply fill the cuvette with the sample, wipe down the sides and place it in the sample holder. Close the sample compartment door and read the absorbance at your chosen wavelength from the display.

You can continue to make measurements at this wavelength without re-blanking the instrument (although it is a good idea to re-blank the instrument from time to time to account for any baseline drift). However, if you change the wavelength then you will need to re-blank the spectrometer at the new wavelength.

Shutdown

To shut the instrument down, simply remove the sample from the sample holder and turn the power switch to the "off" position.

Melting Point Determination Apparatus

J. M. McCormick



Figure 1. A melting point determination apparatus.

The melting point determination apparatus used in many of our laboratories is shown in Fig. 1 above. It is commonly known by the brand name of one manufacturer, Mel-Temp. A rheostat (power level knob) controls the rate at which the sample is heated and the temperature that is ultimately reached. The sample sits in the sample stage (obscured by the blue protective housing in the figure above) and can be viewed through a magnifying glass mounted on the side of the sample stage. A thermometer is used to monitor the sample's temperature (not shown above for clarity).

The rate at which the apparatus heats the sample, and the temperature which is ultimately achieved, is set by the power level control knob, as shown in Fig. 2. The trick to using a Mel-Temp is to match the power setting to the melting point of the material so that the heating is slow enough to give a precise and accurate melting point in a reasonable amount of time. For example, you would not want to use a power level of "10" to determine the melting point of a substance whose melting point is 100 °C; the temperature would be rising so fast that it would be almost impossible to measure the melting point. Ideally, you would want to choose a power level that levels off just above the substance's melting point, so the melting takes place over a longer time. It is not always practical to take this approach, so we often choose a power level that will achieve the desired temperature in about 10 to 15 min, which is usually sufficient. A rough rule of thumb is to set the power level such that the temperature changes 1 °C per minute near the substance's melting point.

It is also possible to change the power level during the course of melting point determination. For example, choose a higher power setting initially and then decrease it to pass slowly through the melting point. When the melting point of the substance is unknown, it takes some degree of trial and error to determine the melting point accurately.




Figure 2. Heating curves and maximum achievable temperature as a function of the power level for a Mel-Temp melting point determination apparatus.

This page intentionally left blank.

Data Analysis, Recordkeeping and Reporting Results

What You need to do Before Coming to Laboratory

1. Before your first lab meeting, make sure that you have written your contact information in the front of your notebook, and have prepared a “Table of Contents”. See [The Laboratory Notebook](#) section for more information on the laboratory notebook.
2. Be sure that the entries (calculations, graphs and conclusions) for any previous experiments that have been completed are ready to be presented to the instructor. If your class has post-lab questions with its laboratory exercises, don’t forget to answer them and be prepared to turn them in, as your instructor directs you.
3. Answer any pre-laboratory questions for the current exercise. Your instructor will give you further information on how to complete these questions.
4. Begin the new exercise on the first blank page. Write the title of the lab, your lab partner’s name, and the date on the top of the page. Write the lab’s title and starting page number in the “Table of Contents”.

Helpful Hints: read the “Background” and “Procedural Outline” sections of the [The Laboratory Notebook](#) section before proceeding. Use a highlighter to note important steps and materials (this is key to writing the “Background” and “Procedural Outline”). The first few times you may want to write out steps 5 through 7 on loose-leaf paper first (not in the notebook). In this way you can edit what you’ve written before you transfer the final version to your notebook. Do **not**, however, get in the habit of doing it this way; it is too time-consuming. You must quickly learn what you need to have in your notebook to perform effectively in the laboratory. As you are preparing your notebook, always remember that once something is written in the notebook, you must never recopy it (you may edit it).

5. Complete the “Statement of Purpose.” Remember that this is a concise (two sentences maximum) description of what is to be done and how it is to be done.
6. Complete any required “Background” information. This section is usually omitted in CHEM 130 and CHEM 131 exercises, but check with your instructor to be sure. If you are required to have background information, check that all the necessary information is present before writing the next section.
7. Write the “Procedural Outline” (see [The Laboratory Notebook](#) section for formatting). The key here is to summarize the procedure in your own words so that you can do the exercise without the laboratory manual or handout.
8. Reserve space for your “Observations” and “Calculations”. Unless you wish to do some of the calculations before lab, leave these spaces blank.

Your instructor may want to check your notebook, or may require you to turn in the duplicate pages from your notebook, before you begin the laboratory exercise. Be sure to read the course’s syllabus and ask the instructor, if you are unclear on his or her requirements.

The Laboratory Notebook

B. D. Lamp, D. L. McCurdy, A. E. Moody, M. C. Nagan and J. M. McCormick

Introduction

The laboratory notebook is perhaps the single most important piece of laboratory equipment. A scientist's notebook may be directly entered as evidence in court, and as such may be worth millions to a company in patent litigation. While you may never be in a situation where your notebook is worth a million dollars, good record keeping is essential in all scientific research. In an academic laboratory, the consequences of poor record keeping are wasted time as you repeat the experiment, or simply failing the exercise. In an industrial laboratory, inadequate lab records ultimately cost the company money, either in the cost of time and materials or as the result of legal action. In either case, the cost to the responsible employee is their job and all possible future employment. Thus, adequate record keeping will be stressed in all chemistry laboratories at Truman.

There are many different sets of rules for keeping a laboratory notebook,^{1,2} which range from the very elaborate rules followed by industrial chemists to the simplified rules listed below. Not all of the points given here will apply to all courses; your instructor will point out modifications to these procedures in his or her syllabus or in the laboratory. No matter what guidelines you use, the goal is to produce a record of a scientific endeavor that is understandable to a knowledgeable reader and which can be used to repeat the experiment and, presumably, get the same results.

Notebook Format and Rules

Laboratory records are to be kept in a bound notebook (i. e., secured with glue), not a spiral notebook or a loose-leaf binder. The pages are to be consecutively numbered. No pages are ever to be removed (except for the copies produced by duplicating notebooks).

All entries are to be made directly in the notebook in black or blue ink. Everything related to the laboratory work must be recorded in the notebook in an organized and neat manner (if it cannot be easily read, it is not adequately recorded). It is critical that the material is intelligible and understandable to the notebook author and any trained chemist who reads the records, attempts to reproduce these results, or endeavors to finish an incomplete analysis. This concept is often known as "traceable" in the industrial world.

It is **unacceptable** under all circumstances to rewrite (or "copy over") an experiment in the notebook outside of lab. It is also unacceptable to type up portions of the laboratory notebook in a word processor and then attach the printout to your notebook. Plan your activities in the laboratory so that all information is properly entered into the notebook while you are in the laboratory.

Include in the notebook a complete description of the work performed, all reference materials consulted, and ideas that you have related to the work. There should be no loose scraps of paper in the notebook. Graphs, charts, spectra, or spreadsheet analyses should be affixed to the pages of the notebook with tape or glue (to both the original and duplicate pages of duplicating notebooks). Label the space where this material is to go with a description of the item and the results it contained. This way, if it is removed, there will be a record of it. Make no notes on the inserted material.

On the first page of your notebook are written the name of the class, your laboratory section and your name. It is also a good idea to put contact information (e. g., your phone number or email address) here, in case your notebook is lost. The next two pages are reserved for a table of contents (some notebooks come with a table of contents section on the inside front cover, or as a separate, removable sheet; these should not be used). The words "Table of Contents" are to be written at the top of these pages. The first entry is to be the table of contents itself. An entry is made in the *Table of Contents* for each experiment when it is begun. This entry includes the title of each experiment and the page number on which the experiment begins.

You may wish to dedicate one page to a "Preface" in which you describe yourself and the contents of the notebook. Another item that is sometimes included is a page titled "Abbreviations and Other Useful Information". These items must be completed before the first laboratory session.

If you are using a non-duplicating notebook, one usually uses only the right side pages. The left pages are reserved for minor calculations, notes of no consequence to the experiment and notes that refer to material found elsewhere in the notebook. In academic laboratories, especially in teaching laboratories, this rule is relaxed as a cost-saving measure for the student. Please check with your instructor on which protocol he/she wishes you to use.

If a page is skipped, a large "X" must be drawn across it. The page is then initialed and dated. While generally frowned upon, you may skip a line as needed to separate sections. There should be no unused empty space on a page, except for the printed margins. Treat large blocks of blank space like a blank page (this assures the reader that no information was added later).

At the top of each page write the title of the experiment that matches that in the table of contents. At the bottom of the page place the date that the last entry was made on that page, your printed name and signature (or initials).

If an error is made, draw a single bold horizontal line through the error so that it can still be read. Write the correct information to the right of the incorrect entry and have a short accompanying explanation of the reason for exclusion. Never use whiteout or completely obliterate the incorrect entry.

Do not copy any information from the notebooks of former or current students. The only exception is when working in a group, and only one member of the group recorded the data during the experiment. In this case, you must indicate in your notebook that the results were copied from the other person's notebook. Write the recorder's name and the page number from which the data were copied next to the copied data.

In general, the notebook should be arranged in chronological order, so that when one experiment ends the next one begins. In an undergraduate laboratory this is very easy to do, but as you progress in your study of chemistry, things are not always so well-ordered. If you must start a new experiment before another is finished, you simply note on the last page of the unfinished experiment the page on which it will be continued.

Arrangement of an Experiment in the Notebook

Each experiment's record includes the following sections: *Title*, *Statement of Purpose*, *Background*, *Procedural Outline*, *Results*, *Calculations*, *Discussion of Conclusions and Error Analysis*, and *Summary of Results*. Each section should be clearly labeled with the underlined words indicated below. Sign and date each page as it is completed. The *Title*, *Statement of Purpose*, *Background* and *Procedural Outline* sections must be prepared prior to the laboratory period (click here for a checklist of what to do before lab).

Title: This should include the experiment's title, your name, the name(s) of your lab partner(s), and the date the experiment was begun.

Statement of Purpose: Clearly and concisely (two or three complete sentences) describe the purpose of the experiment, including the general method that will be used and anticipated results. Do not begin a *Statement of Purpose* with the phrase "The purpose of this lab is to. .". Don't resort to stock phrases; be somewhat creative. The pedagogical purpose of an exercise is not the same as the *Statement of Purpose*. For the "Determination of Density" exercise in CHEM 130, the pedagogical purpose is to learn about precision and accuracy, and the statistical treatment of data. But your statement of purpose might read "The density of a copper block will be determined by two methods: (1) from its dimensions and mass, and (2) from its mass and volume, as measured by water displacement."

Background: This section contains more information on the goals of the experiment, the methods used and the procedure followed. The content of the *Background* section varies with the type of experiment being performed and the requirements of each laboratory course. Check with your instructor about what to include, but in general the *Background* section must include:

- 1) reference(s) to the procedure that you are using following American Chemical Society guidelines. This reference should contain the full title of the article, or the title of the book and the name of the experiment.
- 2) balanced chemical equations for any chemical reactions that you will be performing. Mechanisms are to be included, when appropriate.
- 3) a table of important physical properties of all the materials (starting materials, solvents, and products) with which you will be working. Be sure that you have thoroughly read the experiment before preparing this table so that it includes all the chemicals that you will use. The following information must be in this table: the name of the compound, its molecular structure, and its molar mass. Other properties that may be important are melting points, boiling points, density, optical rotation, etc, depending on the particular laboratory exercise. Textbooks, laboratory manuals and library references (such as the *CRC Handbook of Chemistry and Physics*, the *Merck Index*, and the *Aldrich Catalog of Fine Chemicals*) are good sources of information on chemicals and their properties. There are some internet resources that also contain the same material. Care should be taken in consulting internet sources because there is often no independent scrutiny of these sites.
- 4) record any hazardous properties (flammability, toxicity, etc.) of the substances that you will encounter in the exercise. The *Merck Index* and the *Safety Data Sheet* (SDS) for a

chemical are excellent sources of this information. Both are available from the library or the stockroom.

Procedural Outline: This section is a brief (this section should not be more than one or two pages long, at most), but complete, description of the steps taken to carry out the experiment. It is not a rewrite of the source material (e. g., laboratory manual, textbook or journal article); use your own words. You may use a bulleted list for the steps. At your instructor's discretion, you may not be allowed to bring the source material to the laboratory. So, be sure that your procedure is complete.

Before beginning the *Procedural Outline*, divide the pages that will contain the procedure into two parts by drawing a vertical line on the page, approximately 3/5 of the way across the page from the left-hand margin (many notebooks already have this line drawn for you). Record the procedure on the left-hand side, and any modifications or procedural notes on the right-hand side. You do not record your results on the right-hand side! Results are recorded in the Results section.

Read the experimental section for the exercise before recording any part of the procedural outline in your notebook. This will make writing the outline much easier and minimize errors. As you read, think through the manipulations that are required and re-read sections that indicate particularly hazardous or important steps (usually denoted by "CAUTION!"). Once you are sure of what you are going to do, go back and write out a step-by-step procedure in your notebook.

Results: This section does not need to be completed before you come to the lab, but you may want to prepare blank tables for recording data. Include in this section a listing of the reduced data (e. g., tables), all graphs, spreadsheet results, and spectra. Unlike the *Procedural Outline*, this and all following sections may use the entire right-hand page. A common error is to forget to leave space for the graphs (a hand-drawn graph should take up most, if not all, of the page so as to maximize the results' precision) and other items (e. g., spreadsheet output) that will be prepared as part of the exercise.

All data should be recorded in this section in chronological order. Include all measurements made (with proper units and correct number of significant figures) and any important observations noted when performing the work. When observations are recorded in the laboratory notebook, they are always written in the passive past tense. So instead of "I saw the solution turn green," one writes, "The solution turned green." In general, personal pronouns (e. g., "I," "we") are not used in scientific writing (the overuse of personal pronouns is taken as a sign of arrogance and the passive is thought to sound more objective). The observations are always written in complete sentences.

When possible, set up tables for repetitive data before coming to the lab. Thinking carefully about the data that will be taken should allow you to prepare a data table, which, although difficult to accomplish for the first few experiments, will save time and space in your lab notebook. The use of tables will make it much easier for the reader to assess your methods and results.

Information on the chemicals and instruments used in the experiment are also included in the Results. For a chemical, the name of its manufacturer, its purity, and the lot number of the chemical are recorded, if this information is available (look the bottle's label). It is easiest to record this information when a chemical is first mentioned. For example, "A saline solution was prepared by dissolving 5.00 g NaCl (99.999%, Aldrich, Lot # 56390-BX) in 500 mL of distilled water." The identity of all instruments used must be recorded, preferably including serial number,

model, manufacturer, and any information on the calibration or settings used. Remember that you want to have enough information in your notebook so that you can easily repeat this measurement, if and when necessary (e. g., you find a mistake). If the instrumental data were saved on disk, include the filename(s) with the data (More Info).

Calculations: An example of each calculation performed to reach the final reported answers should be shown with the units clearly shown at each step. For most exercises in a teaching laboratory, only one example of each different calculation needs to be included. Be sure to label each calculation and parallel the order in which the calculations appear in the procedure. You may want to set up the calculations before coming to lab to maximize your laboratory efficiency.

It is sometimes acceptable to include calculations in the Results section as needed. This is usually done in research situations where you need to make a calculation that you did not anticipate at the start of the experiment, but is sometimes allowed in upper-level courses where the laboratory exercises are not “cookbook.” Check your instructor’s syllabus, or ask him/her, for the format that you are to follow in your course.

If you made more than one measurement on the same phenomenon, calculate the average and standard deviation. Perform other statistical analyses as instructed. When an accepted or theoretical value is available, calculate a percent error. Include the output from any programs used to perform these calculations, and the filename under which the data were saved.

Discussion of Conclusions and Error Analysis: Summarize your results paralleling what you set forth in the Statement of Purpose, compare them to the expected results and try to place them in context (either in the larger field of chemistry or what you have done in class). This is not a long section; it may only be two or three pages long in the notebook. The key to a good discussion section is to concisely cover the important points.

Do not write things like “I liked this lab,” “This lab went well” or “This lab was successfully completed,” and do not use personal pronouns. Take your time and put some thought into your conclusions.

The discussion should try to pinpoint various specific sources of error encountered from the standpoint of the most likely determinate and indeterminate errors in the procedure. Once you have identified a source of error in your measurement, evaluate how it affected the result, and then suggest how this error could be minimized or eliminated. Simply attributing everything to “human error” is insufficient, and will be graded accordingly. Some labs won’t have numerical results to discuss, but you can still indicate sources of uncertainty and how they could be, or were, minimized. To help you learn how to organize your discussion, brief outlines for the three types of exercises usually encountered in undergraduate chemistry laboratory exercises are given below. The types are: 1) exercises with a primary focus on measurement, 2) those which are focused on the synthesis of a compound, and 3) those exercises which require you to observe and report on physical phenomena.

Many aspects of the discussion section are the same in all three, but there are subtle differences that you should appreciate. These outlines are meant only as guides; you will need to adapt them to each particular experiment. Some experiments may incorporate components of each of these

three broad categories. In that case, you will need to write a conclusion that combines the three types of discussions.

Outline for Measurement Experiments

A. State the results and associated statistics.

- If an accepted value is known, assess whether the result is accurate (use a calculated percent error).
- Identify what factors lead to a decrease or increase in accuracy.
- Discuss how the accuracy could have been improved.
- Is there evidence for systematic or gross error? What is the source of that error?
- Is the result precise? (use the standard deviation and/or confidence limits)
- State what factors limited the precision (use propagation of error results).
- State what experimental methods or practices maximized the precision.
- Suggest ways that the precision could be improved.
- In the absence of a true value, discuss whether the precision allows you to have any confidence in the accuracy of the result. You may be able to qualitatively assess the accuracy of your results (e. g., if wood floats, then its density must be less than water's, do your data support this conclusion?).

B. Evaluate the experimental procedure.

- Was the procedure sufficient to provide an accurate and precise result?
- If not, how could it be improved?

C. Discuss whether the result(s) is(are) reasonable in comparison to known values or in the context of similar measurements.

Outline for Synthesis Experiments

A. Report the properties of the prepared material and what methods were used in the characterization.

- Do your results match published properties?
- Compare the published and the experimental properties to assess purity.
- In the absence of published values, evaluate the purity based on the material's properties.

B. Report the percent yield.

- Is the yield reasonable? Compare to the literature or others in class.
- Evaluate the factors and experimental techniques that gave a less than, or better than, average yield.

C. Evaluate the utility of the synthetic procedure.

- Does the reaction give the product in high enough yield?
- Is the material sufficiently pure?
- Is it not too complicated or lengthy?
- Suggest improvements.

Outline for Reporting on Physical Phenomena

A. Describe the system that observed and how you probed its properties. Use sufficient detail so that reader can clearly picture the experiment, but avoid being overly verbose.

B. Describe what you saw.

- Did you see what you expected to see, or were there differences?
- Are the results reasonable, based on what you know about chemistry?

C. Try to explain any differences between what you observed and what you expected.

- Was the difference because of your experimental procedure? If so, how could you modify the procedure to change the result?
- Was the difference a result of reactions or other things that you didn't consider initially? If so, how will you need to change your assumptions to accommodate the new data?

Summary of Results: For measurement and synthetic exercises you will need to include a final table summarizing the results of your experiment. For a measurement exercise this table should include each individual value used in the establishment of the average (check with your instructor if you have more than three or four individual values), the standard deviation, and the confidence limit. For a synthetic exercise your summary table should include the percent yield and the measured physical properties of the new substance. Once the Summary of Results has been recorded, sign and date the experiment.

Labels for Products: If you prepared a substance in the exercise, you must place it in a properly labeled bottle and give it to your instructor. The bottle label must include:

Your Name, Instructor's Initials
CHEM xxx, Date Prepared
Name of Chemical
Formula of Compound (Structure is Optional)
Percent Yield
Physical Properties used in Characterization

The other physical properties used in characterization may include: experimental melting point, boiling point, melting point, NMR chemical shifts, IR peaks, etc. Only one or two physical properties need to be listed (check with your instructor as to which one, or ones, to include). The corresponding literature value for each property should also be listed, if known.

References

1. Kanare, H. M. *Writing the Laboratory Notebook*; The American Chemical Society: Washington, D. C., 1985.

2. Eisenberg, A. *J. Chem. Educ.* **1982**, 59, 1045-1046. This article is available as a PDF file at <http://pubs.acs.org/doi/pdfplus/10.1021/ed059p1045> for Truman addresses and *J. Chem. Educ.* subscribers.

Introduction to Statistics in Chemistry

B. D. Lamp, D. L. McCurdy, V. M. Pultz and J. M. McCormick*

Introduction

“Figures often beguile me, particularly when I have the arranging of them myself; in which case the remark attributed to Disraeli would often apply with justice and force: ‘There are three kinds of lies: lies, damned lies and statistics.’” -Mark Twain

Twain’s dyspeptic view of figures and statistics notwithstanding, the predictive power of science derives from quantifying (placing numbers on) physical properties, and then writing mathematical equations to describe the Universe’s behavior (laws and theories).¹ In this way, we can predict the outcomes of experiments that we have not yet done, and thereby test our models of the Universe. Once we have actually done the experiments we will also need a mathematical way of assessing the reliability of our results. It is in the latter application that scientists rely on statistics. The following is intended to provide you with a basic, working understanding of statistical analysis in chemistry. For a more complete treatment of statistics you will want to take STATS 190/2903 and CHEM 222.²⁻⁴

There are three terms that are used by scientists in relation to their data’s reliability. They are *accuracy*, *precision* and *error*. Accuracy is how close a measured value is to the true, or accepted, value, while precision is how carefully a single measurement was made or how reproducible measurements in a series are. The terms accuracy and precision are not synonymous, but they are related, as we will see. Error is anything that lessens a measurement’s accuracy or its precision.

To beginning science students the scientific meaning of “error” is very confusing, because it does not exactly match the common usage. In everyday usage “error” means a mistake, but in science an “error” is anything that contributes to a measured value being different than the “true” value. The term “error” in science is synonymous with “mistake” when we speak of *gross errors* (also known as *illegitimate errors*). Gross errors are easy to deal with, once they are found. Some gross errors are correctable (a mistake in a calculation, for example), while some are not (using the wrong amount of a reactant in a chemical reaction). When met with uncorrectable gross errors, it is usually best to discard that result and start again.

The other types of “errors” that are encountered in science might be better referred to as *uncertainties*. They are not necessarily mistakes, but they place limits on our ability to be perfectly quantitative in our measurements because they result from the extension of a measurement tool to its maximum limits. These uncertainties fall into two groups: *systematic errors* (or *determinate errors*) and *random errors* (or *indeterminate errors*).

A systematic error is a non-random bias in the data and its greatest impact is on a measurement’s accuracy. A systematic error can be recognized from multiple measurements of the same quantity, if the true value is known. For example, if you made three measurements of copper’s density and got values of 9.54, 9.55 and 9.56 g/cm³, you would not be able to determine whether a systematic error was present, unless you knew that the accepted value of copper’s density is 8.96 g/cm³. You might then suspect a systematic error because all of the measured values are consistently too high (although the closeness of the data to each other implies some level of confidence). Often in science one needs to assess the accuracy of a measurement without prior knowledge of the true value. In this case the same experiment is performed with samples where the quantity to be

measured is known. These *standards*, or *knowns*, can reveal systematic errors in a procedure before measurements are made on unknowns, and give the experimenter confidence that they are getting accurate results.

The last type of uncertainty is random error. As the name suggests, these uncertainties arise from random events that are not necessarily under the control of the experimentalist. Random errors can be thought of as background noise. The noise restricts our ability to make an exact measurement by limiting the precision of the measurement. Because indeterminate errors are random, they may be treated statistically.

Assessing Accuracy

Accuracy can be expressed as a percent error, defined by Eqn. 1, if the true value is known. Note that the percent error has a sign associated with it ('+' if the measured value is larger than the true value, and '-' if it is less than the true value). Using the copper density data from above

$$\text{percent error} = \frac{\text{measured value} - \text{true value}}{\text{true value}} \times 100 \quad (1)$$

and Eqn. 1, we can calculate a percent error for each data point of approximately +6.5%. This suggests the presence of a systematic error because, if there were no systematic error, we would expect the percent error for each member of the data set to be very small and that there would be both positive and negative values. When the true value is not known, no conclusion about accuracy may be made using a percent error. In this case, standards must be run or other statistical methods based on the precision can be used. However, the latter can be used only to assess the accuracy of a group of measurements.

In the absence of systematic errors, the average, \bar{x} , of a set of measurements (Eqn. 2) should approximate the true value, as the number of measurements, N , becomes very large (i. e., there are many individual data points, x_i). But if a systematic error is present, then making more

$$\bar{x} = \frac{x_1 + x_2 + x_3 + \cdots + x_n}{N} = \frac{1}{N}(x_1 + x_2 + x_3 + \cdots + x_n) = \frac{1}{N} \sum_{i=1}^N x_i \quad (2)$$

measurements will not make the average approach the true value (as is the case for the copper data we have been discussing). So to make the most accurate measurements (smallest percent error), all systematic errors must be eliminated. Note that the percent error for a set of measurements can be made using the average. The average value of copper's density, using the data that we have been discussing, is 9.55 g/cm³, which has +6.6% error.

Assessing Precision

The range is the simplest, and crudest, measure of the precision for a set of measurements. The range is simply the highest value minus the lowest value, and can be used to get a rough idea of the spread in the data, but not much more. Sometimes you will see a range reported in the form \pm (range/2), which should not be confused with the confidence limits discussed below. To avoid confusion, chemists almost never report the range in this way.

A better measurement of precision for a data set is the standard deviation (σ) which may be calculated using Eqn. 3 for data sets that have more than about 20 points.

$$\sigma = \sqrt{\frac{1}{N} \sum_{i=1}^N (x_i - \mu)^2} \quad (3)$$

In Eqn. 3 μ is the true mean (what the average becomes when N is large). Since it is rare in chemistry to have more than three to five replicate experiments, the estimated standard deviation, S , is used instead (Eqn. 4). In either case, a smaller S or σ indicates higher precision.

$$S = \sqrt{\frac{1}{N-1} \sum_{i=1}^N (x_i - \bar{x})^2} \quad (4)$$

Note the dependence of both S and σ on the number of data points. If the difference terms are all about the same, then the precision should increase (S and σ decrease) as N increases. So, it is statistically advantageous to make more measurements, although this must be balanced with practical considerations. No one wants to do a ten-day experiment 30 times just to get better statistics!

The standard deviation is related to another estimate of precision known as the confidence limit or the confidence interval. The confidence interval is a range of values, based on the mean and the standard deviation of the data set, where there is a known probability of finding the “true” value. A confidence limit is written as $\bar{x} \pm \Delta$ at the given confidence level. For example, a volume expressed as $2.16 \pm 0.05 \text{ cm}^3$ at the 95% confidence level means that there is at least a 95% probability of finding the “true” value in the range 2.11 cm^3 to 2.21 cm^3 (in other words, within $\pm 0.05 \text{ cm}^3$ of the average, 2.16 cm^3). It does not mean that only 95% percent of the time we are confident of the result! To some extent precision is separate from accuracy. However, if enough precise measurements are made in the absence of systematic error, we have increased confidence that our average is a good approximation to the true value, even though we do not know the true value. So, a confidence limit also expresses a level of certainty that the true value lies within $\pm \Delta$ of the average, in the absence of systematic error.

To determine a confidence limit, the uncertainty, Δ , must first be calculated from the estimated standard deviation using Eqn. 5. The value of t (called Student’s t) in Eqn. 5 may be calculated

$$\Delta = t \frac{S}{\sqrt{N}} \quad (5)$$

in Excel using the TINV function, or may be taken from a table such as Table 1, which gives the value of t for various degrees of freedom (usually the number of data points minus one, i.e., $N - 1$) at the 95% confidence level. Note that as the precision of a set of measurements increases, Δ will decrease at a set confidence level. Higher confidence levels also reflect higher precision in the data set.

Table 1. Values of t at the 95% confidence level for various degrees of freedom.

Degrees of Freedom	t	Degrees of Freedom	t
1	12.7	7	2.36
2	4.30	8	2.31
3	3.18	9	2.26
4	2.78	10	2.23
5	2.57	15	2.13
6	2.45	∞	1.96

Precision and Significant Figures

In lecture and on exams and quizzes when we write a number, we assume that the precision is ± 1 in the last number written (for example, the number 31.778 would be assumed to have a precision of ± 0.001). We do this for simplicity. Because when we make this assumption we only need to concern ourselves with significant figures and we can ignore statistics and the propagation of error (see the [Propagation of Uncertainty](#) section for a more thorough discussion of error propagation). In real life we are not so lucky and we must worry about significant figures, statistics and the propagation of error. However, significant figures are always our first step in analyzing our data.

The uncertainty in a number tells us directly how many significant figures our result has. This is because the uncertainty tells us in what place the first uncertain digit is (or you could say it is the first digit where certainty ends). For example, if you had a result 15.678 ± 0.035 kJ/mole at the 95% confidence level, then you could tell from the uncertainty that the first digit that has any uncertainty in it is the tenths place. We know the 1, the 5 and the 6 (and are confident that we know them), but the 7 we have some doubt about. We only really know this digit to ± 3 at 95% confidence and the hundredths place is not known with any certainty. How we show this is discussed below.

Reporting Results

There are three ways in which the statistical information that accompanies a measurement (average, standard deviation, and confidence limit) can be stated. If, for example, five replicate measurements of a solid's density were made, and the average was 1.015 g/cm^3 with an estimated standard deviation of 0.006, then the results of this experiment could be reported in any of the following ways:

- The average density is 1.015 g/cm^3 with an estimated standard deviation of 0.006 g/cm^3 .
- The density is $1.015(6) \text{ g/cm}^3$.
- The density is $1.015 \pm 0.007 \text{ g/cm}^3$ at the 95% confidence limit.

In this example the density has four significant figures, and the uncertainty is in the last decimal place. Sometimes the uncertainty and the number of significant figures in the measurement do not match. This means that each individual measurement was measured more exactly than the reproducibility within the group. If the standard deviation in the density experiment had instead been 0.010 g/cm^3 , then the results might be reported as:

- The average density is 1.02 g/cm^3 with an estimated standard deviation of 0.01.
- The density is $1.02(1) \text{ g/cm}^3$.
- The density is $1.02 \pm 0.01 \text{ g/cm}^3$ at the 95% confidence limit.

The results have been rounded off because the number of significant figures does not reflect the precision of the data set. In other words, the statistical analysis shows us that the first digit where uncertainty begins is the $1/100^{\text{th}}$ place, even though each measurement was made to the $1/1000^{\text{th}}$ place. The last significant figure is in the $1/100^{\text{th}}$ place, so this is where rounding occurs. Sometimes the average and the uncertainties are quoted to the maximum number of significant figures (i. e., $1.015(10) \text{ g/cm}^3$). In this way the precision of each individual measurement and the precision of the set of measurements are shown.

Using Statistics to Identify Hidden Gross Error

Another way in which statistics can be used is in the evaluation of suspect data by the Q -test. The Q -test is used to identify outlying (“bad”) data points in a data set for which there is no obvious gross error. The Q -test involves applying statistics to examine the overall scatter of the data. This is accomplished by comparing the gap between the suspect point (outlier) and its nearest neighbor with the range, as shown in Eqn. 6. The calculated Q is then compared to the critical Q values, Q_c , at given confidence level, like those in Table 2. If the measured Q is greater than Q_c , then that data point can be excluded on the basis of the Q -test.

$$Q = \frac{|\text{suspect value} - \text{closest value}|}{\text{highest value} - \text{lowest value}} \quad (6)$$

Table 2. Critical Q (Q_c) values at the 90% confidence limit for a small number of data points, N .

N	3	4	5	6	7	8	9	10
Q_c	0.94	0.76	0.64	0.56	0.51	0.47	0.44	0.41

For large data sets ($N > 10$) a data point that lies more than 2.6 times S (or σ) from the average may be excluded. Although for medium-sized data sets (between 11 and 15 data points), there is an alternative treatment that is usually sufficient. In these cases, we can use Q_c for $N = 10$, but in doing so, a higher criterion is placed on the data for exclusion of a point than is required by statistics. So, an outlying point that could have been discarded is retained and the precision is quoted as being less than it actually is. But again, it is better to err on the side of caution in our data treatment.

In any case, only one data point per data set may be excluded on the basis of the Q -test. More than one point may be tested, but only one may be discarded. For example, you have measured the density of copper as 9.43, 8.95, 8.97, 8.96 and 8.93 g/cm^3 ; can any of these points be excluded?

First, we must remember that the Q -test is only valid at the extremes, not in the middle of the data set. So before performing a Q -test, it is best to sort the data (as already been done with the data that we are considering). Now look at the extremes and see whether either of the points look odd.

In this case, the low value (8.93 g/cm³) is not that much different than the values in the middle of the set, while the high value (9.43 g/cm³) looks to be suspect.

Having decided that the 9.43 g/cm³ value is suspect, we can calculate Q using Eqn. 6, (suspect value = 9.43, closest value = 8.97, highest value = 9.43 and lowest value = 8.93). This gives $Q = 0.92$ for this point. Since this exceeds Q_c for five data points (for $N = 5$, $Q_c = 0.64$ in Table 2), this point may be excluded on the basis of the Q -test. The Q -test may **not** be repeated on the remaining data to exclude more points.

One last important thing about the Q -test is that it cannot be performed on identical data points. For example, if our data set had been 9.43, 9.43, 8.95, 8.97, 8.96 and 8.93 g/cm³, we would not have been able to use the Q -test on the 9.43 g/cm³ values.

Propagation of Uncertainty

So, now we have an average and an associated uncertainty at given confidence level for a data set. What happens if we use this result in a calculation? The simple answer is that the uncertainty carries through the calculation and affects the uncertainty of the final answer. This carrying over of uncertainty is called propagation of error, or propagation of uncertainty, and it represents the minimum uncertainty in the calculated value due entirely to the uncertainty in the original measurement(s). The equations that describe how the uncertainty is propagated depend on the calculation being done, and can be derived using calculus.⁵ If you would like to learn more about how Eqn. 7 was derived, please see the [Propagation of Uncertainty](#) section. The following example demonstrates how a propagation of uncertainty analysis is done.

The dimensions of a regular rectangular wood block are 15.12 cm, 3.14 cm and 1.01 cm, all measured to the nearest 0.01 cm. What is the volume and the confidence limits on the volume based on this single measurement? The equation for the uncertainty in the volume is given in Eqn. 7, where ΔV , Δx , Δy and Δz are the uncertainties in the *volume* and the x , y and z

$$\Delta V = \pm V \sqrt{\left(\frac{\Delta x}{x}\right)^2 + \left(\frac{\Delta y}{y}\right)^2 + \left(\frac{\Delta z}{z}\right)^2} \quad (7)$$

dimensions, respectively. Do not be confused by the notation! The Δ represents the uncertainty, not a change, in these parameters. Since each measurement was made to the nearest 0.01 cm, $\Delta x = \Delta y = \Delta z = \pm 0.01$ cm. First we calculate the volume, being careful with our significant figures (note the extra “insignificant” figures from the calculator output, shown as subscripts, carried along in the calculation for rounding purposes).

Substituting the known values of V , x , y , z , Δx , Δy and Δz into Eqn. 7 gives

$$\Delta V = \pm 47.9_{52} \text{ cm}^3 \sqrt{\left(\frac{\pm 0.01 \text{ cm}}{15.12 \text{ cm}}\right)^2 + \left(\frac{\pm 0.01 \text{ cm}}{3.14 \text{ cm}}\right)^2 + \left(\frac{\pm 0.01 \text{ cm}}{1.01 \text{ cm}}\right)^2}$$

$$\Delta V = \pm 47.9_{52} \text{ cm}^3 \sqrt{(4.374 \times 10^{-7}) + (1.014 \times 10^{-5}) + (9.803 \times 10^{-5})}$$

$$\Delta V = \pm 47.9_{52} \text{ cm}^3 \sqrt{1.086 \times 10^{-4}}$$

$$\Delta V = \pm 47.9_{52} \text{ cm}^3 (1.042 \times 10^{-2}) = \pm 0.5 \text{ cm}^3$$

So, the volume would be reported as $48.0 \pm 0.5 \text{ cm}^3$ for the single measurement, and this represents a minimum uncertainty in the volume based on the uncertainties in the block's dimensions. Note that the propagated uncertainty usually has only one significant figure.

To see how the propagated uncertainty differs from an uncertainty for a population (data set), imagine that we did this measurement three times and got volumes of 48.1, 47.8 and 48.3 cm^3 . Each individual measurement has an uncertainty of $\pm 0.5 \text{ cm}^3$, from the propagation of uncertainty analysis, but the uncertainty for the set of measurements is $\pm 0.7 \text{ cm}^3$. This was calculated with $S = 0.3 \text{ cm}^3$ (determined using Eqn. 4) and the value of t taken from Table 1 (for $N - 1 = 2$) by substitution into Eqn. 5. Thus, the volume would be reported as $48.1 \pm 0.7 \text{ cm}^3$ at the 95% confidence limit. Notice that the uncertainty in the population is not the same as the uncertainty in each individual measurement. They are not required to be the same, nor are they often the same. In this example, the propagated uncertainty is less than that for a series of volume measurements, indicating another source of uncertainty besides that arising from the uncertainty in the block's dimensions. This is often the case, and in your conclusions to an exercise or experiment you should try to identify its source and discuss its impact on your result.

Regression Analysis

Once we have data from an experiment, the challenge is to determine the mathematical expression that relates one measured quantity to another. The problems that confront us when we attempt to mathematically describe our data are 1) how to establish the mathematical formula that connects the measured quantities and 2) how to determine the other parameters in the equation. The process by which a mathematical formula is extracted from a data set is called fitting, or regressing, the data.

A linear relationship is the simplest, and most useful, mathematical formula relating two measured quantities, x (the independent variable) and y (the dependent variable). This means that the equation takes the form $y = m \cdot x + b$, where m is the slope of the line and b is the intercept. It is possible to relate two quantities with other equations, but unless there is a good theoretical basis for using another function, a line is always your best initial choice. For a linear relationship the values of m and b must be found from the data (x and y values), which is done through a linear least squares regression (or fit). The mathematics behind the fitting algorithm is not relevant at this time, but it is important to know that the least-squares procedure assumes that the uncertainty in the x values is less than the uncertainty in the y values. This means that, if we want to get a meaningful slope and intercept from our fit, we must make the measured quantity with the smallest uncertainty be the independent variable.

Some of pitfalls that you may encounter when performing a regression analysis (and why it is always a good idea to graph your data) have been discussed by Anscombe.⁶

References

1. Paul, R. L.; Elder, L. *A Miniature Guide for Students and Faculty to Scientific Thinking*; Foundation for Scientific Thinking: Dillon Beach, CA, 2003.

2. Devore, J.; Peck, R. *Statistics: the Exploration and Analysis of Data*; Duxbury Press: Belmont, CA, 1993.
3. Skoog, D. A.; West, D. M.; Holler, F. J.; Crouch, S. R. *Analytical Chemistry: an Introduction*; Harcourt College Publishers: New York, 2000.
4. Harris, D. C. *Quantitative Chemical Analysis*; W. H. Freeman: New York, 2003.
5. Andraos, J. J. *Chem. Educ.* **1996**, 73, 150-154. This article is available as a PDF file at <http://pubs.acs.org/doi/pdf/10.1021/ed073p150> for Truman addresses and *J. Chem. Educ.* subscribers.
6. Anscombe, F. J. *The American Statistician* **1973**, 27, 17-21. This article is available as a PDF file on JSTOR (http://www.jstor.org/stable/2682899?seq=1#page_scan_tab_contents). The Wikipedia entry for Anscombe's quartet summarizes the results nicely (available online at https://en.wikipedia.org/wiki/Anscombe%27s_quartet).

Preparing Graphs

J. M. McCormick

General Considerations¹

When we prepare a graph the independent variable is always on the x-axis, and the dependent variable is always on the y-axis. We indicate which variable is which by saying “as a function of” or “versus”, with the dependent variable coming first, and the independent variable coming second. So if someone says, “volume was plotted as a function of mass” or “the volume is plotted versus mass,” it means that mass was on the x-axis and volume was on the y-axis. Watch what you say/write as there is only one correct usage!

All graphs must have axis labels. The axis labels have two parts: the first is the name of the parameter, and the second is the unit. The axis labels can have one of three formats: parameter name (unit); parameter name, unit; parameter name/unit. For example, if the parameter was temperature and it was measured in Kelvin, then the axis label could be Temperature (K), or Temperature, K or Temperature/K.

Devise a scale for each axis so that the tick mark labels end in a 0 or a 5. If this is not possible, use a scale so that the last digit in the tick mark labels is an even number. Tick mark labels with the last digit being an odd number, other than 5, are almost never used. The axes do not need to start at zero. For example, if all the x values occur between 400 and 600 nm, a graph of these data could start at 400 nm. The key to preparing good graphs is selecting a scale that shows all of the data and minimizes large regions of blank space. The one exception is when you need to extrapolate back to a certain value, but the data are not necessarily close to that value. In this case, you would want to have the value to which you will extrapolate shown on the graph, even though there may be some blank space.

The number of significant figures in the tick marks is usually less than that in the original data. For example, if our temperature axis has values between 200 and 320 K with each individual point measured to the nearest 0.01 K, we would not write the tick marks as 200.00, 220.00, etc., but rather 200, 220, etc. (see Fig. 1 and Fig. 2).

Graphs that will appear in a notebook must also have a title. The title should concisely tell the reader what is in the graph. Avoid the obvious use of vs. or versus or the word “plotted” in the title. Graphs that will appear as a figure in a publication or in a formal laboratory report will not have a title (the information is given in the figure caption).

Hand-Drawn Graphs

For hand-drawn graphs in the notebook choose a scale so that the graph fills most, if not all of the page. Use a small, dark dot for each data point. Draw a circle around the point to help the reader locate the actual data point. The radius of the circle usually approximates the uncertainty in the point unless this gives a circle that is too large. With hand-drawn graphs, one usually does a linear regression “by eye”, which means that a ruler is used to put a line through the data such that all points lie as close as possible to the line. The slope and intercept can then be determined from the line itself using points on the line at the extremes of the graph to maximize the precision.




Figure 1. A properly executed hand-drawn graph.

Spreadsheet Graphs

With a spreadsheet it is much easier to prepare graphs, but it is also much easier to produce a poor quality graph. To assure that your graphs are correctly prepared (e. g., look good and are easy to understand by the reader), follow these standard procedures:

- Most of the graphs that you will prepare in a chemistry class are called *XY Scatter* plots in Excel. The other formats will be used sparingly because they are generally not useful in chemistry.
- When creating the graph, make the chart occupy a new sheet; do not create it in the worksheet containing the data.
- Make all lines (the axes, the tick marks, the box around the plot area and any data lines) black with medium intensity (third choice from the top under *Weight* in the format boxes).
- The plot area should have no fill color (not the default setting) and no box around it.
- Turn off all gridlines (not the default setting), unless told otherwise by your instructor.
- Legends showing what a given symbol means are useful, but not always needed (e. g., in a formal laboratory report this information is given in the caption that accompanies the figure). Be sure that your selection of lines and legend titles clearly distinguish between multiple data sets and fits.

-A title should be placed at the top of the graph if the graph is to be placed in the laboratory notebook. This helps the reader immediately know what the graph is. The title should be a concise description of what is being graphed (e. g., “Pressure as a Function of Temperature for Nitrogen”). Usually you do not need to describe in the title the units used in the graph, but there are some instances where this is necessary. Avoid using the words “versus”, “vs.” and “plotted” in the title. The title is **omitted** when the graph will be used as a figure in a publication or formal laboratory report. This is because the information normally put in the title will be included in the figure caption.

-Watch your tick marks. The spreadsheet automatically selects the spacing, which may not be appropriate for your graph (see *General Considerations*).

- Be careful with the use of color to distinguish different data sets on the same graph. Yellow and light blue do not show up very well when printed either on color or black and white printers. Do not distinguish different data sets by color if you do not have a color printer. In general, it is best to dispense with color entirely and make all lines and symbols black (or at least a uniform dark color).

-Data are always shown as symbols and fits to the data are shown as lines or curves. Do not connect the data points with lines. This will confuse the reader as to whether these lines represent a fit, or not.

-For graphs that will be placed in a notebook, you can include the equation of a best-fit line and the R^2 value for the fit in a legend (but remember that this information should also be written in the notebook as part of the graph’s description, in case the graph is removed). However, for graphs that will be submitted for publication or used in a formal laboratory report, this information is not shown on the graph itself. Rather, it is placed either in the figure caption or in the body of the text itself.

-Sometimes a fit is not required, or cannot be done, but you still want to show a trend in the data. In this case a smooth line that passes through the data as an “aid to the eye” is used, and is so indicated in the text accompanying the graph.

-A 16-point font for axis labels and graph titles and a 14-point font for the tick mark labels are better choices for general work than the default settings. The font used for the title should match that used for the axis labels.

-Depending on the final size of the graph, you might want to increase the size of the symbols to 9 point.

-Most spreadsheets automatically select the number of decimal places on the tick labels from the data, but not always correctly. For example, if your data were between 50 and 100, you would want tick labels of 50, 60, 70, 80, 90, and 100. But if the data in the spreadsheet are set to two decimal places, most spreadsheets would make the labels 50.00, 60.00, 70.00, 80.00, 90.00 and 100.00. Although this is not necessarily incorrect, the convention is to show the minimum number of decimal places on the tick labels.

Remember that the default settings in a spreadsheet may not fit your needs. Therefore, it is critical that you know how to change the settings (see the [Guide to Excel](#) section).




Figure 2. A graph that was properly prepared for a laboratory notebook using a spreadsheet.

References

1. Dodd, J. S., Ed. *The ACS Style Guide: a Manual for Authors and Editors*; The American Chemical Society: Washington, DC, 1986.

Guide to Excel

K. N. Carter, B. D. Lamp, D. L. McCurdy and J. M. McCormick*

Updates by S. Warner and V. Pultz

Introduction

The purpose of this section is to give the average student a brief introduction to Excel and provide a quick reference to those students who are more familiar with the program. It supplements, but does not replace, hands-on experience. This tutorial was written for Excel 2010, but other versions are similar, although where you find functions on the menu bar are different. If you encounter trouble, use the online help, or ask a more experienced user.

Cell References

All cells in a spreadsheet are referred to by their column (letter) and row (number). For example, cell A1 is the first cell in the upper left-hand corner of the worksheet. Many times it is easier to enter a formula once and then copy and paste it into other cells. When we do this, Excel will automatically change the row and column references. If we don't want the reference to a cell to change, we need to use an *absolute reference*. An absolute reference is denoted with the "\$" symbol, and may precede the column or row designation, or both. For example, if we wanted to subtract the value in cell B13 from that in cell A2, but keeping the reference to cell B13 constant, we might enter in cell C2 the formula " $=A2-\$B\13 ". If we then copy cell C2 and paste it into cell C3, cell C3 would contain the formula " $=A3-\$B\13 ". But if we pasted it into cell D2, cell D2 would contain the formula " $=B2-\$B\13 ".

Basic Mathematical Operations

Formulas may be entered with any combination of numbers and cell references using the standard symbols shown in Table 1. Note that all formulas begin with the equals sign. So, entering " $=A1+B1$ " in cell C1 would add the contents of the two cells and place that result in C1, while entering " $=A1+5$ " in C1 adds 5 to the value in A1 and places it in C1. The order of operation is important! Multiplication and division are performed before addition and subtraction. Use parentheses to ensure proper calculation. Convince yourself that the order of operations is important by entering " $=(2+3)/4$ " in one cell and " $=2+3/4$ " in another. These calculations do not give the same result!

Table 1. Symbols for basic mathematical operations in Excel.

=	Equals (starts all formulas entered in Excel)
+	Addition
-	Subtraction
/	Division
*	Multiplication
^	Power

Built in Functions

Pushing the "Insert Function" button in the "Formula" tab will give access to all of the built-in functions. These functions all require an *argument* (i. e., the *arg* in the parentheses), which can be either a list of numbers, cell references separated by commas, or two cell references separated by a

colon to denote a range of cells. Some useful functions are shown in Table 2, and examples of how they are used are given in Table 3.

Table 2. Some useful built-in Excel functions.

AVERAGE(<i>arg</i>)	Calculates the average of a set of numbers.
COS(<i>arg</i>)	Calculates the cosine of the argument (in radians).
COUNT(<i>arg</i>)	Counts the number of cells that contain numbers in the specified range.
COUNTA(<i>arg</i>)	Counts the number of cells that are not empty in the specified range.
EXP(<i>arg</i>)	Raises e to the power of the argument.
LN(<i>arg</i>)	Returns the natural logarithm (ln) of argument.
LOG(<i>arg</i>)	Returns the logarithm of argument; defaults to base 10 if no base is specified
LOG10(<i>arg</i>)	Returns the base 10 logarithm (log) of argument.
PI()	Returns the number pi (this function has no argument).
SQRT(<i>arg</i>)	Takes the square root of the argument.
STDEV(<i>arg</i>)	Calculates the standard deviation of a set of numbers.
SIN(<i>arg</i>)	Calculates the sine of the argument (in radians).
SUM(<i>arg</i>)	Adds all the numbers in a range of cells.
TINV(α , $N-1$)	Calculates the value of t at a given confidence limit and degree of freedom ($N - 1$). Note that $(1 - \alpha)$ is the fractional confidence limit (<i>e. g.</i> , for the 95% confidence limit, $\alpha = 0.05$).

Starting with Excel 2010 Microsoft introduced T.INV.2T to replace TINV. The new function may provide improved accuracy. The TINV function is available in Excel 2010 for backward compatibility, but it may not be available in future versions of Excel.

Table 3. Examples of some Excel functions as they would appear in a spreadsheet cell.

=AVERAGE (A2:A5)	Calculates the average of cells A2, A3, A4 and A5.
=AVERAGE (A2,A4,A5)	Calculates the average of cells A2, A4 and A5.
=SUM(A2:A5)	Adds all the numbers in cells A2, A3, A4, and A5.
=EXP(-1*A2)	Raises e to the -A2 power (note that writing -A2 instead of -1*A2 may not work in older versions of Excel).
=TINV(0.05,COUNT(A2:A10)-1)	Calculates the value of t at the 95% confidence limit for the data contained in cells A2 to A10, inclusive. Note the use of the COUNT function to determine the degrees of freedom.

Selecting Cells

To select a single cell, simply click (note that unless otherwise noted “click” means to momentarily depress the left mouse button) on it. To select more than one cell in the same row or column, click on the first cell, and while still holding the left mouse button down, drag the mouse to the last cell. The selected cells will be highlighted. You can also select a block of cells by clicking on a corner cell and dragging the mouse to the diagonal corner. An alternate way to select cells is to click on the first cell, hold the *Shift* key down and then select the last cell containing data. Once the cells are selected; let up the *Shift* key. If you want to select an entire column or row, click on the column or row header. To select multiple columns or row, simply













click on the header and drag. This last method of selecting cells is very useful if you have a large amount of data that needs to be formatted.

Non-adjacent cells can be selected by holding the *Ctrl* key while clicking on each cell. Highlighting non-adjacent rows or columns is a combination of the steps for selecting nonadjacent cells and rows/columns. Once the first column (or row) has been highlighted, the second column is selected by pressing and holding the *Ctrl* key while you click with the left mouse button in the first cell. While holding the *Ctrl* key and the left mouse button, drag to select the values in this column. The other way to select the second column is to press and hold the *Ctrl* key while clicking in the first cell, then let up the *Ctrl* key and press and hold the *Shift* key. Click on the last cell in the column that contains data and all the cells in between will be highlighted.

Formatting and Toolbar Functions

Many formatting functions are depicted as icons in the “Home” tab along the top of the spreadsheet. The formatting functions are grouped in a category (such as “Font” or “Alignment”). If the function you want is not displayed, there is an expansion button in the bottom right corner of the category which, when clicked, gives more formatting functions. For the “Font,” “Alignment,” and “Number” categories, this expansion button brings opens the “Format Cells” window. (see Fig. 1: *Format Cells* dialog box). For the “Clipboard” category, it opens the clipboard. If your Excel window is not maximized to full-screen, the function categories become buttons instead of displaying any of the functions. To access the functions, click the category’s button. Some of the more commonly encountered formatting and toolbar functions are given in Table 4. You will note that some of these tools are the same as in word processing programs and that some have menu bar buttons or shortcut keys (e. g., to make

Table 4. Common tool bar function buttons.

Button	Function	Shortcut	Menu Bar Location
	Align left	<i>Ctrl-l</i>	<i>Home, Alignment</i>
	Align center	<i>Ctrl-e</i>	<i>Home, Alignment</i>
	Align right	<i>Ctrl-r</i>	<i>Home, Alignment</i>
	Bold	<i>Ctrl-b</i>	<i>Home, Font</i>
	Italic	<i>Ctrl-i</i>	<i>Home, Font</i>
	Underline	<i>Ctrl-u</i>	<i>Home, Font</i>
	Increase number of decimal places		<i>Home, Number</i>
	Decrease number of decimal places		<i>Home, Number</i>
	Undo last action(s)	<i>Ctrl-z</i>	<i>Window Title, Undo</i>
	Redo last action(s)	<i>Ctrl-y</i>	<i>Window Title, Redo</i>
	Sort cells in ascending order		<i>Data, Sort & Filter</i>
	Sort cells in descending order		<i>Data, Sort & Filter</i>

something bold you can either click the bold icon in the use the key sequence *Ctrl-b*). To perform any one of the various formatting options, select the object or cell that you wish to format, then

use the appropriate shortcut or click the desired formatting button in the “Home” tab. Several formatting functions are also available in the menu displayed when a cell is right clicked. It is often more efficient to select multiple objects or cells to format before performing the actual formatting. To access the right-click menu with multiple cells selected, select the desired cells then right click one of the selected cells.

Clicking on any object, and selecting the appropriate format from the respective tab and category, allows you to change the formatting of the object. Right clicking on a selected object and choosing the format option will also bring up the format dialog box. To get the superscript (or a subscript), select the cell you wish to edit. The contents of the cell will appear in the formula bar. Use the mouse to highlight the text to be superscripted (or subscripted), right click, then select “Format Cells...” and click on the superscript (or subscript) check box. Another way to do this same formatting is to highlight the desired text and click the *Format Font* dialog box launcher (see Fig. 1) and click the superscript (or subscript) check box.







Figure 1. “Format Cells” dialog Box (font tab selected)




To put symbols in axis labels or titles, use the same procedure, except change the font to Symbol font. The correspondence between the Greek letters (Symbol font) and the Roman alphabet (Normal font) for selected symbols is given in Table 5. Symbols that do not appear in Table 5 (such as \pm , \bullet , \neq , \leq , \geq , Å, Â, Ä, Æ, Ñ, ..., etc.) can be added by selecting “Symbol” (Ω) from the Insert tab. Clicking on the appropriate symbol and then clicking the Insert button will complete the task. When finished inserting symbols, close the *Symbol* window.

Table 5. The relationship between symbol and Roman fonts.


Symbol	Roman	Symbol	Roman	Symbol	Roman	Symbol	Roman
α	a	ν	n	A	A	N	N
β	b	o	o	B	B	O	O
χ	c	π	p	X	C	Π	P
δ	d	θ	q	Δ	D	Θ	Q
ε	e	ρ	r	E	E	P	R
ϕ	f	σ	s	Φ	F	Σ	S
γ	g	τ	t	Γ	G	T	T
η	h	u	u	H	H	Y	U
ι	i	ϖ	v	I	I	ζ	V
φ	j	ω	w	ϑ	J	Ω	W
κ	k	ξ	x	K	K	Ξ	X
λ	l	ψ	y	Λ	L	Ψ	Y
μ	m	ζ	z	M	M	Z	Z


The number of significant figures displayed can be changed using the increase/decrease decimal place button ( ) in the *Home* tab and *Number* category. One can also select the *Format Font* dialog box launcher (, see Fig. 1) and the *Number* tab to increase/decrease the number of decimal places displayed. The latter method also allows you to change all aspects of how the numbers are displayed. Note that Excel uses the format “0.00E+00” for scientific notation, not “0.00 x 10+00” or “0.00 x 1000”. If you enter a number in the latter formats, Excel will treat it as text.

Sorting

To sort data you must first select the cells to sort. Go to the *Data* tab and click on the *Sort* button () in the *Sort & Filter* category. If you have selected one column and there are data in adjacent cells Excel will give a warning message. If you Continue with the current selection, the selected cells will be sorted and the adjacent cells will **not** be sorted. Do **not** select this option if your data are (x, y) data pairs, because this will ruin the pairs. If you want to sort data pairs, select the cells containing both values and then sort. When several rows or columns have been selected, using the toolbar buttons such as  causes the left-most column or top-most row to be sorted. If you use , you will be prompted for which row/column to sort by.


Graphing

To graph data in the spreadsheet, first select the data. Excel assumes that the *x*-values are in the first column/row selected and that the *y*-values are in the second and any subsequent column/row selected. Now select the *Insert* tab. In the *Charts* category, you will be able to select which type of chart is most appropriate. For most cases in chemistry, select the “Scatter”() drop-down menu and “Scatter with only Markers” as the graph type.¹ The chart will be generated automatically and you will then be able to use the *Chart Tools* tabs. See the [Preparing Graphs](#) section for more information on graphing.¹

The first *Chart Tools* tab is *Design* where you can revise your chart type, manipulate the data selected, and modify the layout of your chart. The *Layout* tab is where you can insert photos/shapes/text, add chart labels (chart title, axis titles, and legend labels), and turn gridlines on or off. The *Format* tab is where you can change colors of data selections and chart components and change the size of your chart. Your charts should be on separate Excel worksheets from your data. To do this, select your chart and go to the *Design* tab and *Location* category. There you will find a “Move Chart” option. With your chart selected, click on “Move Chart” (); a “Move Chart” dialog box will appear. Select *New Sheet* for the chart location.

Once a graph has been created it can be edited by simply right-clicking the chart object (axis, data series, etc.), and choosing *Format* from the drop-down menu.

A trendline can be added to a graph by selecting the *Layout* tab and in the *Analysis* category click the *Trendline* button. In chemistry we usually want a *Linear Trendline*. If there is more than one data set on the graph, you can choose which data set requires a trendline. If you click on *More Trendline Options...* you can check the boxes labeled “Display Equation on chart” and “Display R-squared value on chart”. If you have a trendline but do not have its equation displayed, right-click the trendline and click “Format Trendline”. From the “Format Trendline” dialog box, select the “Display Equation on chart” check box. Generally we will not force the line to go through a particular value on the y-axis, so leave the “Set Intercept” box unchecked. Once you close the dialog box, Excel will perform a linear least squares fit of the selected data. Note that while a trendline will give you the equation of the best-fit line through the data, it is usually insufficient for the data analysis that is often required in chemistry. Please see the *Regression Package* section below for a more powerful means of fitting data.

It is possible to copy a chart from Excel and paste it into another program. Make the chart that you wish to copy be the active chart, and while holding the shift key, click on *Copy* () in the *Home* tab and *Clipboard* category. When in the other program, paste the object by clicking the *Paste* button in the *Home* tab of Word, right clicking and selecting a *Paste* option, or use keyboard shortcut *Ctrl+v*. Once in Word, the image size can be modified to fit the page and further edited by selecting the *Format Picture* tab on the ribbon.

Regression Package

In addition to the trendline, Excel also has a more complete regression package, which gives much more information than the simple trendline treatment. This package is part of the *Analysis ToolPak* add-in that can be activated in Excel. Depending on your computer, you may need to install the *Analysis ToolPak* before using the regression tools described here. The installation varies depending on the version of Excel that you are using. Version-specific instructions are below. **NOTE:** The *Data Analysis* option will only appear if the active window is a worksheet; it will not appear if the active window is a chart.

•If you are using Excel 2010 or Excel 2007, select the *Data* tab. If *Data Analysis* appears as an option in the *Analysis* category, the *Analysis ToolPak* is already installed and the regression package can be accessed by clicking on *Data Analysis* and in the pop-up window selecting *Regression* and clicking *OK*. If *Data Analysis* is not an option on the *Data* tab, go to <http://office.microsoft.com/en-us/excel-help/load-the-analysis-toolpak-HP010342659.aspx?CTT=1> (last visited January 10, 2012) for Excel 2010 or to

<http://office.microsoft.com/en-us/excel/HP100215691033.aspx> (last visited January 10, 2012) for Excel 2007 to install the *Analysis Toolpak*.

•If you are using Excel 2003, select the *Tools* menu. If *Data Analysis* appears as an option, the *Analysis ToolPak* is already installed. To use the regression package, select *Tools*, *Data Analysis* from the menu bar and in the pop-up window select *Regression* and click *OK*. If these statistical tools are not present, go <http://office.microsoft.com/en-us/excel-help/load-the-analysis-toolpak-HP001127724.aspx?CTT=1> (last visited January 10, 2012) to install the *Analysis Toolpak*.

•Mac Users: The most recent version of Excel for the Mac does not have the *Analysis Toolpak* built in. Microsoft suggests the following alternative, which functions nearly identically to the *Analysis ToolPak*. You can use *StatPlus:mac LE* to perform many of the functions that were previously available in the *Analysis ToolPak*, such as regressions, histograms, analysis of variance (ANOVA), and *t*-tests. First, visit the *AnalystSoft* web site (<http://www.analystsoft.com/en/products/statplumacle/>, last visited January 10, 2012), and then follow the instructions on the download page. After you have downloaded and installed *StatPlus:mac LE*, open the workbook that contains the data that you want to analyze. Then open *StatPlus:mac LE*. The functions are located on the *StatPlus:mac LE* menus.





Figure 2. Pop-up window for configuring the regression routine in Excel.

When the regression package is selected, the window shown in Fig. 2 will appear. Input the *Y Range* and *X Range* (if you push the  button in the regression window, the window will minimize and you can choose the desired range by clicking and dragging). Because we do not generally want to force a fit to go through zero, leave the check box for *Constant is Zero* unselected. By clicking on the check box next to *Confidence Level* and entering a new number in the input box, one can select a new confidence limit. For our purposes this is best left at 95%. Leave the “*Output Options*” set for a *New Worksheet Ply*. Press *OK*, and Excel will perform the regression and output the results to a new worksheet, an example of which is shown as Fig. 3. The new worksheet contains a good deal of statistical information, but for a basic analysis we need be

concerned with only a few items. The first is the value of R^2 (R Square), which is a measurement of how related the two measured quantities are ($R^2 = 1$ is a perfect linear relationship). For many of the fits that you will encounter in chemistry $R^2 > 0.999$ are common, and so the use of this parameter to determine how good the fit is or to compare two fits, is somewhat limited. In this case, the standard deviation about the regression, which is given as the *Standard Error* under the *Regression Statistics* heading, is more useful.

The values of the intercept and slope (*X Variable 1*) are under the heading of *Coefficients*, while the estimated standard deviation of each parameter is given under the heading *Standard Error*. The *Lower 95%* and *Upper 95%* entries give the confidence interval for a parameter at 95% confidence. If you selected a confidence level other than 95% the last two columns (the second set of columns in Fig. 3 labeled *Lower 95.0%* and *Upper 95.0%*) will contain the confidence interval at the level. You can calculate the uncertainty (at the $xx\%$ confidence level) by subtracting the entry under *Lower $xx\%$* from the *Upper $xx\%$* and dividing by 2. For the example given in Fig. 3, the uncertainty in the slope is $(2.727447482 - 1.662430005)/2$ which is rounded to one significant figure, giving ± 0.5 as the answer.

Important! If you are performing a propagation of error analysis (see also the [Introduction to Statistics in Chemistry](#) and [Propagation of Uncertainty](#) sections of this document) on a value calculated from a regression analysis, you do not use the uncertainty! Instead, you must use the standard error (i. e., the standard deviation), not the confidence interval. For example, if we used the slope shown in Fig. 3 in a calculation and we wanted to see how its uncertainty affected the calculated value, we would use ± 0.2 (the standard error in the slope rounded to one significant figure) not ± 0.5 ! If an intercept of zero was specified in the regression window, or if an error occurred, the error “#N/A” will appear for some entries.

SUMMARY OUTPUT								
Regression Statistics								
Multiple R	0.999509156							
R Square	0.999018554							
Adjusted R Square	0.998992726							
Standard Error	0.005716936							
Observations	40							
ANOVA								
	df	SS	MS	F	Significance F			
Regression	1	1.264204688	1.264204688	38680.37433	9.01275E-59			
Residual	38	0.001241968	3.26834E-05					
Total	39	1.265446655						
	Coefficients	Standard Error	t Stat	P-value	Lower 95%	Upper 95%	Lower 95.0%	Upper 95.0%
Intercept	-0.329484975	0.002183057	-150.9282695	2.08155E-54	-0.333904342	-0.325065608	-0.333904342	-0.325065608
X Variable 1	-0.00049891	2.53675E-06	-196.6732679	9.01275E-59	-0.000504045	-0.000493775	-0.000504045	-0.000493775

Figure 3. Typical output of the Excel regression package.

Excel's regression package also has a number of advanced features that you may read about, if you wish, in the [Advanced Regression with Microsoft Excel](#) section.


Printing


The appearance of hardcopy can be viewed by going to the File tab and selecting *Print*. On the right hand of the screen now is a preview of what your printed page would look like. To adjust

page settings, such as margins, go to the *Page Layout* tab. You will need to experiment with combinations of page orientation (portrait or landscape), scaling, margins, etc., until you have the page that matches what you want. If possible, a table should fit completely on a single page. When satisfied with the page settings, return to the *Print* menu under the *File* tab, and push the large “Print” button.

Miscellaneous Features

Excel has a “fill” feature that is useful when one wants to construct a set of data where the x values are all equally spaced between some starting and ending value. For example, if we wanted a graph of the function e^{-x} , we could create the x values for the graph using the fill function. To fill a series first type the starting value in the first cell of the series and select it by clicking on that cell.

Next, in the *Home* tab in the *Editing* category, select the *Fill* button (). Select whether you want the created series to occupy *Rows* or *Columns* under the “Series in” heading. Under the heading “Type”, select *Linear*, which is the most commonly encountered way to fill a data series. Finally, select the *Step* value and the *Stop* value and hit *Okay*. Excel will fill in the column or row starting at the value initially by adding the step value to it until the stop value is reached.

Error bars, which show the uncertainty in a given point, may be added to a graph in Excel, by first going to the *Layout* tab, *Analysis* category, and selecting *Error Bars*. There are three default options. For more options, click “More Error Bar Options...” You will then be able to select among various options for how the error bars appear (on both sides, one side, or none) and how large to make them. The *Fixed* value option is good when all of the data have approximately the same uncertainty, while *Percentage* assigns the uncertainty as a certain percent of the y values. This author has never found the *Standard deviation* option to give any useful depiction of the uncertainty, but it may be useful in certain situations. The *Standard error* assigns the uncertainty based on the standard error, which seems to be very similar to the fixed value option. The final option of having *Custom* error bars is very useful when each point has a different uncertainty. You can always simply enter whatever value you want in the boxes, but by pushing the  button you can specify a range of cells in your spreadsheet that contains your uncertainties. For example, if you have a set of volumes in cells B3 through B10 and their estimated uncertainties from a propagation of error analysis in cells C3 through C10, you would only need to enter the cell references for the uncertainties in the boxes, and each point has its own error bar.

References

1. Dodd, J. S., Ed. *The ACS Style Guide: a Manual for Authors and Editors*; the American Chemical Society: Washington, DC, 1986.

Advanced Regression with Microsoft Excel

J. M. McCormick

It is possible to have Excel perform a non-linear least square regression. One simple trick is to create columns each containing the variable of interest to the requisite power. For example, if we wanted to fit a set of data to a third order polynomial (i. e., $y = ax + bx^2 + cx^3 + d$, where a , b , c and d are constants that we need to find), then we would create columns containing the independent variable to the desired powers, as shown in Fig. 1. Note that in addition to x^1 , x^2 and x^3 there must be a column containing x^0 , which contains only ones. Once the spreadsheet is set up as shown below, select Tools, Data Analysis from the menu bar, and scroll down to *Regression*, select it and click OK.

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O
	x^0	x^1	x^2	x^3	y										
1	1	-13	169	-2197	2571.7										
2	1	-12	144	-1728	2598.8										
3	1	-11	121	-1331	2625.5										
4	1	-10	100	-1000	2651.7										
5	1	-9	81	-729	2677.5										
6	1	-8	64	-512	2702.8										
7	1	-7	49	-343	2727.5										
8	1	-6	36	-216	2751.6										
9	1	-5	25	-125	2775.3										
10	1	-4	16	-64	2798.4										
11	1	-3	9	-27	2821										
12	1	-2	4	-8	2843.1										
13	1	-1	1	-1	2864.7										
14	1	1	1	1	2905.9										
15	1	2	4	8	2925.4										
16	1	3	9	27	2944.4										
17	1	4	16	64	2962.7										
18	1	5	25	125	2980.5										
19	1	6	36	216	2997.6										
20	1	7	49	343	3014.1										
21	1	8	64	512	3029.8										
22	1	9	81	729	3044.8										
23	1	10	100	1000	3059.1										
24	1	11	121	1331	3072.6										
25	1	12	144	1728	3085.4										
26	1	13	169	2197	3097.5										
27	1	14	196	2744	3108.8										

Figure 1. Sample spreadsheet that is ready to be fit to the cubic expression $y = ax + bx^2 + cx^3 + d$ using Excel's regression package.

The regression pop up will appear (Fig. 2). In the *Input Y Range* box type in the cell addresses that contains your y values (or click and drag to select them from the worksheet). In Fig. 1 this would be the cells F2:F28 (enter F2:F28 in the *Input Y Range* box). While entering the y values is no different than for a linear regression, inputting the x values is very different in that the *Input X Range* box must contain all of the columns containing a power of x. In Fig. 1 this would be columns B through E (cells B2:E28). Once the x and y ranges are set, you can set any other parameters as desired, and click OK.




Figure 2. Excel's regression set-up window

If you have selected to have the results of the fit put in a new worksheet ply, then you will see something like what is shown in Fig. 3 (note that the columns in this spreadsheet have been modified to clearly show their contents). Each of the rows labeled *X Variable*, except *X Variable 1*, correspond to one of the constants in the polynomial. The row labeled *Intercept* will always be the constant in the polynomial that is not multiplied by x . In this case where the polynomial is $y = ax + bx^2 + cx^3 + d$, the *Intercept* is d , *X Variable 2* is a , *X Variable 3* is b and *X Variable 4* is c . *X Variable 1* will always be meaningless. The values under the *Standard Error* and *Upper* and *Lower 95%* for each variable have the same meaning as they do in a simple regression.

SUMMARY OUTPUT									
Regression Statistics									
Multiple R	0.99999992								
R Square	0.99999984								
Adjusted R Square	0.956521558								
Standard Error	0.071913988								
Observations	27								
ANOVA									
	df	SS	MS	F	Significance F				
Regression	4	741557.0877	185389.2719	47796553.98	5.62092E-76				
Residual	23	0.118947299	0.005171622						
Total	27	741557.2067							
Coefficients									
Intercept	2885.475246	0.021356602	135109.2833	7.4211E-104	2885.431066	2885.519425	2885.431066	2885.519425	
X Variable 1	0	0	65535	#NUM!	0	0	0	0	
X Variable 2	20.57469078	0.004209556	4887.615914	1.06412E-70	20.56598265	20.58339891	20.56598265	20.58339891	
X Variable 3	-0.300961361	0.000244281	-1232.028513	6.19266E-57	-0.301466696	-0.300456027	-0.301466696	-0.300456027	
X Variable 4	-0.002082423	3.30678E-05	-62.97437398	2.93583E-27	-0.002150829	-0.002014017	-0.002150829	-0.002014017	

Figure 3. Excel worksheet showing the results for the fit of the data in Fig. 1 to a cubic function. Note that the values of R and R^2 are not going to be particularly useful in defining how well this expression fit the data, or in comparing this fit to a fit using another expression. For this one must use the standard error about the regression (listed as the *Standard Error* under the *Regression*

Statistics heading). Also note that we could set up a column of calculated y values based on this fit (for example, in column G of the spreadsheet shown in Fig. 1) by entering the function in the appropriate column and using absolute references to the values in the output spreadsheet (Fig. 3).




Figure 4. Excel spreadsheet set up to fit a data set to any expression in the variable x using the *Solver* tool.

Another way to have Excel perform a non-linear least squares fit is to use the *Solver* tool. Before getting started, you will need to set up the spreadsheet, but in a different way than for the *Regression* package. Again let us assume that we want to fit a set of y values to a cubic expression in x (i. e., $y = ax + bx^2 + cx^3 + d$). We only need a single column for the x values, but we will need to have two columns for the y values (see Fig. 4). The first column of y values (column C in Fig. 4) will contain the y values that we actually measured. It is labeled y_{obs} (for observed values) to avoid confusion with the calculated y values (column D in Fig. 4), which we label y_{calc} (for calculated values). The y_{calc} are determined using the values of a , b , c and d . Thus, to find the values of y_{calc} we will have to have cells containing trial values of a , b , c and d (these trial values are sometimes referred to as seed values). These can be placed in any convenient place in the spreadsheet. In Fig. 4 they are in cells J2 through J5 (note that to minimize confusion the cells I2 through I5 contain labels for the adjacent cells). We could type in any number that we want as trial values, but to maximize the chance of quickly obtaining a successful fit, the trial values should be reasonable guesses. The last column that we will need is one that contains the square of each y_{obs} minus the corresponding y_{calc} (hence the name least squares). In the current example, this is column E in Fig. 4. The final step is to create a cell containing the sum of the values in column E, which is cell H6 in Fig. 4.

Once the spreadsheet is ready, select *Tools, Solver* from the menu bar. The *Solver* window will pop up (Fig. 5). Under *Set Target Cell*, enter the cell containing the value of $S(y_{obs} - y_{calc})^2$, which in our example is cell H6. Select the *Min* (minimum) radio button under the *Equal To* heading

(note that if you select *Value Of* and enter zero, you may not get a reasonable answer). In the *By Changing Cells* box, enter the cell(s) containing the trial values that you want to change. Since we want to fit the data by adjusting all of the values in cells H2 to H5, we would enter H2:H5 in this box (or by highlighting all of the cells in the spreadsheet). Note that you can choose to vary all or only some of the trial values by which cells you select at this point. Sometimes it is useful to vary only a subset of all the parameters, while holding others constant, but this is highly dependent on the situation. You may also apply additional constraints to the fit in the *Subject to Constraints* box by clicking the *Add* button. For example, to constrain the value in cell H2 to be positive, one would click *Add* and then in the new pop-up enter H2 in the *Cell Reference* box, select the "> =" in the middle pull-down box, and then enter "0" in the *Constraint* box.




Figure 5. Solver set-up window

Once everything is set, click *Solve*. *Solver* will then attempt change the cells that you told it to change, subject to any constraints that you specified, to minimize the difference between the observed and calculated y values. If it is successful it will display a window stating that it found a solution and asking whether you want to keep it (note that the best fit values of your variables and the function that you were minimizing will be displayed in the spreadsheet behind the pop-up). If you want to keep the fit, click *OK*, if not click *CANCEL*. If a satisfactory fit was not found, a similar message will be displayed and you can either keep or discard the result.

Sometimes you will need to restart the *Solver* tool (that is reinitiate the procedure with the values from a previous fit) before a satisfactory fit is obtained. In any case, it is good practice to perform the fit with different values of the initial parameters to be sure that the procedure converges to a unique solution. It may also be necessary to adjust the fitting procedure itself (done by clicking the *Options* button in the main *Solver* window) to get a good result. However, if this occurs it may be advisable to find a more robust fitting algorithm or to find a simpler model.

Unlike the *Regression* package, *Solver* does not give uncertainties in the fit parameters. However, there is a way to estimate the uncertainties,¹ but it is somewhat time consuming.

References

1. Harris, D. C. *J. Chem. Educ.* **1998**, 75, 119-121. This article may be viewed as a PDF file at <http://pubs.acs.org/doi/pdf/10.1021/ed075p119> for Truman addresses and *J. Chem. Educ.* subscribers only.

Propagation of Uncertainty

J. M. McCormick

Introduction

Every measurement that we make in the laboratory has some degree of uncertainty associated with it simply because no measuring device is perfect. If a desired quantity can be found directly from a single measurement, then the uncertainty in the quantity is completely determined by the precision of the measurement. It is not so simple, however, when a quantity must be calculated from two or more measurements, each with their own uncertainty. In this case the precision of the final result depends on the uncertainties in each of the measurements that went into calculating it. In other words, uncertainty is always present and a measurement's uncertainty is always carried through all calculations that use it.

Fundamental Equations

One might think that all we need to do is perform the calculation at the extreme of each variable's confidence interval, and the result reflecting the uncertainty in the calculated quantity. Although this works in some instances, it usually fails, because we need to account for the distribution of possible values in all of the measured variables and how that affects the distribution of values in the calculated quantity. Although this seems like a daunting task, the problem is solvable, and it has been solved, but the proof will not be given here. The result is a general equation for the propagation of uncertainty that is given as Eqn. 1.¹ In Eqn. 1 f is a function in several variables, x_i , each with their own uncertainty, Δx_i .

$$\Delta f = \pm \sqrt{\sum_{i=1}^n \left(\frac{\partial f}{\partial x_i} \right)^2 (\Delta x_i)^2} \quad (1)$$

From Eqn. 1, it is possible to calculate the uncertainty in the function, Δf , if we know the uncertainties in each variable and the functional form of f (so we can calculate the partial derivatives with respect to each variable). It is easier to understand how this all works by doing several examples.

Example 1: $f = x + y$ (the result is the same for $f = x - y$).

Let the uncertainty in x and y be Δx and Δy , respectively. Taking the partial derivatives with respect to each variable gives: $\frac{\partial f}{\partial x} = 1$ and $\frac{\partial f}{\partial y} = 1$. The uncertainty in f is then

$$\Delta f = \pm \sqrt{(1)^2(\Delta x)^2 + (1)^2(\Delta y)^2} = \pm \sqrt{(\Delta x)^2 + (\Delta y)^2} \quad (2)$$

Example 2: $f = x \cdot y$ (also works for $f = x/y$)

Again let the uncertainty in x and y again be Δx and Δy , respectively. Taking the partial derivatives with respect to each variable gives: $\frac{\partial f}{\partial x} = y$ and $\frac{\partial f}{\partial y} = x$. The uncertainty in f is then $\Delta f = \pm \sqrt{(y)^2(\Delta x)^2 + (x)^2(\Delta y)^2}$.

This result is more commonly written by dividing both sides by $f = x \cdot y$ to give

$$\frac{\Delta f}{f} = \pm \sqrt{\left(\frac{\Delta x}{x}\right)^2 + \left(\frac{\Delta y}{y}\right)^2} \quad (3)$$

Although the idea of error propagation may seem intimidating, you have already been using it since your first chemistry class when you applied the rules for significant figures in calculations. These rules are simplified versions of Eqn. 2 and Eqn. 3, assuming that Δx and Δy are both 1 in the last decimal place quoted. The formal mathematical proof of this is well beyond this short introduction, but two examples may convince you.

If we add 15.11 and 0.021, the answer is 15.13 according to the rules of significant figures. This assumed that $\Delta x = 0.01$ ($x = 15.11$) and $\Delta y = 0.001$ ($y = 0.021$), substituting these values into Eqn. 2, we get $\Delta f = \pm \sqrt{(0.01)^2 + (0.001)^2} = \pm 0.010050$. Remembering our basic statistics (see [Introduction to Statistics in Chemistry](#)), we know that the uncertainty begins in the first non-zero decimal place, which in this case this means that the last significant figure in the sum is the 1/100th place. According to the rules for propagation of error the result of our calculation is 15.13 \pm 0.01, exactly what the significant figure rules gave us.

If we had multiplied the numbers together, instead of adding them, our result would have been 0.32 according to the rules of significant figures. Again assuming $\Delta x = 0.01$ and $\Delta y = 0.001$, and using Eqn. 3, we can determine Δf as follows.

$$\Delta f = \pm f \sqrt{\left(\frac{\Delta x}{x}\right)^2 + \left(\frac{\Delta y}{y}\right)^2} = \pm (0.31731) \sqrt{\left(\frac{0.01}{15.11}\right)^2 + \left(\frac{0.001}{0.021}\right)^2} = \pm 0.015111$$

Once again we see that the uncertainty begins in the second decimal place, which gives the same result as the significant figures result gave.

The significant figure rules are important to know and use in all chemistry calculations, but they are limited in that they assume an uncertainty in the measured quantities. So while the significant figure rules are always to be used in any calculation, when precision matters a propagation of error analysis must also be performed to obtain an accurate prediction of the uncertainty arising from the precision of the measured quantities.

Worked Examples

Problem 1

In CHEM 130, you have measured the dimensions of a copper block (assumed to be a regular rectangular box) and calculated the box's volume from the dimensions. In that exercise you were given an equation that allowed you to calculate the minimum uncertainty that could be expected in the box's volume based solely on the uncertainties in the measured dimensions, now derive that equation using the procedure given above.

Solution

Let x , y and z be the box's length, width and height, respectively, and the uncertainties be Δx , Δy , Δz . Since $V = x \cdot y \cdot z$, we can use Eqn. 1 to determine the uncertainty in the volume (ΔV), which results in Eqn. 4. We know that $\frac{\partial V}{\partial x} = yz$, $\frac{\partial V}{\partial y} = xz$ and $\frac{\partial V}{\partial z} = xy$, and can then make these substitutions in Eqn. 4 to give Eqn. 5.

$$\Delta V = \pm \sqrt{\left(\frac{\partial V}{\partial x}\right)^2 (\Delta x)^2 + \left(\frac{\partial V}{\partial y}\right)^2 (\Delta y)^2 + \left(\frac{\partial V}{\partial z}\right)^2 (\Delta z)^2} \quad (4)$$

$$\Delta V = \pm \sqrt{(yz)^2 (\Delta x)^2 + (xz)^2 (\Delta y)^2 + (xy)^2 (\Delta z)^2} \quad (5)$$

Dividing both sides by V gives Eqn. 6 and simplifying gives Eqn. 7 (which you probably could have guessed from the form of Eqn. 1 and Eqn. 3). Multiplying both sides by V then gives the equation used in the CHEM 130 *Determination of Density* exercise.

$$\frac{\Delta V}{V} = \pm \sqrt{\frac{(yz)^2 (\Delta x)^2 + (xz)^2 (\Delta y)^2 + (xy)^2 (\Delta z)^2}{(xyz)^2}} \quad (6)$$

$$\frac{\Delta V}{V} = \pm \sqrt{\left(\frac{\Delta x}{x}\right)^2 + \left(\frac{\Delta y}{y}\right)^2 + \left(\frac{\Delta z}{z}\right)^2} \quad (7)$$

Note that there are several implications of Eqn. 7. First, if one side has a large uncertainty relative to the length of that side (such as when one side is very short), then this side will dominate the uncertainty. Second, when the volume is large and the uncertainty in measuring a dimension is small compared to the uncertainty in the measurement, then the uncertainty in the volume will be small. The experimental implication of this is that, if you want the smallest uncertainty in a box's volume, make sure it is a big box, with no unusually short side and use the most precise measurement tool possible.

Problem 2

You have measured the volume and mass of a set of regular wooden blocks and have fit a graph of their volume as a function of their mass to a straight line using the regression package in Excel. What is the predicted uncertainty in the density of the wood (Δd) given the uncertainty in the

slope, s , of the best fit line is Δs and the uncertainty in the intercept is Δb ? Note that you have also seen this equation before in the

[Determination](#) of Density exercise, but now you can derive it.

Solution

The relationship between volume and mass is $V = \frac{1}{d}m$. This is a linear equation ($y = s \cdot x + b$) where $= \frac{1}{s}$. Note that b does not affect the value of d and so Δb has no effect on Δd . The relationship between Δs and Δd can be calculated by simply substituting d in place of f and s in place of x in Eqn. 3 to give $\Delta d = \pm d \sqrt{\left(\frac{\Delta s}{s}\right)^2} = \pm d \left(\frac{\Delta s}{s}\right)$.

We could have also have used Eqn. 1. First we need to find the first derivative of the density with respect to the slope, which is $\left(\frac{\partial d}{\partial s}\right) = -\frac{1}{s^2}$. Substituting this into Eqn. 1 gives

$\Delta d = \pm \sqrt{\left(-\frac{1}{s^2}\right)^2 (\Delta s)^2}$, which rearranges to $\Delta d = \pm \sqrt{\left(\frac{\Delta s}{s}\right)^2 \left(\frac{1}{s^2}\right)}$. Recognizing the relationship between s and d , this simplifies to $\Delta d = \pm d \left(\frac{\Delta s}{s}\right)$.

This problem is the simplest example of how one determines the uncertainty in a quantity extracted from a best-fit line. In general you will have the uncertainty in the slope and intercept and the relationship between each of these to the desired quantities. It is then a simple process to apply Eqn. 1, where f is either the slope or intercept.

Propagation of Uncertainty through a Calibration Curve

A situation that is often encountered in chemistry is the use of a calibration curve to determine a value of some quantity from another, measured quantity, for example, in the [Preparation and Analysis of Alum](#) exercise you used a calibration curve in the determination of the percent by mass of aluminum in alum. In that exercise, we did not propagate the uncertainty associated with the absorbance measurement through the calibration curve to the percent by mass. However, in most quantitative measurements, it is necessary to propagate the uncertainty in a measured value through a calibration curve to the final value being sought. The general procedure is quite straight-forward, and is covered in detail in [CHEM 222](#). Therefore, only a very basic review of the fundamental equations and how to implement them in Excel will be presented here. You are referred to any analytical chemistry textbook for more details.²

For a linear least squares analysis we need to define several parameters. We will assume that the equation of a straight line takes the form $y = mx + b$ (where m is the slope and b the intercept) and that the x values are known precisely. Let there be N individual data points (so there are N ordered pairs x_i, y_i) in the calibration curve. Further, let y_{meas} be the average response of our unknown sample based on M replicate measurements, and let S_{meas} be the standard deviation of the result from the calibration curve. Note that S_{meas} is the standard deviation associated with the x value (x_{meas}) corresponding to y_{meas} , and should not be confused with S_r , the standard deviation about the regression. We can then draw up the following table to summarize the equations that we need to calculate the parameters that we are most interested in (x_{meas} and S_{meas}).

Table 1. Relationships between standard equations encountered in a linear least squares analysis and the Excel regression package output and Excel commands. Note that *arg* in the Excel command refers to a range of cells over which the command is to be calculated (e. g., E5:E10).

Equation	Location on Regression Output Worksheet	Excel Command
$S_{xx} = \sum (x_i - x_{avg})^2$	-----	DEVSQ(<i>arg</i>)
$S_{yy} = \sum (y_i - y_{avg})^2$	Under the <i>ANOVA</i> heading it is the entry in the row labeled <i>Total</i> in the <i>SS</i> column.	DEVSQ(<i>arg</i>)
$S_{xy} = \sum (x_i - x_{avg})(y_i - y_{avg})$	-----	-----
$x_{avg} = \frac{\sum x_i}{N}$	-----	AVERAGE(<i>arg</i>)
$y_{avg} = \frac{\sum y_i}{N}$	-----	AVERAGE(<i>arg</i>)
$m = \frac{S_{xy}}{S_{xx}}$	Coefficient listed under <i>X Variable 1</i> .	SLOPE(known y's, known x's)
$b = y_{avg} - mx_{avg}$	Coefficient listed under <i>Intercept</i> .	INTERCEPT(known y's, known x's)
$S_r = \sqrt{\frac{S_{yy} - m^2 S_{xx}}{N - 2}}$	<i>Standard Error</i> under the <i>Regression Statistics</i> heading.	STEYX(known y's, known x's)
$x_{meas} = \frac{y_{meas} - b}{m}$	-----	-----
$S_{meas} = \left(\frac{S_r}{m}\right) \sqrt{\frac{1}{M} + \frac{1}{N} + \frac{(y_{meas} - y_{avg})^2}{m^2}}$	-----	See below.

Although one could enter formulas in various cells to calculate all of the intermediate parameters needed to determine S_{meas} , it is not necessary. One only needs to have a cell in which to enter the number of replicate measurements on the unknown (M) and then it is possible to calculate S_{meas} using only the STEYX, SLOPE, INTERCEPT, COUNT, DEVSQ and SQRT Excel functions. For example, in the spreadsheet shown in Fig. 1, cell D16 contains the formula

“=(STEYX(D3:D13,C3:C13)/SLOPE(D3:D13,C3:C13))*SQRT((1/D15)+(1/COUNT(D3:D13))+((D18-AVERAGE(D2:D13))^2/(SLOPE(D3:D13,C3:C13)^2*DEVSQ(C2:C13))))”

which calculates S_{meas} directly from the potential as a function of temperature data. Adding a cell that will contain y_{meas} (cell D17 in Fig. 1), allows calculation of x_{meas} value (cell D18) and its uncertainty at 95% confidence (cell D19). To review how this is done using S_{meas} and Student's t ,

please see the [Introduction to Statistics in Chemistry](#) section). Note that instead of using N in the calculation of the uncertainty from S_{meas} , one must use $N - 2$ because two degrees of freedom have been used to find the slope and the intercept.




Figure 1. An example of an Excel spreadsheet that may be used to calculate an x value (temperature, in this case) from a measured y value (potential) along with the uncertainty in the measured x value at 95% confidence.

Propagation of Uncertainty of Two Lines to their Intersection

Sometimes it is necessary to determine the uncertainty in the intersection of two lines. This problem is not trivial and the reader is referred to the literature for more details.³

References

1. Andraos, J. J. *Chem. Educ.* **1996**, 73, 150-154. This article is available as a PDF file at <http://pubs.acs.org/doi/pdf/10.1021/ed073p150> for Truman addresses and *J. Chem. Educ.* subscribers.
2. Skoog, D. A.; West, D. M. and Holler, F. J. *Fundamentals of Analytical Chemistry*, 5th Ed.; Saunders College Publishing: New York, 1988; p. 39-42.
3. Carter Jr., K. N.; Scott, D. M.; Salmon, J. K. and Zarcone, G. S. *Anal. Chem.* **1991**, 63, 1270-1270. This article is available as a PDF file at <http://pubs.acs.org/doi/pdf/10.1021/ac00013a017> for Truman addresses and *Anal. Chem.* subscribers.

The Laboratory Report

M. C. Nagan and J. M. McCormick

Introduction

The research paper is the primary means of communication in science. The research paper presents the results of the experiment and interpretation of the data, describes the rationale and design of the experiment, provides a context for the results in terms of previous findings and assesses the overall success of the experiment(s). Scientists working in industrial laboratories do not write as many journal articles as their colleagues in academia, but they routinely write progress reports, which take the same form as a journal article. So no matter what your career goals are, it is important that you become familiar with this style of writing.

There are set rules for preparing a journal article (or a laboratory report). The style requirements vary only slightly from journal to journal, but there are far more similarities than differences in the scientific writing style. If you are writing an article for publication in a particular journal (or preparing a laboratory report in the style of a particular journal) you should consult the *Instructions to Authors* section of the journal's website (this information is also included in the journal's first issue of each year).

There are several style guides^{1,2} and articles³ to help scientists and students prepare their manuscripts. The most useful of these to chemists is the American Chemical Society's (ACS) *ACS Style Guide*,¹ which may be found in the Truman library or may be purchased from the ACS web site. Because of the variation in journal styles, and the requirements for a specific course, your instructor will inform you of specific style requirements for his or her class. This guide is based on the Journal of the American Chemical Society style,⁴ and is meant to provide a good starting point for writing a laboratory report. It is not meant to be the definitive style guide; you must adjust your style to your audience and the journal in which your results will be published.

General Editorial Issues

Although we shouldn't, all of us are swayed by first impressions. How your paper appears to the journal editor or reviewer is their first impression of your science, and it will color their impression of your results, if you let it. Nothing is worse than a sloppily prepared paper with no page numbers, a font that can't be read or which is full of grammatical errors. Remember that everyone will assume that if you did not take the time to write your paper carefully, you did not take the time to do your science carefully.

The following are some general editorial guidelines to follow that will leave a good first impression with your readers.

- 1) Double-space your paper throughout (including figure captions and tables, too).
- 2) Use a reasonably sized font such as 12 point Times.
- 3) For figures, you may choose to use a sans-serif font for better graphics quality such as Arial or Helvetica.
- 4) Use at least 1" margins on all sides.

5) Number the pages. Place the page numbers in the top, right-hand corner, or centered on the bottom of the page. Either style is acceptable and whichever one you choose remain consistent in your numbering scheme throughout the paper.

6) Do not start sentences with symbols or numbers; rather, spell out the full name of the symbol if it is used at the beginning of the sentence. For example, write “Alpha-lactalbumin” instead of “ α -Lactalbumin” when beginning a sentence. Also spell out symbols or numbers in a title, except when part of a chemical name (e.g. 2-hexanol).

7) Spell check the document thoroughly. Have someone, who will give you an honest and complete critique of your paper, read the paper. Revise, revise, revise!

General Stylistic Issues

Uniformity of style is the key to scientific communication. The journal editors, the referees who review a manuscript and the journal readers who are interested in the results presented in a paper all expect certain things to be present in a manuscript and that they are in a certain order. Just like the sloppy-looking paper, a paper that does not adhere to the expected style reflects poorly on the author, no matter how good the science is.

1) The paper should be written in a third person, passive voice. Occasionally, but rarely, it is appropriate to use “we” when describing the intention of the authors. It generally depends upon the intended subject of the sentence. Consider the two sentences below:

- a) Calcium solid (5 g) was poured into a beaker.
- b) We poured calcium solid (5 g) into a beaker.

In the first sentence (a), which is passive, the subject is the calcium solid. In the second sentence, the subject is the experimenters. In scientific articles, the subject is most often the science and not the experimenters.

2) Use the past tense in general (e. g., what was or has been done). However, use present tense when describing properties of molecules or organisms because they still have these properties.

3) Unless directed otherwise, assume the reader of your laboratory report is your peer, the average chemistry student, not the chemistry professor. Therefore, everything should be explained as if the reader knows some chemistry, but is not an expert in the subject of the paper. By no means does the reader know what you are doing, or why you are conducting your experiment. Think about what you would want to know about the subject if you were the reader.

4) Avoid repetition in language. Try not to start each sentence with the same construction and words.

5) Do not use quotes. Unlike humanities or literature papers, quotations are rarely found in scientific articles. However, it is appropriate to paraphrase other authors.

6) Explain technical terms.

Example

“Hemoglobin has a Hill constant, a value that describes the degree of cooperative ligand binding, of 2.8.”

7) Define abbreviations.

Example

“The official colors of Truman State University (TSU) are purple and white.”

8) Place a space between a number and a unit.

Example

“Sephadex (10 g) was combined with deionized H₂O (100 mL) at 25 °C.”

9) Do not start a sentence with a number or “Figure 1” or “Table 1”, etc..

Correct: Milk samples (50 µL) were analyzed by high performance liquid chromatography under three different buffer conditions (Figure 1).

Incorrect: Figure 1 shows the high performance liquid chromatography chromatograms for the sample run under three different buffer conditions.

Incorrect: 50 µL of milk was analyzed by high performance liquid chromatography using three different buffer conditions.

10) There are three ways to refer to a paper in the text.

For example, the citation of the work authored by Jackson, A. K.; Wilson, R. S.; Houk, K. L.*, could appear in the text in any of the following ways. (Note that *et al.* is an abbreviation for *et alia* and that it is italicized because it is not English.⁵)

- a) Jackson et al.
- b) Jackson and coworkers
- c) Houk and coworkers

In the last example we assumed that the author whose name is starred is the principle investigator on the project, and gave them more credit for the work. Note that it is an American convention to list the principle investigators last, while many European and Japanese journals place them first.

Often there are two principle investigators, and in this case both should be mentioned. For example, the work by Jackson, A. K.; Wilson, R. S.*; Houk, K. L.* should be referred to, in the format given in example (c) above, as “Wilson, Houk and coworkers”. If there are more than two principle investigators, it is best to use either of the formats given in example (a) or (b), or to use some other wording to avoid this construction entirely.

Organization/Components

Sections

Sections should appear in your paper in the order described below. All sections but the title have the section explicitly labeled, usually in bold letters to differentiate it from the rest of the text, and left aligned on the page. A blank line should appear after the last word of the section to separate the various sections, but a line should not be placed after the section title.

- 1) Title/Title page
- 2) Abstract
- 3) Introduction
- 4) Experimental (Materials and Methods in some journals)
- 5) Results
- 6) Discussion
- 7) Conclusions
- 8) Acknowledgements
- 9) References
- 10) Tables
- 11) Schemes
- 12) Figure Legends
- 13) Figures
- 14) Supporting Information

Please note that you should not physically assemble your paper in this order. Instead, it is suggested that you compose: a) Materials and Methods, b) Figures, Figure Legends and Tables, c) Results, d) Discussion, e) Conclusions, f) Introduction and Schemes, g) Abstract, and h) Title. Then put all the sections together in the final paper in the order outlined above.

A template is available to help you organize your report at http://chemlab.truman.edu/LabReports_files/FormalTemplate.doc.

Subsections

It may be helpful to organize sections further into subsections. These subsections should have their own titles that are italicized and followed by a period.

Description of Paper Components

Title/Title Page

A title reflects the emphasis and contents of the paper. It tells the reader the paper's topic and it also entices the reader to continue reading further. Therefore, it is not uncommon for the title to reveal the results or major conclusions of the experiment. Examples are given below. The title should be on its own page (the title page), left-aligned at the top of the page, in bold letters. Note that in some journals the title's font size is 2 points larger than the text (i. e., 14-point, if the rest of the paper is in a standard 12-point font). However, this is not standardized and you should check with your instructor for which format he/she wants you to follow.

The title must be brief (2 lines maximum) and grammatically correct. Under the title, write your name and your professional address in italics (*Department of Chemistry, Truman State University, 100 East Normal, Kirksville, MO 63501*).

Example Titles

- 1) Determination of the Differential Fluidity of Water and Benzene by Viscosity Measurements
- 2) Purification of *Alpha*-Lactalbumin from Bovine Skim Milk by Immobilized Metal Ion Affinity Chromatography
- 3) Synthesis and Characterization of Potassium Tris(oxalato)ferrate(III)
- 4) Ionic Composition of Drinking Water Influenced by Pipe Materials: An Atomic Absorption Spectroscopic Analysis

Abstract

The abstract is a one-paragraph summary of the paper that is written in the present tense. As the abstract is the only part of the paper that is entered into article databases, it should be able to stand alone, separate from the paper. The first one to three sentences of the abstract should briefly introduce the reader to the problem studied. Next, the scientific approach, major results and primary significance of the findings should be presented. The abstract is generally 150-200 words (less for shorter papers). This section is normally written after the body of the paper. Because the abstract is separate from the paper, all abbreviations should be written out, or defined, and any references should be written out in full. An example of how a reference might appear in an abstract is

Inhaled fumes from permanent markers have been shown to cause brain damage (Johnson, A. J. *Permanent markers and the brain. J. Am. Brain. Res.* **2004**, 18, 215–218).

Note that in some journals that inclusion of the title in a reference is not required, or accepted practice (*vide infra*).

Introduction

The introduction should present the scientific problem at hand to the reader. Explain to the reader why the experiment was conducted, how it was designed and perhaps, if appropriate, what was found. Literature that is relevant should be incorporated and will help the reader understand the context of your study. A good rule of thumb is to start at the most general topic and progressively move towards the specific.

Here is a general outline for an introduction:

- I. Broad significance of the topic to the chemistry discipline and society in general
- II. Introduction to the topic within chemistry
- III. Description of the specific problem
- IV. General goals and significance of the experiment or research topic

In this section, consider including figures, schemes and equations that complement the text.

While this is similar to the information that you should have written your notebook, the introduction to a paper is different than the background that you included for an experiment (or experiments) in your notebook. Remember that you are trying to reach a larger, more general audience with your paper, and the introduction must be structured to draw the reader in and help them focus on your important results.

Experimental

The *Experimental* section of your paper should be a logical, coherent recount of the experiment(s) conducted. This section should be complete enough for a trained scientist to pick up your report and replicate your experiment. The experimental section in a laboratory report is more concise than the corresponding section in the laboratory notebook. It should not be a step-by-step procedure of the activities carried out during the laboratory period.

The first paragraph of the *Experimental* section contains information on key chemicals used in the procedure. When the chemicals are used as received, there will usually be a statement to that effect and further details are not usually necessary. You will list the chemical supplier's name and the substance's purity will be noted in cases where the chemical is hard to find, it is of a special purity or if there is only one supplier. Do not list lot numbers. If a starting material was synthesized according to a literature procedure, then state this in the opening paragraph and reference the procedure. If purification or drying of the compounds is required, it is described here, also.

The first paragraph often will also list the instruments used to characterize the newly synthesized substances. All instruments and equipment should be specified including the model number of the instrument and the name of the manufacturer (serial numbers are not included). When a spectroscopic or physical method is the focus of the report, it will be described in its own subsection. You are not required to write the experimental in this fashion.

For common techniques, laboratory textbooks or material provided by your instructor should be referenced. However, if a previously published procedure was modified, then this is stated and only the modifications performed are included. If the procedure is your own, then outline the procedure with the main points, including details that are critical to replicating the experiment. These might include the type and size of your HPLC column, the buffer or the concentrations of chemicals.

When the syntheses of substances are reported, the synthetic procedure used to make each substance is described in its own separate paragraph. The paragraph begins with the name of substance, or its abbreviation (if the abbreviation was defined earlier in the paper), in bold face. If numbers are assigned to the compounds, these are also included (in parentheses). Often the synthesis will be written out, even when a literature procedure was followed. The mass and percent yields must be reported. Some of the new compound's characteristics are included at the end of the paragraph describing its synthesis. These include: melting point range (and literature value, if known), elemental analysis (both calculated and found), selected peaks from the mass spectrum (with assignments), selected IR peaks (also with assignments), and any NMR peaks with their chemical shift, multiplicity and integration (you will often find the observed coupling quoted and the assignment of the peaks). The following is an example of how to report a compound's synthesis.

Tris-(2-pyridylmethyl)amine: To a stirred solution containing 10.11 g 2-pyridylmethyl chloride hydrochloride and 3.20 ml 2-pyridylmethyl amine in 20 ml H₂O was added in a slow, drop-wise manner (~1 drop every 25 sec) a solution containing 5.03 g NaOH in 12 ml H₂O so that all of the solution was added in about 1.5 hr. Upon complete addition of the NaOH, the reaction mixture was heated on a heating mantle to 70 °C for 20 min. The cooled reaction mixture was then extracted four times with 50 ml CH₂Cl₂. The combined extracts were dried over Na₂SO₄ and the CH₂Cl₂ was removed using a rotary evaporator. The resulting red oil solidified upon standing. The red solid was then dissolved in a minimum of hot hexane. The yellow solution was decanted from a red oil which did not dissolve and filtered hot. Upon cooling the product crystallizes in large needles, which were recovered by filtration and air-dried. Recrystallization from hexane gave 2.08 g of the product (23% yield). The melting point of product is 85 °C, sharp (literature 87 – 89 °C).¹⁰ ¹H NMR (CDCl₃, ppm): 3.89 (s, 6 H, methylene), 7.14 (m, J = 1.3, 6.1 Hz, 3 H, pyridyl), 7.58 (d, J = 7.8 Hz, 3 H, pyridyl), 7.63 (m, J = 1.8, 7.6 Hz, 3 H, pyridyl), 8.15 (m, J = 0.9, 4.9 Hz, 6 H, pyridyl). ¹³C NMR (CDCl₃, ppm): 60.13 (methylene), 122.01, 122.97, 136.48, 149.06, 159.25 (pyridyl).

The *Experimental* section has two quirky wrinkles on the general scientific style. These are:

- 1) when citing previously published procedures, authors' names are generally not included,

Correct "Purification of the bovine brain isolate was performed according to previously published procedures."^{ref}

Incorrect "The previously published procedure of Jackson *et al.*^{ref} was followed with modifications outlined below."

- 2) when citing the use of a kit, pre-packaged-assay or other commercial equipment with directions, include just the company's name in parentheses; it should not be a full reference.

Example

"The Bradford assay (Sigma) was carried out to determine the total protein concentration of the five protein isolates.

Results

In the *Results* section, the results are presented and summarized in a reader-friendly form. Raw data are not presented here. For instance, it is appropriate to include the average calculated concentration of a solution but not the original absorbance values that were collected from the spectrophotometer; that information is best left in your laboratory notebook.

Graphs and tables often make the data more interpretable and understandable (see the [Preparing Graphs](#) section to review graph preparation). A graph is presented in the paper as a figure. In general, a graph or table is an appropriate representation of the data when more than 2 or 3 numbers are presented. Data that are presented in the form of a graph or table should be referred to but should not be repeated verbatim in the text as this defeats the purpose of a graph. More information on figures and tables is presented later.

The *Results* section also reports comparable literature values for the properties obtained and/or calculated in the paper. Observation of trends in the numerical data is acceptable. However, interpretation of the trend should be saved for the *Discussion* section.

Remember; do not simply report your numerical results. The *Results* section must have a narrative that describes your results. This narrative can include a description of the data (such as spectra or data in graphs), what problems were encountered during data acquisition (and how they were resolved, or not) and a general description of how the raw data were processed to give the final results (not a step-by-step description of everything you did). The reader wants to know what you did, how you did it, what problems you encountered and finally what your results were. Each of these topics must be addressed in the *Results* section in a way that is clear, yet concise.

Discussion

This is the section where the results are interpreted. This section of the paper is analogous to a debate. You need to present your data, convince the reader of your data's reliability and present evidence for your convictions. First, evaluate your data. Do you have good, mediocre, terrible, or un-interpretable data? Evaluate your results by comparing to literature values or other precedents. Explain what results should have been obtained and whether you obtained these expected values. Note that even if expected results were not obtained, you did not fail. Unexpected results are often the most interesting. Perhaps your hypothesis was not correct. Why is this? What new hypothesis do your data suggest? If you feel that your results are not reliable, you need to explain why. Use statistical analysis or chemical principles to support your claims. Was there a systematic error? Is the error due to the limitations of your apparatus? Does your data look the same to within a standard deviation? Evaluate the statistical significance of your data (click here to review the statistical treatment of data). After validating your data, you should interpret your results; state what you believe your results mean. How do your results help us understand the scientific problem? What do your results mean in the context of the bigger picture of chemistry, or of science? How do your results relate to the concepts outlined in the introduction? Do not assume that your experiment failed or was successful. You need to prove to the reader, with logical arguments and supporting evidence, the value of your study.

The conclusions that you wrote in your laboratory notebook are a good starting point from which to organize your thoughts. Your paper's discussion section is structured very similarly to the conclusions section in your notebook, and it might be good idea to review that now (see the [The Laboratory Notebook](#) section to review the structure of the conclusions in the laboratory notebook).

Conclusions

The *Conclusions* section is typically a one-paragraph summary of your laboratory report. Here you summarize the goal(s) of your experiment, state whether you reached that goal, and describe briefly the implications of your study. Note that in some chemistry sub-disciplines it is acceptable to combine the *Discussion and Conclusions* sections. Consult your course syllabus or check with your instructor on the specific format to be used in your class.

Acknowledgements

The *Acknowledgements* section is where you thank anyone who helped you significantly with the project or with the manuscript. For instance, you would thank your laboratory partners if they're not authors on the paper, anyone who helped with the design of the experiment or the preparation of the paper. You might also include funding sources such as a Truman State University summer scholarship or a National Institutes of Health grant.

References

Most of the ideas presented in your paper are probably not exclusively yours. Therefore, you should cite other people's work wherever appropriate. However, you do not need to cite information that is common knowledge or is exclusively your idea. The References section is a compilation of all citations made within the paper. It is not a bibliography and therefore should not list sources that are not directly referred to in the text.

References Format

The format of references varies amongst journals. For your chemistry laboratory reports, you should follow, by default, the ACS guidelines as outlined in The ACS Style Guide and Journal of the American Chemical Society, JACS (all examples given in this document conform to JACS format). If your professor requires you to conform to a specific journal's format, look at articles from that journal or refer to the journal's "Instructions to Authors." The specifications for most ACS journals are:

- 1) References should be compiled at the end of the paper in the References section.
- 2) References should be numbered in the order that they appear in the paper. For citations in the narrative, numbers should be superscripted and appear after the punctuation mark.
- 3) No empty lines should be inserted between reference entries.
- 4) This section should be double spaced just like the rest of your paper.
- 5) A reference is only listed once in the References section. If multiple citations of the reference are made in the text, then the number corresponding to that reference is placed in the text each time. The common abbreviations used in footnotes and references (e. g., *op. cit.*, *ibid.*) are not generally used in scientific writing.

Types of References

Articles. Journal articles are the primary source found in laboratory reports. An example is given below. Notice that the authors' initials are given instead of the first and middle names. Also, there is no "and" before the last author's name. Some journals require that the article's title be included in the reference (check with your instructor to see if he/she wants you to use this style). When included, the article's title should start with a capital letter but the other words in the title, unless they are proper nouns, should not be capitalized (see below). The journal title is abbreviated (click [here](#) for a list of the ACS abbreviations for common journals). Also, the year and the comma after the year are in bold. Lastly, the reference has inclusive pagination (first and last pages are given).

The following are examples of the same journal article with the first given in style where the article's title is included in the reference, while the second is in the style where the article's title is omitted.

(1) Bergmann, U.; Glatzel, P.; deGroot, F.; Cramer, S. P. High resolution K capture X-Ray fluorescence spectroscopy: a new tool for chemical characterization. *J. Am. Chem. Soc.* **1999**, *121*, 4926-4927.

(1) Bergmann, U.; Glatzel, P.; deGroot, F.; Cramer, S. P. *J. Am. Chem. Soc.* **1999**, *121*, 4926-4927.

Books. Books should be cited in the following manner:

(2) Brünger, A. T. *X-PLOR Manual, Version 3.1: A System for X-ray Crystallography and NMR*; Yale University: New Haven, CT, 1990; pp 187-206.

(3) Cheatham, T. E., III; Kollman, P. A. *In Structure, Motion, Interaction, and Expression of Biological Macromolecules*; Sarma, R. H. and Sarma, M. H., Eds, Adenine: New York, 1998; p. 99.

Computer Programs. Citations for computer programs vary. If a person in academia wrote the program, there is often a journal-article source. In other cases, the program is simply distributed by a company.

Journal Article

(4) Humphrey, W.; Dalke, A.; Schulten, K. VMD: Visual Molecular Dynamics. *J. Mol. Graph.* **1996**, *14*, 33-38.

Company Distribution

(5) Case, D. A.; Pearlman, D. A.; Cladwell, J. W.; Cheatham, T. E.; Ross, W. S.; Simmerling, C. L.; Darden, T. A.; Merz, K. M.; Stanton, R. V.; Cheng, A. L.; Vincent, J. J.; Crowley, M.; Ferguson, D. M.; Radmer, R. J.; Seibel, G. L.; Singh, U. C.; Weiner, P. K.; Kollman, P. A. AMBER version 5.0; University of California: San Francisco, 1997.

(6) Insight II; San Diego, CA: Molecular Simulations, 2000.

Websites. Journal articles are much preferred over websites. Websites are dynamic and are usually not peer reviewed. One of the only instances when a website is an acceptable reference is when it is referring to a database (however, an article is usually associated with the creation of the database). If you must use a website, the reference should include a title for the site, the author(s), year of last update and URL. It is unacceptable to use a website as a reference for scientific data or explanations of chemical processes.

(7) Cheatham, T. E., III Simulation Protocol for Polynucleotides; 1998,
<http://www.amber.ucsf.edu/amber/tutorial/polyA-polyT>.

Tables, Schemes and Figures

Tables, schemes and figures are all concise ways to convey your message. As you prepare these items for your report, remember to think of your reader. You want them to derive the maximum amount of information with the minimum amount of work. Pretend to be the reader and ask yourself, “Does this enhance my understanding?”, “Can I find everything?”, “Can I read it without being distracted?” Poorly prepared tables, schemes and figures will reflect badly on your science and you as a scientist, so think carefully about these items as you prepare your report.

Tables

A table is a way to summarize data or ideas in a coherent, grid-like fashion. This is not simply output from a spreadsheet! You should prepare the table in a word-processor so that its formatting matches the rest of your report. In general, tables have no more than ten rows and columns to avoid overwhelming the reader. One common exception is in review articles (such as in Chemical Reviews) where an author is attempting to summarize results from an entire field. Another common exception is in the reporting of X-ray crystallography data. These tables have their own special formatting rules, and will not be discussed here.

Tables are referred to in the text as “Table #”. Tables, schemes and figures are labeled separately, with Arabic numbers, in the order they are referred to in the paper. Tables have a table caption, which in some journals appears above the table, while in others it appears below. In either case, the table caption is always on the same page as the table.

Don’t use lines or boxes in your table except where absolutely necessary. Use spaces between your columns instead (helpful hint: it is better to use your word processor’s table formatting tools than trying to get the columns to line up using tabs or spaces). All column or row headings should have clear subtitles and units if needed (usually in parentheses). Any numbers that are presented should have proper significant figures, and an indication of the error should be shown (please see the [Introduction to Statistics in Chemistry](#) for information on how to report uncertainty in one’s data). An example table is given below.

Table 1. Aminoacylation efficiency of duplexAla substrates containing base pair substitutions at the 2:71 position.

2:71 Base Pair	k_{cat}/K_M (relative) ^a	Fold decrease ^b	$-\Delta\Delta G^\ddagger$ (kcal/mol) ^c
G:C (wild-type)	1		0
<u>Watson-Crick Pur:Pyr Base Pairs</u>			
I:C	0.51	1.9	0.39
G:4HC	0.25	3.9	0.81
2AA:U	0.23	4.3	0.86
2AP:U	0.18	5.6	1.0

^aValues reported are averages of at least three determinations with average standard deviations of $\pm 26\%$.

^bFold decrease in k_{cat}/K_M is given relative to wild-type duplexAla.

^c $\Delta\Delta G^\ddagger$ is defined as $RT\ln[(k_{cat}/K_M)^{\text{variant}}/(k_{cat}/K_M)^{\text{wild-type}}]$, where $R=1.98272$ cal/mol·K and $T=298$ K.

Schemes

A scheme is usually a sequence of two or more chemical reactions that together summarize a synthesis. A scheme may also show the steps in a purification with each step or reaction giving the reactants, products, catalysts, and yields. A scheme that shows a chemical reaction may also


show possible intermediates. Note that mechanisms are not usually conveyed using a scheme because they are more complicated and illustrate where electrons are proposed to move. Mechanisms are most often placed in a figure.

It is a common convention in a scheme to write a bold number underneath chemical species referred to in the text. Note that for the first occurrence of the bold number in the text, the chemical's name is given, but after that only the bold number is used to identify it. This method of defining abbreviations for compounds can also be done in the experimental section, if there is no scheme. This is very useful when a compound's name is long or complicated.

The one-step yield is usually written to the right of the equation, although it is also proper to write the yield under the arrow. Note also how the reaction conditions can be summarized (i. e., the first step below), which saves the reader from flipping to the experimental section for these details.

Each scheme also has a caption, which is included under the scheme. The caption should briefly summarize what is in the scheme. If the scheme is from another source, the reference to this source should appear at the end of the caption.

The following is an example of a scheme that might appear in a synthetic paper. The text below it shows how the scheme could be referred to in the body of the paper.



Scheme 1. Synthesis of benzoyl chloride (3).

Benzamide (1) was refluxed under aqueous acidic conditions for 1 hour to yield benzoic acid (2). Acid (2) was then refluxed with SOCl_2 to yield benzoyl chloride (3).

Sometimes a scheme may be used to illustrate a non-chemical process or how an instrument's components are connected. These could also be presented as figures, and there is no definitive rule that will tell you when to use a scheme and when to use a figure. When in doubt, think of the reader and use the method that conveys the most information in the most easily understood format

Figures

Figures fall into two broad categories; those that are pictorial representations of concepts that are presented in the text, and those which summarize data. Again, it is critical to your report that your figures are clear, concise and readable, and that they support the arguments that you are making. Remember that you must refer to and discuss every figure in the text! If a figure is not mentioned, you don't need it!

Figures that are pictorial representations of concepts usually appear in the Introduction, but it is also appropriate to include them in the Discussion. Use this type of figure to make your writing more concise (remember the conversion factor: 1 picture = 1 keyword). Remember, humans are very visually oriented and we can grasp complex concepts presented as picture more easily than when they are presented in words or as mathematical formulae. Some examples of concept figures include:

- 1) An illustration of the deposition of metals onto a silicon wafer.
- 2) A diagram of the HIV life cycle.
- 3) A depiction of microwaves exciting water molecules.
- 4) A diagram illustrating the Frank-Condon principle.
- 5) A proposed organic mechanism.

Graphs are figures that present data. You use a graph when you have more data than will fit in a table. The general rules for preparing good figures for your notebook also apply in a laboratory report (click here to review graph preparation). Formatting tips: do not use colored backgrounds or gridlines, and do not draw a box around the graph.

You may find it more concise to combine all your data into one graph. For example, it may be appropriate to put six lines with absorbance as a function of time, with varying concentrations of a reactant on the same graph rather than constructing six different graphs. However, when doing this, be careful not to over-clutter the graph.

Standard curves should not be included in this section unless that was the primary goal of the experiment. They should be put in the Supporting Information.

Figures have figure captions compiled in the *Figure Legend* section, located on a separate page at the end of the paper. Journals chose this format because of typographical issues, and it has been retained despite its inconvenience to the reader. Each figure should appear on its own page in the order it is discussed in the text. Figure captions appear in the *Figure Legends* section and do not appear on the same page as the figure. However, in the bottom, right-hand corner of the page the following identifying text appears:

"First author's last name et al., Figure number"

Figure Legends

All figure legends (captions) should be found in the section entitled "Figure Legends". The format for a figure legend is usually: "Figure number" (italics and bold), a short title (followed by a

period) and then a description of what is in the figure. All figure legends are compiled on the same page separated by a blank line. Be sure to define in the caption any symbols used in the figure, and note whether lines that pass through data points are fits, or “guides to the eye”.

Example Figure Caption

Figure 1. Nucleic acid bases. The chemical structures of (a) adenine, (b) guanine, (c) cytosine, and (d) thymine.

Supporting Information

This section (also known as Supplemental Material) is where you can include information that may be helpful, but not essential, for evaluation of your data. Items in this section may include calibration curves, and spectra (from which you extracted only one absorbance value for your analysis). Figures or tables of data whose contents were summarized in the text, or which were not critical to the conclusions, are also to be placed in the supporting information. An example of this type of material is the table of atom positions generated in an X-ray crystal structure.

References

1. *The ACS Style Guide*; 2nd Ed.; Dodd, J. S., Ed.; American Chemical Society: Washington, D.C., 1997.
2. Booth, W. C.; Colomb, G. G.; Williams, J. M. *The Craft of Research*; The University of Chicago Press: Chicago, IL, 1995.
3. Spector, T. J. *J. Chem. Educ.* **1994**, 71, 47-50. This article is available as a PDF file at <http://pubs.acs.org/doi/pdfplus/10.1021/ed071p47> for Truman addresses and *J. Chem. Educ.* subscribers.
4. The instructions for authors for several American Chemical Society journals may be found at <http://pubs.acs.org/page/jacsat/submission/authors.html> (Journal of the American Chemical Society); <http://pubs.acs.org/page/inocaj/submission/authors.html> (Inorganic Chemistry) and <http://pubs.acs.org/page/chreay/submission/authors.html> (Chemical Reviews), all last accessed December 10, 2015.
5. Any non-English word should be italicized. This includes Greek and German words, and their abbreviations, that appear as part of chemical names (e. g., *ortho*-, *meta*-, *para*-, *cis*-, *trans*-, *E*-, *Z*-, *alpha*-, *beta*-, etc.). Also italicized are the condensed forms of secondary (*sec*-), tertiary (*tert*-), etc. The primary exception to the rule for italicizing non-English words are the Greek and Latin prefixes that denote numbers in chemical names (e. g., mono-, bi-, tri-, etc.). Some common Latin phrases that appear in scientific writing are *vide infra* (“see later”), *vide supra* (“see earlier”), *et al.* (abbreviation of *et alia*, Latin for “and others”), e. g. (from Latin *exempli gratia*, “for example”, not usually italicized) and *i. e.* (from Latin *id est*, “it is”, also not usually italicized). Other Latin phrases and abbreviations commonly used in footnotes and references (e. g., *op. cit.*) are not used in scientific writing.

Example Formal Report

Synthesis and Characterization of Luminol

Maria C. Nagan, Eric V. Patterson and James M. McCormick

Science Division, Truman State University, 100 East Normal, Kirksville, MO 63501

Abstract

Luminol (5-amino-2,3,-dihydrophthalazine-1,4-dione) was synthesized and its chemiluminescent properties were explored. Both spectral and physical data indicate luminol was successfully synthesized but in low yield, presumably due to entropically unfavorable reactions. A more reactive starting material such as an acid chloride may help overcome entropic barriers.

Chemiluminescence of luminol, the excitation and emission of light, was observed in the presence of H_2O_2 and an oxidant enhancer. The chemiluminescent quantum yield of luminol, which gives an indication of luminescent intensity, was determined to be moderate (0.10) when compared to other chemiluminescent compounds.

Introduction

The chemiluminescence of luminol, 5-amino-2,3-dihydrophthalazine-1,4-dione, was first discovered by Albrecht in 1928.¹ Since then, the luminescent properties of luminol and other hydrazines have been thoroughly investigated^{2,3} and have found applications in such diverse areas as the detection of ion concentrations in aqueous solutions,⁴⁻⁶ monitoring the progress of H₂O₂-dependent reactions⁷ and even the detection of blood at crime scenes.⁸ In the current study, luminol, one of the more widely studied chemiluminescent molecules, is synthesized and its spectroscopic properties are examined.

The mechanism of luminol light emission and how it can be enhanced is still not well understood. However, a few pioneering studies have explored the chemiluminescent properties of luminol in both aprotic and protic solvents. In aprotic solvents such as dimethylsulfoxide, luminol has been shown to chemiluminesce via an excited 3-aminophthalate (**2**) dicarboxylate anion (Scheme 1).^{9,10} The reaction requires two equivalents of base, such as NaOH, to deprotonate the amide nitrogens before subsequent reaction with O₂ gives **2**. In protic solvents, the reaction is more complex and the chemiluminescent reaction involves H₂O₂, whose exact role is undetermined, and some additional oxidant such as potassium ferricyanide in addition to luminol.² The emission of light requires both H₂O₂ and the oxidant enhancer to produce **2**.^{2,11} A free-radical mechanism has been proposed,¹² but the mechanism of oxidant enhancement, increased luminescence in the presence of the oxidant, is still not well understood. Enhanced luminescence can be achieved with a number of different ferric enzymes including horseradish peroxidase and luciferin.^{13,14}

Chemiluminescence is distinct from fluorescence and phosphorescence in that a chemical is used to excite a molecule and an external source of photons is not needed. In chemiluminescence, a molecule is excited through a set of chemical reactions. The excited molecule, usually in the lowest excited singlet state,¹⁰ can spontaneously lose its electronic energy in the form of a photon. If the electronic spin ($\Delta S = 0$) is preserved, this process is referred to as fluorescence.¹⁵

$$I_{CL} = \phi_{CL} \frac{dC}{dt} \quad (1)$$

The radiant intensity of chemiluminescence, I_{CL} , can be expressed as¹⁶ and is dependent upon both the rate of chemical reaction, dC/dt , and the chemiluminescent quantum yield ϕ_{CL} . The chemiluminescent quantum yield, the number of photons emitted per chemiluminescing molecule, provides a means to quantify the efficiency and degree of chemiluminescence.

In this study, luminol was synthesized and the chemiluminescent properties of it were explored. The luminescence was first qualitatively observed in a hydrogen peroxide-dependent reaction, using a Cu(II) species as the oxidant. Then the chemiluminescent emission spectrum of luminol in aqueous solution was obtained to determine the chemiluminescent quantum yield of luminol in water and the extent of chemiluminescence.

Experimental

Melting point data was obtained on a Thomas-Hoover apparatus. Infrared (IR) spectra were collected on a Beckman Coulter IR spectrophotometer. Proton and carbon nuclear magnetic resonance (NMR) spectra were recorded on a Bruker Avance 400 MHz spectrometer. All chemicals were obtained from Sigma (St. Louis, MO) and used without further purification.

5-amino-2,3-dihydrophthalazine-1,4-dione. Luminol was synthesized from 3-nitrophthalic acid (Scheme 2) according to previously published procedures.¹⁷ In a test tube, 3-nitrophthalic acid (1.0024 g, 4.748 mmol) and 8 % aqueous hydrazine (2.05 mL) were combined. Triethyleneglycol (3.02 mL) was then added to the mixture and excess water was distilled away by heating. The mixture was refluxed at 215 °C for 2 minutes, and quenched with hot water (15 mL). Cooling the reaction mixture in an ice bath afforded a precipitate, which was collected by vacuum filtration. To the isolated solid was added 3 M sodium hydroxide (5 mL) and then sodium hydrosulfite (2.5100 g). The mixture was refluxed for 5 minutes at 180 °C and quenched with glacial acetic acid (2 mL). Upon cooling on ice, a yellow solid was collected and recrystallized from 95% ethanol: 35% yield; m.p. 326-334 °C; IR (KBr) 1496 (s), 1603 (s), 1667 (s), 3029 (s-broad), 3442 (m), 3464 (m) cm⁻¹; ¹H NMR (400 MHz, DMSO-d₆) δ 6.939 (m;1), 6.999 (d;1), 7.360 (d; 2), 7.490 (m; 1), 11.200 (s;2) ppm.

Luminescence. To a solution (250 mL) containing 0.3 M sodium bicarbonate, 4 mM anhydrous ammonium carbonate, 1.8 mM copper (II) sulfate, 37.5 mM sodium carbonate, product (0.1002 g) from the above synthesis was added, followed by an equal volume to a 0.15% aqueous hydrogen peroxide.

Luminescence Spectroscopy. The chemiluminescence emission spectrum of 99.5% pure luminol (Sigma) was obtained using a Perkin-Elmer luminescence spectrometer (model LS-5B). The previously determined maximum emission wavelength of luminol in water (424 nm) was used in the present study.¹⁰ The chemiluminescence spectra for luminol were obtained for three different samples (25 mM) in water (Figure 1).

Results

Luminescence. A blue fluorescence was observed upon combination of the above solutions.

Luminescence Spectroscopy. From luminescent spectra of luminol (Fig. 1), the chemiluminescent quantum yield was determined to be 0.10 ± 0.03 .

Discussion

Both physical and spectral data indicate that luminol was indeed synthesized. The IR spectrum (See Supporting Information) shows the characteristic amide vibrational stretch at about 3400 cm⁻¹. No carboxylic acid stretch (1700 cm⁻¹) or nitro group stretch (1550 cm⁻¹) was present in the IR spectrum but instead an amine at 1603 cm⁻¹, indicating successful conversion of the starting material to the desired product. The proton NMR spectrum (See Supporting Information) indicates that both aromatic amine (δ = 7.360 ppm; 2) and amide protons (δ = 11.200 ppm; 2) are present in addition to three aromatic protons (δ = 6-8 ppm). Lastly, the yellow solid's melting point of 326-334 °C, although a little low and the range rather broad, is similar to the literature melting point of luminol (332 °C).¹⁷

The 35 % overall yield for the reaction was low. This may be due to an intrinsically low reaction-efficiency, which is expected since the overall reaction is not entropically favorable. The reaction requires two consecutive condensations, the first intermolecular and the second an intramolecular cyclization, both of which are likely to have $\Delta S < 0$. With a high reaction temperature and a net negative ΔS , the overall free energy of reaction is therefore less favorable. A modification of the

synthesis, by replacing the carboxylic acid groups of 3-nitrophthalic acid with more reactive group (e.g. acid chlorides) may increase overall yield. Lastly, the acidity of the solution is a critical component because the carbonyls need to be activated in the first step of the reaction (Scheme 1) but too many protons may lead to protonation of hydrazine. As this reaction is conducted in one pot, careful control of the pH in future experiments may increase yield.

The chemiluminescence of luminol was examined both in a qualitative and quantitative manner. Luminol chemiluminesced blue in the presence of H₂O₂ and Cu(II). In addition, the chemiluminescent quantum yield of luminol was determined (0.10 ± 0.03) and shown to be moderate relative to most chemiluminescent compounds which have quantum yields in the range of 0.01 to 0.2.¹⁶ Other derivatives of luminol have been found to luminesce more intensely,² specifically, changing the 5-amino substituent to an alkoxide or an alkyl amino group has been shown to increase light emission.^{18,19} As well, enhancers such as ferric enzymes greatly increase luminescence intensity.^{13,14} It is therefore not surprising that many applications of chemiluminescence employ either a luminol derivative or an oxidant enhancer.

Conclusion

Luminol has been successfully synthesized in pure form. The yield of 35% was rather low but expected since the reaction is not entropically favorable. The chemiluminescent quantum yield was determined to be 0.10, which is indicative of a rather moderate luminescence intensity.

Supporting Information Available


Full synthetic details and characterization are available in this section (PDF). The material is available free of charge via the Internet at <http://pubs.acs.org>.

References

- (1) Albrecht, H. O. *Z. Phys. Chem. (Leipzig)* **1928**, *136*, 321.
- (2) White, E. H.; Roswell, D. F. *Acc. Chem. Res.* **1970**, *3*, 54-62.
- (3) Gundermann, K. D. *Chimia* **1971**, *25*, 261.
- (4) Kurihara, M.; Hasebe, T. and Kawashima, T. *Bunseki Kagaku* **2002**, *51*, 205-233.
- (5) Meseguer-Lloret, S.; Campins-Falco, P.; Tortajada-Genaro, L. A. and Blasco-Gomez, F. *Int. J. Env. Anal. Chem.* **2003**, *83*, 405-416.
- (6) Nieman, T. In *Chemiluminescence and Photochemical Reaction Detection in Chromatography*; Birks, J. W., Ed.; VCH: New York, 1989, pp 99-123.
- (7) Yuan, J. and Shiller, A. M. *Anal. Chem.* **1999**, *71*, 1975-1980.
- (8) Kent, E. J. M.; Elliot, D. A. and Miskelly, G. M. *J. Foren. Sci.* **2003**, *48*, 64-67.
- (9) White, E. H.; Zafirliou, O.; Kagi, H. H. and Hill, J. H. M. *J. Am. Chem. Soc.* **1964**, *86*, 940-941.
- (10) White, E. H. and Bursey, M. M. *J. Am. Chem. Soc.* **1964**, *86*, 941-942.

- (11) Lind, J. S.; Merényi, G. and Eriksen, T. E. *J. Am. Chem. Soc.* **1983**, *105*, 7655-7661.
- (12) White, E. H. In *Light and Life*; McElroy, W. D. and Glass, B., Eds.; Johns Hopkins Press: Baltimore, MD, 1961, p 183.
- (13) Whitehead, T. P.; Thorpe, G. H. G.; Carter, T. J. N.; Groucutt, C. and Kricka, L. J. *Nature* **1986**, *305*, 185-189.
- (14) Easton, P. M.; Simmonds, A. C.; Rakishev, A.; Egorov, A. M. and Candeias, L. P. *J. Am. Chem. Soc.* **1996**, *118*, 6619-6624.
- (15) Levine, I. N. *Physical Chemistry*, 3rd Ed.; McGraw-Hill: New York, 1988.
- (16) Skoog, D. A.; Holler, F. J. and Nieman, T. A. *Principles of Instrumental Analysis*, 5th Ed.; Thomson Learning, 1998.
- (17) Williamson, K. L. In *Macroscale and Microscale Organic Experiments*, 2nd Ed.; D.C. Heath: Lexington, MA, 1994, pp 153-162.
- (18) White, E. H. and Bursey, M. M. *J. Org. Chem.* **1966**, *31*, 1912.
- (19) Gundermann, K. D.; Horstmann, W. and Bergmann, G. *Ann.* **1965**, *684*, 127.

Scheme 1. Chemiluminescence of luminol (5-amino-2,3-dihydrophthalazine) (**1**) in aprotic solvents.



Scheme 2. Preparation of luminol (5-amino-2,3-dihydrophthalazine-1,4-dione) (**1**) from 3-nitrophthalic acid (**3**) via a nitrohydrazine intermediate (**4**).




Figure Legends

Figure 1. Chemiluminescence emission spectra of luminol. Spectra of 25 mM luminol in water were obtained on a Perkin-Elmer luminescence spectrometer (LS-5B). Spectra were obtained in triplicate and each spectrum is indicated by a separate line (—, - - - or ...).

Luminescence Spectrum goes here.
All axes would be appropriately labeled with units.

Nagan *et al.*, Figure 1